# C. DIFF QUIK CHEK®

A rapid test for detection of *C. difficile* glutamate dehydrogenase in fecal specimens

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International Symbol Key:

Catalog Number

*In Vitro* Diagnostic Medical Device

Lot Information

Contains sufficient reagents  
for <n> tests

Temperature Limitation

Use By/Expiration Date

CE Symbol

Caution, consult  
accompanying documents

## C. DIFF QUIK CHEK®

### INTENDED USE

The *C. DIFF QUIK CHEK*® test is a rapid membrane enzyme immunoassay for use as a screening test to detect *Clostridium difficile* antigen, glutamate dehydrogenase, in fecal specimens from persons suspected of having *C. difficile* disease. The test does not distinguish toxigenic from nontoxigenic strains of *C. difficile*. With the use of additional tests that detect *C. difficile* toxins, the test is to be used as an aid in the diagnosis of *C. difficile* disease. As with other *C. difficile* tests, results should be considered in conjunction with the patient history.  
FOR *IN VITRO* DIAGNOSTIC USE.

### EXPLANATION

After treatment with antibiotics, many patients develop gastrointestinal problems ranging from mild diarrhea to severe pseudomembranous colitis. Many cases of the milder forms of gastrointestinal illness and most cases of pseudomembranous colitis are caused by toxigenic strains of *Clostridium difficile* (1). This organism is an opportunistic anaerobic bacterium that grows in the intestine once the normal flora has been altered by the antibiotic. Toxigenic strains of *C. difficile* carry the genes encoding the toxins while non-toxigenic strains do not carry the toxin genes. The disease results from the toxins that the toxigenic organism produces. The clinical symptoms associated with the disease are believed to be primarily due to toxin A, which is a tissue-damaging enterotoxin (2,3). *C. difficile* also produces a second toxin, designated toxin B. Toxin B, which has been referred to as the cytotoxin of the organism, is the toxin detected by the tissue culture assay currently used by some laboratories. Toxigenic *C. difficile* strains produce both toxins or only toxin B (4-7). The glutamate dehydrogenase of *C. difficile* is a good antigen marker for the organism in feces because it is produced in high amounts by all strains, toxigenic or non-toxigenic (8-10). The antigen can be detected in fecal specimens by using the *C. DIFF QUIK CHEK*® test. A positive result in the test, which is specific for the glutamate dehydrogenase of *C. difficile*, confirms the presence of this organism in a fecal specimen; a negative result indicates the absence of the organism. A positive result should be followed by a toxin-specific test to confirm the presence of toxigenic *C. difficile*.

### PRINCIPLE OF THE TEST

The *C. DIFF QUIK CHEK*® test uses antibodies specific for the glutamate dehydrogenase of *C. difficile*. The device contains a *Reaction Window* with two vertical lines of immobilized antibodies (Fig. 1a). The test line ("T") contains antibodies against *C. difficile* glutamate dehydrogenase. The control line ("C") contains anti-IgG antibodies. The *Conjugate* consists of antibodies to glutamate dehydrogenase coupled to horseradish peroxidase. To perform the test, the sample is added to a tube containing a mixture of *Diluent* and *Conjugate*. The diluted sample-conjugate mixture is added to the *Sample Well* and the device is allowed to incubate at room temperature for 15 minutes. During the incubation, any glutamate dehydrogenase in the sample binds to the antibody-peroxidase conjugate. The antigen-antibody-conjugate complexes migrate through a filter pad to a membrane where they are captured by the immobilized glutamate dehydrogenase-specific antibodies in the line. The *Reaction Window* is subsequently washed with *Wash Buffer*, followed by the addition of *Substrate*. After a 10 minute incubation period, the "T" reaction is examined visually for the appearance of a vertical blue line on the "T" side of the *Reaction Window*. A blue line indicates a positive test. A positive "C" reaction, indicated by a vertical blue line on the "C" side of the *Reaction Window*, confirms that the test is working properly and the results are valid.

### MATERIALS PROVIDED

**MEM | DEV** *Membrane Devices* – 25 pouches, each containing 1 device

DIL   SPE	<b>Diluent (14 mL)</b> – Buffered protein solution with graduated dropper assembly
WASH   REAG	<b>Wash Buffer (12 mL)</b> – A buffered solution with graduated dropper assembly
SUBS   REAG	<b>Substrate (3.5 mL)</b> – Solution containing tetramethylbenzidine
CONJ   ENZ	<b>Conjugate (2.0 mL)</b> – Mouse monoclonal antibody specific for glutamate dehydrogenase coupled to horseradish peroxidase in a buffered protein solution
CONTROL   +	<b>Positive Control (1mL)</b> – Antigen in a buffered protein solution

**Disposable plastic transfer pipettes**, 50 (graduated at 25  $\mu$ L, 400  $\mu$ L and 500  $\mu$ L)

#### MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

<i>Small test tubes (e.g., plastic Eppendorf tubes)</i>	<i>Applicator sticks</i>	<i>Timer</i>
<i>Disposable gloves for handling fecal samples</i>	<i>Pipettor and tips</i>	<i>Vortex mixer</i>

#### SHELF LIFE AND STORAGE

The expiration date of the kit is given on the label. Expiration dates for each component are listed on the individual labels. The kit should be stored between 2° and 8°C.

#### PRECAUTIONS

1. Reagents from different kits should not be mixed. Do not use a kit past the expiration date.
2. Bring all components to ROOM TEMPERATURE BEFORE USE!
3. Caps, tips and dropper assemblies are color-coded; do NOT mix or interchange!
4. Do not freeze the reagents. The kit should be stored between 2° and 8°C.
5. The pouch containing the *Membrane Device* should be at room temperature before opening, and opened just before use. Keep the membrane devices dry before use.
6. Use fecal specimens within 72 hours of collection to obtain optimal results. Specimens that are frozen may lose activity due to freezing and thawing.
7. Specimens that have been preserved in 10% Formalin, merthiolate formalin, sodium acetate formalin or polyvinyl alcohol cannot be used.
8. Specimens in transport media such as Cary Blair and C&S can be used as specified in the specimen preparation protocol.
9. Hold reagent bottles vertically to dispense reagents to ensure consistent drop size.
10. Specimens and membrane devices should be handled and disposed of as potential biohazards after use. Wear disposable gloves when doing the test.
11. Membrane devices cannot be reused.
12. The test has been optimized for sensitivity and specificity. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test. Do not deviate from the specified procedure.
13. Microbial contamination of reagents may decrease the accuracy of the assay. Avoid microbial contamination of reagents by using sterile disposable pipettes if removing aliquots from reagent bottles.
14. Be attentive to the total assay time when testing more than one fecal specimen. Add *Diluent* first, then add the *Conjugate* to each tube of *Diluent*. Then add specimen to the tube of *Diluent/Conjugate*. Thoroughly mix all of the diluted specimens, and then transfer to the *Membrane Device*. The 15-minute incubation step begins after the last diluted sample-conjugate mixture has been transferred to the final *Membrane Device*.
15. If the *Substrate* reagent changes to a dark blue/violet color call technical services for replacement.

#### COLLECTION AND HANDLING OF FECAL SPECIMENS


1. Standard collection and handling procedures used in-house for fecal specimens are appropriate. Specimens should be stored between 2° and 8°C; test specimens that are less than 24 hours old, whenever possible.

2. Store specimens frozen ( $\leq -10^{\circ}\text{C}$ ) if the test cannot be performed within 72 hours of collection, but note that freezing and thawing of the specimen may result in loss of activity due to degradation of the antigen. Do not freeze and thaw multiple times.
3. Make sure that specimens are thoroughly mixed PRIOR to performing the assay.
4. Storing fecal specimens in the *Diluent* is NOT recommended.
5. Do not allow the fecal specimens to remain in the *Diluent* and *Conjugate* for more than 8 hours.

### SAMPLE PREPARATION

1. Bring all reagents and the required number of devices to room temperature before use.
2. Set up and label one small tube for each specimen, and optional external controls as necessary.
3. Add 500  $\mu\text{L}$  *Diluent* to each tube for fecal specimens using the graduated black dropper assembly (or equivalent). For specimens in transport media such as Cary Blair or C&S, add 425  $\mu\text{L}$  of *Diluent* to the tube.
4. Holding the dropper bottle in a vertical position, add one drop of *Conjugate* (red capped bottle) to each tube.
5. Obtain one disposable plastic transfer pipette (supplied with the kit) for each sample the pipettes have raised graduations at 25  $\mu\text{L}$ , 400  $\mu\text{L}$  and 500  $\mu\text{L}$ .

#### Graduated Transfer Pipette:

6. Mix all specimens thoroughly regardless of consistency- it is essential that the specimens be evenly suspended before transferring.

**Liquid/Semi-solid specimens** – pipette 25  $\mu\text{L}$  of specimen with a transfer pipette (graduated at 25  $\mu\text{L}$ , 400  $\mu\text{L}$  and 500  $\mu\text{L}$ ) and dispense into the *Diluent/Conjugate* mixture. Use the same transfer pipette to mix the diluted specimen.

**Formed/Solid specimens** – Care must be taken to add the correct amount of formed fecal specimen to the sample mixture. Mix the specimen thoroughly using a wooden applicator stick and transfer a small portion (approximately 2 mm diameter, the equivalent of 25  $\mu\text{L}$ ) of the specimen into the *Diluent/Conjugate* mixture. Emulsify the specimen using the applicator stick.

**Fecal specimens in Cary Blair or C&S transport media** - pipette 100  $\mu\text{L}$  of sample into the *Diluent/Conjugate* mixture.

7. **Optional External Control Samples:**

**External Positive Control** – holding the dropper bottle in a vertical position, add one drop of *Positive Control* (gray-capped bottle) to the appropriate test tube.

**External Negative Control** - add 25  $\mu\text{L}$  *Diluent* to the appropriate test tube using a graduated transfer pipette.

**NOTE:** Transferring too little specimen, or failure to mix and completely suspend the specimen in the *Diluent* mixture, may result in a false-negative test result. The addition of too much fecal specimen may cause invalid results due to restricted sample flow.

### TEST PROCEDURE

1. Obtain one *Membrane Device* per specimen, and one device per optional external positive or negative control as necessary. The foil bags containing the devices should be brought to room temperature before opening. Label each device appropriately and orient it on a flat surface so the letter “C” on the device is on the left, the letter “T” is on the right, and the small *Sample Well* is located in the top right corner of the device (Fig. 1a).

2. Close each tube of diluted specimen and mix thoroughly. Proper mixing can be achieved by vortexing or inverting the tube. Immediately proceed to Step #3.
3. Using a transfer pipette (graduated at 25  $\mu\text{L}$ , 400  $\mu\text{L}$  and 500  $\mu\text{L}$ ), transfer 400  $\mu\text{L}$  of the diluted sample-conjugate mixture into the **Sample Well** (smaller hole in the top right corner of the device) of a *Membrane Device*, making certain to expel the liquid sample onto the wicking pad inside of the *Membrane Device*.
4. Incubate the device at room temperature for 15 minutes – the sample will wick through the device and a wet area will spread across the *Reaction Window* (larger hole in the middle of the device).

**NOTE FOR SAMPLES THAT FAIL TO MIGRATE:**

*Occasionally, a diluted fecal specimen cannot be tested because it clogs the membrane and the Reaction Window does not wet properly. If the diluted fecal specimen fails to migrate properly within 5 minutes of adding the sample to the Sample Well (i.e. the membrane in the Reaction Window does not appear to be completely wet), then add 100  $\mu\text{L}$  (4 drops) of Diluent to the Sample Well and wait an additional 5 minutes (for a total of 20 minutes).*


5. After the incubation, add 300  $\mu\text{L}$  of *Wash Buffer* to the **Reaction Window** using the graduated white dropper assembly (or equivalent). Allow the *Wash Buffer* to flow through the *Reaction Window* membrane and be absorbed completely.
6. Holding the dropper bottle vertically, add 2 drops of *Substrate* (blue-capped bottle) to the **Reaction Window**. Read and record results visually after 10 minutes.

**INTERPRETATION OF RESULTS**

1. Interpretation of the test is most reliable when the device is read immediately at the end of the reaction period. Read the device at a normal working distance in a well-lit area. View with a line of vision directly over the device.
2. Observe device for the appearance of a blue line on the “C” side of the *Reaction Window* representing the internal positive control line. Observe device for the appearance of a blue line on the “T” side of the *Reaction Window* representing the test line. The lines may appear faint to dark in intensity.
3. **Positive Result:** A positive result may be interpreted at any time between the addition of *Substrate* and the 10-minute read time. Two blue lines are visible, the control line (“C”) and the test line (“T”). The lines may appear faint to dark in intensity. The appearance of a blue line on the “T” side along with a blue control line is interpreted as a positive result. An obvious partial line is interpreted as a positive result. Do not interpret membrane discoloration or shadow as a positive result. A positive result indicates the presence of *C. difficile* antigen.
4. **Negative Result:** A test cannot be interpreted as negative or invalid until 10 minutes following the addition of *Substrate*. A single blue line is visible on the control (“C”) side of the *Reaction Window* and no test line is visible on the “T” side of the *Reaction Window* (Fig. 1b). A negative result indicates *C. difficile* antigen is either absent in the specimen or is below the detection limit of the test.
5. **Invalid Result:** A single line is visible on the test (“T”) side of the *Reaction Window*, or no lines are visible in the *Reaction Window* (Fig. 1c, 1d). The test result is invalid if a control line is not present at the completion of the reaction period.

**Figure 1. C. DIFF QUIK CHEK® INTERPRETATION OF RESULTS**

*Reaction Sample*  
*Window Well*

**Fig. 1a. Positive**

**Fig. 1b. Negative**

**Fig. 1c. Invalid**

**Fig. 1d. Invalid**

## QUALITY CONTROL

**Internal:** A blue control line must be visible on the “C” side of the *Reaction Window* on every *Membrane Device* that is tested. The appearance of the blue control line confirms that the sample and reagents were added correctly, that the reagents were active at the time of performing the assay, and that the sample migrated properly through the *Membrane Device*. A clear background in the result area is considered an internal negative control. If the test has been performed correctly and reagents are working properly, the background will be clear to give a discernible result.

**External:** The reactivity of the *C. DIFF QUIK CHEK*® test should be verified on receipt using the *Positive Control* and negative control (*Diluent*). The *Positive Control* is supplied with the kit (gray-capped bottle). The *Positive Control* confirms the reactivity of the other reagents associated with the assay, and is not intended to ensure precision at the analytical assay cut-off. *Diluent* is used for the negative control. Additional tests can be performed with the controls to meet the requirements of local, state and/or federal regulations and/or accrediting organizations.

## LIMITATIONS

1. The *C. DIFF QUIK CHEK*® test is used to detect *C. difficile* antigen in fecal specimens. The test confirms the presence of *C. difficile* glutamate dehydrogenase in feces and this information should be taken under consideration by the physician in light of the clinical history of the patient and the results of toxin detecting tests. A positive result with the *C. DIFF QUIK CHEK*® test DOES NOT confirm the existence of a TOXIGENIC strain of *C. difficile*. If a positive result is seen using the *C. DIFF QUIK CHEK*® test, it is suggested that additional tests be used to confirm the presence of toxin. Alternatively, tissue culture cytotoxicity assay may be performed to verify the presence of toxin.
2. Fecal specimens are extremely complex. Optimal results with the *C. DIFF QUIK CHEK*® test are obtained with specimens that are less than 24 hours old. Most undiluted specimens can be stored between 2° and 8°C for 72 hours before significant degradation of the antigen is noted. If specimens are not assayed within this time period, they may be frozen and thawed. However, freezing and thawing may result in loss in the immunoreactivity of the antigen.
3. Some specimens may give weak reactions. This may be due to a number of factors such as the presence of low levels of antigen, the presence of antigen-binding substances, or inactivating enzymes in the feces. Under these conditions a fresh specimen should be tested. Additional tests that may be used in conjunction with the *C. DIFF QUIK CHEK*® test include isolation of the organism on selective media, toxin-specific ELISAs, or tissue culture cytotoxicity assays for the detection of *C. difficile* and its toxin. The isolation of the organism does not confirm that the organism is toxigenic. Additional testing of the isolate using a toxin-specific ELISA or tissue culture assay must confirm this.
4. Fecal specimens preserved in 10% formalin, merthiolate formalin, sodium acetate formalin, or polyvinyl alcohol CANNOT be used.
5. The *C. DIFF QUIK CHEK*® test is qualitative. The intensity of the color should not be interpreted quantitatively.
6. Colonization rates of up to 50% have been reported in infants. A high rate has also been reported in cystic fibrosis patients (1,3).
7. The accuracy of the result when testing pediatric patients is uncertain because the performance of *C. DIFF QUIK CHEK*® test of specimens from pediatric patients has been evaluated with a limited number of samples.

## EXPECTED VALUES

*Clostridium difficile* disease is primarily a nosocomial disease of elderly patients, and the frequency of the disease is dependent on factors such as patient population, type of institution and epidemiology. The reported incidence of *C. difficile* disease in

patients with antibiotic-associated diarrhea may range from 5 to 20%, and hospitals may experience rates lower or higher than this range. It is important to consider any test results in conjunction with clinical symptoms because some healthy adults and large numbers of healthy infants (up to 50%) will be positive for *C. difficile* toxin. In addition, *C. difficile* carriage rates of 22 to 32% have been reported in cystic fibrosis patients (1,3). Because the *C. DIFF QUIK CHEK*<sup>®</sup> test detects both toxigenic and non-toxigenic strains of *C. difficile*, the expected values of this test are higher compared to that of any toxin test. A positive result in the *C. DIFF QUIK CHEK*<sup>®</sup> test confirms the presence of *C. difficile* in a fecal specimen; a negative result indicates the absence of the organism. A positive result should be followed by a toxin-specific test to confirm the presence of toxigenic *C. difficile*. The prevalence of a positive *C. DIFF QUIK CHEK*<sup>®</sup> test at an independent study site was 18.2% (N=578).

## PERFORMANCE CHARACTERISTICS

### Clinical evaluation

The *C. DIFF QUIK CHEK*<sup>®</sup> test was compared to the bacterial culture test at two clinical laboratories and in-house at TECHLAB<sup>®</sup>, Inc. Specimens included in the evaluation were submitted to the clinical laboratories for routine testing. The presumptive bacterial culture test was performed according to the in-house procedures. The results are shown in Table 1.

**Table 1. Summary of clinical performance comparing *C. DIFF QUIK CHEK*<sup>®</sup> test to presumptive bacterial culture**

n=979	Presumptive Bacterial Culture positive	Presumptive Bacterial Culture negative
<i>C. DIFF QUIK CHEK</i> <sup>®</sup> positive	206	56
<i>C. DIFF QUIK CHEK</i> <sup>®</sup> negative	16	701
		95% Confidence Limits
Sensitivity	92.8%	88.3% - 95.7%
Specificity	92.6%	90.4% - 94.3%
Predictive Positive Value	78.6%	73.1% - 83.3%
Predictive Negative Value	97.8%	96.3% - 98.7%
Correlation	92.6%	91.7% - 93.4%

Discrepant samples were resolved using commercial tests for *C. difficile* glutamate dehydrogenase in ELISA or membrane test formats, or by a research PCR assay for the detection of *C. difficile* GDH gene *gluD*. Twenty-nine of the 56 apparent false positive samples were positive by another GDH test, and were considered true positives. Twenty-seven remained false positive. Thirteen of the 16 apparent false negative samples were negative by another GDH test, and were considered true negatives. Three remained false negative.

### EFFECT OF FECAL SPECIMEN CONSISTENCY

The reaction of fecal specimens of varying consistencies in the *C. DIFF QUIK CHEK*<sup>®</sup> test is shown in table 2. A total of 305 fecal samples of known consistency was included in the analysis. The percentages of positive reactions using either resolved culture assay or the *C. DIFF QUIK CHEK*<sup>®</sup> test were similar in all three types of fecal

specimens (liquid, semi-solid, and solid). All of the specimens were submitted for *C. difficile* testing. The basis of the submission was the clinical history of the patient and not the consistency of the specimen. The results show the *C. DIFF QUIK CHEK*<sup>®</sup> test performed similarly to the culture assay when testing samples of different consistencies.

**Table 2. Reaction of fecal specimens of varying consistencies in the *C. DIFF QUIK CHEK*<sup>®</sup> test**

Number of Specimens (n=305)	Liquid Specimens (n = 117)	Semi-solid Specimens (n = 142)	Solid Specimens (n = 46)
Positive by resolved bacterial culture assay	19 (16.2%)	37 (26.1%)	18 (39.1%)
Positive by <i>C. DIFF QUIK CHEK</i> <sup>®</sup>	18 (15.4%)	36 (25.4%)	21 (45.7%)
Negative by resolved bacterial culture assay	98 (83.8%)	105 (73.9%)	28 (60.9%)
Negative by <i>C. DIFF QUIK CHEK</i> <sup>®</sup>	99 (84.6%)	106 (74.6%)	25 (54.3%)

#### TRANSPORT MEDIA

A total of 69 fecal specimens diluted in C&S Transport Media was tested in the *C. DIFF QUIK CHEK*<sup>®</sup> test and the results were compared to those obtained by routine testing. The test exhibited an agreement of 95.7% for the detection of *C. difficile* antigen in specimens prepared in Transport Media.

#### REPRODUCIBILITY

A total of eight fecal specimens, 6 positive and 2 negative were coded to prevent identification and were sent to each of three independent laboratories for analysis using the *C. DIFF QUIK CHEK*<sup>®</sup> test. The results from each laboratory were compared with in-house results. The positive specimens were confirmed to be positive and the negative specimens were confirmed to be negative at each site.

#### CROSS REACTIVITY

Fecal specimens inoculated with the following microorganisms to a final concentration of approximately 10<sup>8</sup> or higher organisms per mL did not react in the *C. DIFF QUIK CHEK*<sup>®</sup> test:

**Bacteria:** *Aeromonas hydrophila*, *Bacillus cereus*, *Bacillus subtilis*, *Bacteroides fragilis*, *Campylobacter coli*, *Campylobacter fetus*, *Campylobacter jejuni*, *Candida albicans*, *Clostridium bifermentans*, *Clostridium butyricum*, *Clostridium perfringens* Type A, *Clostridium septicum*, *Clostridium sordellii* (toxigenic and nontoxigenic), *Clostridium sporogenes*, *Enterococcus faecalis*, *Escherichia coli* EIEC, *Escherichia coli*, *Escherichia coli* 0157 H7, *Escherichia coli* ETEC, *Klebsiella pneumoniae*, *Peptostreptococcus anaerobius*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Salmonella typhimurium*, *Shigella dysenteriae*, *Shigella flexneri*, *Shigella sonnei*, *Staphylococcus aureus*, *Staphylococcus aureus* (Cowans), *Staphylococcus epidermidis*, *Vibrio parahaemolyticus*, *Yersinia enterocolitica*

The following viruses of  $10^{3.3}$  to  $10^{7.5}$  TCID units per 0.2mL did not react in the *C. DIFF QUIK CHEK*® test:

**Viruses:** Adenovirus types 1,2,3,5,40,41, Human coronavirus, Coxsackievirus B2,B3,B4,B5, Echovirus 9,11,18,22,33, Enterovirus type 68,69,70,71.

#### **INTERFERING SUBSTANCES**

The following substances had no effect on test results when present in feces in the concentrations indicated: mucin (3.5% w/v), human blood (40% v/v), barium sulfate (5% w/v), Imodium® (5% w/v), Kaopectate® (5 mg/mL), Pepto-Bismol® (5% w/v), steric/palmitic acid (fecal fats, 40% w/v), Metronidazole (0.25% w/v), Vancomycin (0.25% w/v).

#### **ANALYTICAL SENSITIVITY**

The analytical sensitivity of the *C. DIFF QUIK CHEK*® test is 0.8 ng/mL.

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C. DIFF QUIK CHEK® is a trademark of TECHLAB®, Inc.
© 2009 TECHLAB®, Inc. All rights reserved.