



ASSESSMENT OF A NEW PARASITOLOGY SCREENING DIAGNOSTIC ELISA FOR THE DETECTION OF ANTIGENS OF GIARDIA SPP., CRYPTOSPORIDIUM SPP. AND ENTAMOEBIA HISTOLYTICA IN FECAL SPECIMENS

Nathaniel C. Christy¹, Janice D. Hencke², William A. Petri, Jr.¹, Aleya D. Escueta³, Forida Nazib⁴, Heidrun v.Thien⁵, Rashidul Haque⁴, Tomoyoshi Nozaki³, Egbert Tannich⁵, Joel F. Herbein²




¹University of Virginia, Charlottesville, VA, United States, ²TechLab[®], Inc, Blacksburg, VA, United States, ³National Institutes of Infectious Disease, Tokyo, Japan, ⁴International Center for Diarrheal Disease Research, Bangladesh, Dhaka, Bangladesh, ⁵Bernard Nocht Institute for Tropical Medicine, Hamburg, Germany



Introduction:

Episodes of diarrhea are common worldwide and while death from diarrheal infections has decreased, levels of morbidity have not declined in comparison to historical levels and thus remain a significant health problem, especially in the developing world. The three most common causes of protozoan-associated diarrheal infections are *Giardia* spp., *Cryptosporidium* spp. and *Entamoeba histolytica*. *Giardia* is a bi-nucleated parasite with a bi-phasic life cycle that causes infection following the ingestion of cysts. The parasite multiplies in the host as a flagellated trophozoite. Encystation of infectious cysts perpetuates the cycle of infection after exiting the host. *Cryptosporidium* spp. are common throughout both the developed and developing world and cause persistent diarrhea in HIV-infected populations. Infection occurs after exposure to oocysts, and control of this organism can prove difficult due to its resistance to standard disinfection methods (e.g. chlorinated water sources). *E. histolytica* is a single-cell ameba that is the cause of amebiasis. As with *Giardia*, infection occurs after ingestion of a multi-nucleated cyst and the emergence of disease-causing trophozoites that multiply in the intestine. Potentially fatal clinical manifestations include dysentery, toxic megacolon, and perforation of the intestine while some invasive cases progress to the development of deadly liver abscesses. Disease incidence caused by *Cryptosporidium* spp and *E. histolytica*, and possibly *Giardia* spp, is associated with malnutrition. Repeated infections are common and can cause developmental delay in small children. While treatment regimens for these infections are available, the surveillance for and prompt diagnosis of diarrheal illness in a population is critical for both prevention and treatment of disease. There remains a need for rapid and cost-effective diagnostic screening methods.

Background:



- Classic Fecal-Oral cycle of infection
- Giardiasis, Amebiasis, and Cryptosporidiosis are endemic in many third world countries and are common causes of enteric parasitic infection in the developed world.
- Repeated infections can cause malnutrition and developmental issues in children. (Mondal et al, 2009)(Petri et al, 2008)
- Inexpensive and easy-to-use diagnostic tools are needed to effectively screen for these parasites before appropriate treatment regimens can be applied.



Study Focus and Goals:

- o Evaluate the performance of a newly developed ELISA-type screen (The *Tri-Combo* parasite screen, TechLab[®], Inc.) for the presence of *Giardia* spp., *Cryptosporidium* spp., and/or *E. histolytica* antigens in human stool samples under clinical laboratory conditions.
- o The *Tri-Combo* parasite screen was compared to established ELISA-type tests on the market for the individual detection of *Giardia* spp., *Cryptosporidium* spp., or *E. histolytica* antigen in human stool samples

Methods:

- o A total of **618** clinical samples at three international sites:

International Center for Diarrheal Disease Research, Dhaka, Bangladesh:
297

National Institutes for Infectious Disease, Tokyo, Japan:
87

Bernard Nocht Institute for Tropical Medicine, Hamburg, Germany:
234

- o All samples were run on the *Tri-Combo* test and the *Giardia II*, *Cryptosporidium II*, and *E. histolytica II* ELISA tests (TechLab[®]).
- o Discrepant samples were re-run on the *Tri-Combo* and individual tests.

Results(1):

Comparison of the *TRI-COMBO* parasite screen by study site:

Study Site:	No. of Specimens with result			
	True Positive	False Positive	True Negative	False Negative
	NIID: Tokyo, Japan	13	0	73
ICDDR,B: Dhaka, Bangladesh	123	10	164	0
BNI: Hamburg, Germany	47	3	181	3
Combined Panel: All Sites	183	13	418	4

Results(2):

Comparison of the *TRI-COMBO* parasite screen by study site:

Study Site:	% of Specimens with result			
	Sensitivity(%)	Specificity(%)	Positive Predictive	Negative Predictive
			Value(%)	Value(%)
NIID: Tokyo, Japan	92.9	100	100	98.7
ICDDR,B: Dhaka, Bangladesh	100	94.25	92.48	100
BNI: Hamburg, Germany	94	98.4	94	98.4
Combined Panel: All Sites	97.9	97.0	93.4	99.1

Results(3):

Comparison of the *TRI-COMBO* parasite screen by organism:

Reference ELISA Test:	No. of Specimens with result			
	True Positive	False Positive	True Negative	False Negative
	<i>Giardia</i>	96	13 ¹	418
<i>Cryptosporidium</i>	47	13 ¹	418	1
<i>E. histolytica</i>	42	13 ¹	418	3

¹False positive results from the entire panel of samples (a total of 13) were recorded as FP for all 3 parasites.

Conclusions:

- o The *Tri-Combo* parasite screen effectively detected antigens from *Giardia* spp., *Cryptosporidium* spp., and *E. histolytica* parasites in clinical stool samples at three international sites as compared to three diagnostic tests specific for each individual organism.

References and Acknowledgments:

- Pierce, K. K., and B. D. Kirkpatrick. 2009. Update on human infections caused by intestinal protozoa. *Curr. Opin. Gastroenterol* 25:12-17.
 - Mondal, D., R. Haque, R. B. Sack, B. D. Kirkpatrick, and W. A. Petri Jr. 2009. Attribution of malnutrition to cause-specific diarrheal illness: evidence from a prospective study of preschool children in Mirpur, Dhaka, Bangladesh. *Am. J. Trop. Med. Hyg* 80:824-826.
 - Petri, W. A., Jr, M. Miller, H. J. Binder, M. M. Levine, R. Dillingham, and R. L. Guerrant. 2008. Enteric infections, diarrhea, and their impact on function and development. *J. Clin. Invest* 118:1277-1290.
 - Tzipori, S., and H. Ward. 2002. Cryptosporidiosis: biology, pathogenesis and disease. *Microbes and Infection* 4:1047-1056.
 - Carranza, P. G., and H. D. Lujan. 2010. New insights regarding the biology of *Giardia lamblia*. *Microbes and Infection* 12:71-80.
 - WHOPAN/UNESCO report. A consultation with experts on amoebiasis. Mexico City, Mexico 28-29 January, 1997. *Epidemiol Bull* 18:13-14.
- This work was supported by NIH grant AA03596 to W.A.P. and NIH (NIAD) grant R42 AI77171 to J.H. We thank TechLab[®], Inc. for the contribution of the *TRI-COMBO*, *GIARDIA II*, *CRYPTOSPORIDIUM II*, and *E. HISTOLYTICA II* stool antigen detection tests. WP receives royalties from a licensing agreement for amebiasis diagnostics with TechLab, Inc. Those royalties are donated in their entirety to the American Society for Tropical Medicine and Hygiene without benefit to Dr. Petri.