TRI-COMBO PARASITE SCREEN

A Monoclonal ELISA for Detecting *Giardia*, *Cryptosporidium* and
E. histolytica Antigens in Fecal Specimens

Catalog No. T30408 (96 Tests)

IVD In Vitro Diagnostic Medical Device

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TRI-COMBO PARASITE SCREEN

INTENDED USE

The *TRI-COMBO PARASITE SCREEN* test is an enzyme immunoassay for the simultaneous qualitative detection of *Giardia* spp., *Cryptosporidium* spp. and/or *E. histolytica* antigen in human fecal specimens. The test is indicated for use as a screen for fecal specimens from patients with gastrointestinal illness including dysentery to aid in the diagnosis of enteric infection resulting from giardiasis, cryptosporidiosis and amebiasis. The *TRI-COMBO PARASITE SCREEN* test serves as a screen to rule out negative specimens and reduce the number of specimens that require additional follow-up testing.

EXPLANATION

Giardia is a binucleated, flagellated protozoan parasite which exists in two forms: a non-infectious, pear-shaped trophozoite (9 to 21 μm) inhabiting the small intestine and the highly-infectious cyst form, which is elliptical in shape and ranges in size from 8 to 12 μm (reviewed in 1). Survival outside its host varies greatly between the two forms: the trophozoite, which is extremely labile, lasts only a matter of hours outside the body, while the cyst form may survive for several days in an external environment (1). *Giardia* often is responsible for infections due to water contamination, and travelers have been found to contract giardiasis from endemic areas. Transmission also occurs by direct contact, especially with asymptomatic carriers and by food contamination. Giardiasis is known to be a common sexually-transmitted disease (2). Animal fecal contamination, especially of water, is another route of transmission in humans (1).

Clinical manifestations of giardiasis range from asymptomatic carriage with the passing of cysts to chronic debilitating diarrhea, weight loss and malabsorption. High-risk categories include young children, immunocompromised patients and those without previous exposure (1).

Cryptosporidium is a protozoan parasite of vertebrates originally thought to cause diarrhea only in animals (3). In 1976, the first human infection was reported (4). Since then, *Cryptosporidium* spp. has been found to be associated with diarrheal illness in most parts of the world and is a frequent cause of traveler's diarrhea. The disease is transmitted by the thick-walled oocyst form, 2-6 μm in diameter, which is remarkably resistant to common disinfectants and routine chlorination of drinking water. Person-to-person transmission, especially among children, is common. *Cryptosporidium* has little or no host specificity, and animals such as rodents, cattle and domestic pets can serve as a reservoir for zoonotic transmission to humans, transmission occurring either by direct contact or by contamination of water supplies with fecal matter. Cryptosporidiosis is a serious opportunistic infection in patients with acquired immunodeficiency syndrome (AIDS) and is potentially a sexually transmitted disease (5,6).

Clinical manifestations of cryptosporidiosis include cholera-like diarrhea, abdominal pain, nausea, vomiting and weight loss. In normally healthy persons, the infection is usually self-limiting and of short term. In AIDS and other immunocompromised patients, cryptosporidiosis can result in prolonged and life-threatening illness due to excessive fluid loss (7). In these patients the infection may also spread to the respiratory and biliary tracts (5).

The life cycle of *Entamoeba histolytica* is similar to that observed with other amoebae. The organism exists either as a trophozoite, ranging from 10-60 μm in diameter, or as a cyst (12-15 μm) (8). Humans serve as the primary reservoir, with the organism being spread through ingestion of contaminated food and water or through sexual contact. The cyst is very stable. Once it enters the intestine, it begins to divide into trophozoites that bind to the intestinal mucosa via a galactose or N-acetyl-D-galactosamine-binding lectin referred to as the galactose adhesin. Once the trophozoites have attached, they release tissue-damaging enzymes and proteins that lyse the mucosal cells. The galactose adhesins of *E. histolytica* and *E. dispar* cross-react serologically but contain distinct epitopes.

The adhesin is antigenically conserved, but monoclonal antibodies are used to distinguish between the adhesin from pathogenic *E. histolytica* and non-pathogenic *E. dispar* (reviewed in 9).

E. histolytica and *E. dispar* are intestinal parasites that infect approximately half a billion people worldwide annually (10). Only malaria and schistosomiasis are believed to be more prevalent parasitic causes of morbidity and mortality. Of the huge number of persons infected, most are infected with *E. dispar*, which has not been associated with disease. Infection with *E. dispar* rather than *E. histolytica* is believed to explain, at least in part, the low prevalence of disease considering the high rate of infection. Even so, it is estimated that approximately 10% of the half billion people infected each year are infected with *E. histolytica*. These individuals can become symptomatic and develop colitis and liver abscesses, resulting in a mortality rate estimated between 40,000 and 120,000 persons annually. It is necessary to distinguish between the two species because *E. histolytica* is pathogenic, whereas *E. dispar* is not associated with colitis or liver abscess. Inaccurate diagnosis may result in unwarranted and unnecessary treatment with harsh drugs with side effects.

Patients with pathogenic *E. histolytica* may show a wide range of conditions. Some are completely asymptomatic. These persons may shed millions of cysts daily, representing a potential reservoir for dissemination. Some patients may show mild diarrhea that develops into bloody diarrhea with abdominal cramps, eventually resulting in fulminant colitis. Because of the tissue damage that can occur, perforation of the intestine may result and the amoebae may disseminate to other parts of the body. Approximately 10% of persons with invasive amebiasis develop liver abscesses (10).

The most common method used to detect giardiasis, cryptosporidiosis and amebiasis has been diagnosis by O&P microscopy. However, accuracy of O&P results depends upon the skill of the technician and also relies on the presence of intact cysts in the feces, which may not be present in all samples. In addition, it is only rarely possible to distinguish between the pathogenic and non-pathogenic species of *Entamoeba* using microscopy. The success rate of fecal examination by microscopy varies between 50 and 70%, and multiple specimens are usually necessary to establish a diagnosis. When infection is present but parasites are not detected by microscopy, sampling and testing of duodenal fluid may detect trophozoites, but this method is invasive and expensive. Detection of the organism and antigens by ELISA provides an alternative method of diagnosis and is sensitive and specific (9). The ELISA procedure is straightforward to perform and exhibits increased sensitivity compared to microscopic examination. Large numbers of specimens can be tested rapidly and objectively, and the procedure is less labor-intensive than most microscopy methods.

PRINCIPLE OF THE TEST

The **TRI-COMBO PARASITE SCREEN** test uses monoclonal and polyclonal antibodies to cell-surface antigens of *Giardia*, *Cryptosporidium* and *E. histolytica*. The microassay plate in the kit contains immobilized monoclonal antibodies against the antigens, and the *Conjugate* consists of polyclonal antibodies against the antigens. In the assay, an aliquot of a diluted fecal specimen is transferred to a microassay well. The immobilized monoclonal antibodies bind the *Giardia*, *Cryptosporidium* and/or *E. histolytica* antigens if they are present. Upon addition, *Conjugate* then binds to the antigen/antibody complex. Any unbound materials are removed during the washing steps. Following the addition of *Substrate*, a color is detected due to the enzyme-antibody-antigen complexes that formed in the presence of antigens and conjugate.

MATERIALS PROVIDED

- CONJENZ** **Conjugate (7 mL)** – Antibodies against antigens of *Giardia* spp., *Cryptosporidium* spp., and *E. histolytica* in a buffered protein solution containing 0.02% Thimerosal
- DIL SPE** **Diluent (50 mL)** – Buffered protein solution containing 0.02% thimerosal. The *Diluent* is also to be used as the negative control solution (see TEST PROCEDURE)
- H₂SO₄0.6N** **Stop Solution (7 mL)** – 0.6 N sulfuric acid. CAUTION: Avoid contact with skin; flush with water immediately if contact occurs
- CONTROL+** **Giardia Positive Control (3.5 mL)** – *Giardia* antigen in a buffered protein solution with 0.02% thimerosal.

CONTROL+	Cryptosporidium Positive Control (3.5 mL) – <i>Cryptosporidium</i> antigen in a buffered protein solution with 0.02% thimerosal
CONTROL+	E. histolytica Positive Control (3.5 mL) – <i>E. histolytica</i> antigen in a buffered protein solution with 0.02% thimerosal
SUBS/REAG	Substrate (14 mL) – solution containing tetramethylbenzidine and peroxide
WASHBUF/20X	Wash Buffer Concentrate (50 mL) – 20X concentrate containing phosphate buffered saline, detergent, and 0.2% thimerosal
MA/PLT	Microassay Plate – 12 strips, each consisting of 8 wells coated with antibodies to <i>Giardia spp.</i> , <i>Cryptosporidium spp.</i> , and <i>E. histolytica</i> antigen (stored with desiccant)
IVD	<i>In Vitro</i> Diagnostic Medical Device

ACCESSORIES

100 Disposable plastic transfer pipettes 2 Plastic adhesive sheets 1 Wash Solution Label

MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

Squirt bottle for wash reagent	Vortex mixer
950 mL distilled water for diluting wash reagent	Discard container
Absorbent paper	Applicator sticks
ELISA reader capable of reading at 450 nm or 450/620 nm	
Small tubes for dilution of fecal specimens (e.g., microcentrifuge tubes)	

SHELF LIFE AND STORAGE

The expiration date of this kit is given on the box label. Expiration dates for each component are listed on the individual labels. The kit should be stored between 2° and 8°C and should be returned to the refrigerator as soon as possible after use.

PRECAUTIONS

- Each component in the kit should be inspected for any signs of leakage. Upon arrival, the kit should be inspected to ensure that components are not frozen or warm to the touch due to improper shipping conditions.
- Reagents from different kits should not be mixed. Do not use a kit past the assigned expiration date.
- Bring all components to ROOM TEMPERATURE BEFORE USE!
- Do not freeze the reagents. The kit should be stored between 2°C and 8°C.
- Caps and tips are color-coded; do NOT mix or interchange!
- When handling assay wells, avoid scratching the bottom of the wells because this may result in elevated absorbance readings.
- Unused microwells must be placed back inside of the resealable pouch with the desiccant to protect them from moisture. Check the desiccant pack before using the microwells. The color indicator on the desiccant pack should be blue. If the color turns pink, the quality of the microwells may be compromised. Please do not use microwells stored with pink desiccant.
- Hold reagent bottles vertically when dispensing to ensure proper drop size and correct volume.
- Handle specimens and used microassay wells as if capable of transmitting infectious agents. Wear gloves when performing the test.
- Reagents contain 0.02% thimerosal as a preservative and should be handled with normal laboratory caution.
- Microbial contamination of reagents may decrease the accuracy of the assay. Avoid microbial contamination of reagents by using sterile disposable pipettes if removing aliquots from reagent bottles.
- All reagents, with the exception of the *Wash Buffer Concentrate*, are supplied in ready-to-use bottles. Reagents can be dispensed directly from the dropper bottles or decanted for use with multichannel pipettes. If excess reagent has been decanted, the excess should be

discarded. Do not pour back into the bottle. The *Substrate* should be stored in and used from the light-protected bottle in which it is supplied. If an aliquot is removed from the original bottle for any reason, do not return unused *Substrate* to the original bottle.

- Perform the washing procedure as directed to avoid high background reactions.
- The *Stop Solution* contains 0.6 N sulfuric acid. Flush with water immediately if contact with skin or eyes occurs.
- The *Substrate* is light sensitive and should be protected from direct sunlight or UV sources.
- The test has been optimized for sensitivity and specificity. Do not deviate from the specified procedure. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test.
- Fecal specimens may contain potentially infectious agents and should be handled at "Biosafety Level 2" as recommended in the CDC/NIH Manual, "Biosafety in Microbiological and Biomedical Laboratories."

PRELIMINARY PREPARATIONS

- All reagents must be at room temperature prior to use in the assay.**
- Prepare 1X Wash Solution.** The *Wash Buffer Concentrate* is supplied as a 20X concentrate (a precipitate may be noticed). It should be diluted to a total volume of 1 liter by adding 50 mL of the concentrate to 950 mL of distilled water. The 1X *Wash Solution* can be stored between 2° and 8°C.
- Assay Strip Preparation.** Each strip contains 8 wells coated with antibodies to *Giardia spp.*, *Cryptosporidium spp.*, and *E. histolytica* antigen. Each specimen or control will use one of these coated wells. Determine the number of wells to be used. Avoid contact with the base of the wells. Assay wells not used must be returned to the foil pouch and carefully resealed with desiccant.

SPECIMEN HANDLING AND PREPARATION


Standard collection and handling procedures used in-house for fecal specimens are appropriate. No modification of collection methods used for standard microscopic O&P examinations is needed.

Acceptable Sample Types
Fresh Fecal Specimen
Frozen Fecal Specimen (frozen undiluted samples)

Do Not Use
Fecal specimens in Formalin-based fixative (e.g., sodium acetate formalin, 10% formalin)
Fecal specimens in alcohol-based fixative (e.g., polyvinyl alcohol)
Concentrated Fecal Specimens

- Specimens should be kept between 2° and 8° C and tested within 24 hours of collection. Specimens that cannot be tested within this time should be stored at ≤ -10°C until tested. Repeated freeze/thaw cycles should be avoided. If using frozen specimens, thaw at room temperature.
- Set up and label one test tube for each sample as necessary.
Add 400 µL Diluent to each tube.
- Obtain one disposable plastic transfer pipette (supplied with the kit) for each sample.

Transfer Pipette

4. **Mix all specimens and cultures thoroughly regardless of consistency - it is essential that the samples be evenly suspended before sampling.** Using the disposable plastic transfer pipette, add 100 μ L (second graduation mark) of fecal specimen to the tube containing *Diluent* and mix well. If the specimen cannot be pipetted, use an applicator stick to transfer approximately 0.1 g of feces. This is about the size of a small pea (about 4mm in diameter).
5. Close each tube of diluted sample and mix thoroughly. Proper mixing can be achieved by vortexing or inverting the tube several times.
6. If using semi-automated or automated washing equipment, once diluted, specimens must be centrifuged (5000 x g for 10 minutes) to remove any particulate matter from the supernatant before transfer to assay wells.

TEST PROCEDURE

1. Bring all reagents and the required number of test strips to room temperature before use.
2. **Prepare Controls**
 - a. Positive Controls – shake each *Positive Control* bottle for several seconds, then add one drop to each positive control well.
 - b. Negative Controls – add 100 μ L of *Diluent* to an additional well to serve as the negative control.
3. **Using a new transfer pipette for each specimen, transfer 100 μ L of diluted specimen to the assay well of the *Microassay Plate*.**
4. Cut the adhesive plastic sheet to the size necessary to cover the wells. **Cover the wells and incubate them at for 1 hour at room temperature.**
5. Shake out contents of assay wells into a discard pan.
6. **Wash each well using the 1X Wash Solution in a squirt bottle with a fine-tipped nozzle,** directing the *Wash Solution* to the bottom of the well with force. Fill the wells, and then shake the **Wash Solution** out of the well into a discard pan. Slap the inverted plate on a dry paper towel.
Note: If using semi-automated or automated washing equipment, add 350 μ L of 1X Wash Solution to each well. Wash for a total of 5 times
7. Repeat step 6 four additional times using a dry paper towel each time. If any particulate matter is seen in the wells, continue washing until all the particulate matter is removed.
8. After washing, completely remove any residual liquid in the wells by striking the plate onto a dry paper towel until no liquid comes out. Dispose of paper towels and specimen containers properly.
9. **Add 1 drop (50 μ L) of *Conjugate* (red cap) to each well.**
10. Seal with a plastic adhesive sheet and gently tap to mix. **Incubate the wells for 30 minutes at room temperature.**
11. After incubating, **repeat the washing procedure described in steps 5-8.**
12. **Add 2 drops (100 μ L) of *Substrate* (blue cap) to each well.** Gently tap the wells to mix the substrate. Incubate the wells at room temperature for 10 minutes. Gently tap the wells at 5 minutes.
13. **Add 1 drop (50 μ L) of *Stop Solution* (yellow cap) to each well.** Gently tap the wells and wait 2 minutes before reading. The addition of the *Stop Solution* converts the blue color to a yellow color which may be quantitated by measuring the optical density at 450 nm on a microplate ELISA reader. The instrument should be blanked against air. If a dual wavelength reader is used, blank against air at 620 and read at 450 nm. Wipe the underside of each well before measuring the optical density. If an ELISA reader is unavailable, the test may be read visually in good light against a white background. Read within ten minutes after adding *Stop Solution*.

QUALITY CONTROL

1. Positive and negative controls must be run with each series of test specimens. The positive controls demonstrate that the assay is functioning properly for the detection of each parasite antigen in fecal specimens. The negative control demonstrates that the assay is not reacting nonspecifically.
2. Positive and negative controls must fall within their respective ranges or the test results are not valid.
 - a) **Positive Control must be a visible yellow color.** If read on a spectrophotometer, the OD at 450 nm or using dual wavelength at 450/620 nm must be ≥ 0.500 . Any well that gives a positive reading without visible color should be repositioned, wiped on the underside of the well, and read again.
 - b) **Negative Control must be visually clear.** If read on a spectrophotometer, the OD at 450 nm must be < 0.200 . If read at 450/620 nm the absorbance must be < 0.160 . If not, the test is invalid and should be repeated, paying attention to the wash procedure.
3. Visual readings must be taken in good light against a white background.
4. Test results are not valid unless the performance characteristics of the positive and negative controls are met. If these results are not observed, call Technical Services.

INTERPRETATION OF RESULTS

	Spectrophotometric Reading	
	Single Wavelength at 450 nm	Dual Wavelength at 450/620 nm
Negative	OD < 0.200	OD < 0.160
Positive	OD ≥ 0.200	OD ≥ 0.160

Visual Interpretation

The negative control well should be colorless or have only a faint yellow color. The positive control well should give an obvious yellow color. If these results are not observed, call Technical Services. A test sample is considered positive if it has an obvious yellow color when compared to the negative control well. It may be less yellow or more yellow than the color observed in the positive control well. A test sample is considered negative if the reaction is colorless or less yellow than the negative control well.

A positive result indicates that antigen from *Giardia* spp., *Cryptosporidium* spp. and/or *E. histolytica* is present in the specimen. A negative result indicates that *Giardia* spp., *Cryptosporidium* spp. and/or *E. histolytica* antigens are absent or the level is below the detection limit of the test.

LIMITATIONS OF THE PROCEDURE

1. The **TRI-COMBO PARASITE SCREEN** test detects the presence of *Giardia*, *Cryptosporidium* and/or *E. histolytica* antigens in fecal specimens. Test results should be taken under consideration by the physician in light of the clinical history and physical examination of the patient.
2. The **TRI-COMBO PARASITE SCREEN** test is intended for the qualitative detection of *Giardia* spp., *Cryptosporidium* spp., and/or *E. histolytica* antigens in unpreserved fecal specimens. The intensity of the color should not be interpreted quantitatively.
3. Concentrated fecal specimens should not be used in the test and may not give accurate results.

EXPECTED VALUES

Normal healthy individuals should not be infected with *Giardia*, *Cryptosporidium* or *E. histolytica* and should test negative in the *TRI-COMBO PARASITE SCREEN* test. A positive test result in the *TRI-COMBO PARASITE SCREEN* test indicates that the person is shedding detectable amounts of *Giardia*, *Cryptosporidium* and/or *E. histolytica* antigen. The incidence of *Giardia*, *Cryptosporidium* and *E. histolytica* infection varies significantly between populations and geographic regions. Children in daycare settings have exhibited higher rates of infection with *Giardia* than the general population (11). In addition, homosexual men have shown higher rates of infection (9,12). In general, laboratory-confirmed incidence of cryptosporidiosis in developed countries ranges from 1 to 2% overall with a higher incidence in children. It is estimated that *Entamoeba histolytica* infects about 50 million people around the world (9). Roughly 90% of these persons remain asymptomatic, whereas about 10% develop clinical symptoms ranging from gastrointestinal disease to liver abscesses. High risk groups include persons who have traveled abroad, immigrants, immunocompromised persons, migrant workers, and active male homosexuals (9,14). Nonpathogenic strains (*E. dispar*) are predominant among male homosexuals (15). The disease often is transmitted by asymptomatic carriers of *E. histolytica*.

PERFORMANCE CHARACTERISTICS

The performance of the *TRI-COMBO PARASITE SCREEN* test was compared to a panel of single parasite detection predicate devices (for the detection of *Giardia* spp., *Cryptosporidium* spp., and *E. histolytica*), and included 275 fresh and 645 frozen specimens. The following table shows a summary of the clinical performance of the *TRI-COMBO PARASITE SCREEN* test to the predicate panel. The results show that the test exhibited a positive percent agreement of 97.8%, negative percent agreement of 97.5% and overall percent agreement of 97.6%.

n = 920	Predicate Panel Positive	Predicate Panel Negative
TRI-COMBO Positive	307	15
TRI-COMBO Negative	7	591

	95% Confidence Limits	
Positive Agreement	97.8%	95.3 – 99.0%
Negative Agreement	97.5%	95.9 – 98.6%
Overall Agreement	97.6%	97.3 – 97.8%

REPRODUCIBILITY

The reproducibility of the *TRI-COMBO PARASITE SCREEN* test was determined using 20 human fecal specimens coded to prevent their identification during testing. Testing was performed at 2 independent laboratories and on-site at TECHLAB®, Inc. The samples were tested twice a day over a 5-day period by multiple technicians at each site using 2 different kit lots. Positive and negative controls were run with each panel of the masked samples. The results from each laboratory were submitted to TECHLAB®, Inc. and compared with in-house results. The results were consistent among the different locations and exhibited a correlation of 100%. The samples produced the expected results 100% of the time.

CROSS REACTIVITY

The *TRI-COMBO PARASITE SCREEN* test was evaluated for cross-reactivity with the bacterial and viral strains listed below. None of the strains were shown to interfere with the performance of the *TRI-COMBO PARASITE SCREEN* test.

<i>Aeromonas hydrophila</i>	<i>Escherichia coli</i> ETEC
<i>Bacillus cereus</i>	<i>Escherichia coli</i> 0157:H7
<i>Bacillus subtilis</i>	<i>Klebsiella pneumoniae</i>
<i>Bacteroides fragilis</i>	<i>Salmonella enterica typhimurium</i>
<i>Campylobacter coli</i>	<i>Shigella dysenteriae</i>
<i>Campylobacter fetus</i>	<i>Shigella flexneri</i>
<i>Campylobacter jejuni</i>	<i>Shigella sonnei</i>
<i>Candida albicans</i>	<i>Staphylococcus aureus</i>
<i>Clostridium bifermentans</i>	<i>Staphylococcus aureus</i> Cowans
<i>Clostridium difficile</i>	<i>Staphylococcus epidermidis</i>
<i>Enterococcus faecalis</i>	<i>Vibrio cholerae</i>
<i>Escherichia coli</i>	<i>Vibrio parahaemolyticus</i>
<i>Escherichia coli</i> EIEC	<i>Yersinia enterocolitica</i>

Adenovirus, Type 1, 2, 3, 5, 40, and 41	Echovirus, Type 9, 11, 18, 22, and 33
Coronavirus, Human	Enterovirus, Type 68, 70, and 71
Coxsackievirus, Type, B2, B3, B4, and B5	Rotavirus, Human

INTERFERING SUBSTANCES (U.S. FORMULATIONS)

The following substances had no effect on positive or negative test results analyzed at the concentrations indicated: Hog gastric mucin (3.5% w/v), Human blood (40% v/v), Barium sulfate (5% w/v), Imodium® (5% v/v), Kaopectate® (5% v/v), Pepto-Bismol® (5% v/v), Mylanta (5% v/v), Steric/Palmitic Acid (40% w/v), Metronidazole (0.25% w/v), Vancomycin (0.25% w/v).

PRECISION – INTRA-ASSAY

For the determination of intra-assay performance, 6 positive and 5 negative fecal specimens were analyzed. Each specimen was assayed in replicates of eight. All positives remained positive and all negatives remained negative.

PRECISION – INTER-ASSAY

The inter-assay precision of the *TRI-COMBO PARASITE SCREEN* test was determined by using 20 fecal specimens (8 negative, 4 *Giardia* positive, 4 *Cryptosporidium* positive, and 4 *E. histolytica* positive). The samples were tested, twice a day over a 5-day period using 2 different kit lots. Positive and negative controls were run on each day. All positive samples remained positive and all negative samples remained negative.

ANALYTICAL SENSITIVITY

The analytical sensitivity of the test was determined by using purified cyst or trophozoite dilution curves in a sample matrix. The concentration of *Giardia* cysts, *Cryptosporidium* oocysts and *E. histolytica* trophozoites (Pathogenic Zymodeme Amoeba [PZ]) where specimens were positive by the *TRI-COMBO PARASITE SCREEN* test 95% of the time was utilized to describe the assay limit-of-detection (LOD). Test results determined the LOD for the assay to be 8,000 cysts/mL for *Giardia* (equivalent to 160 cysts detected per test), 40,000 oocysts/mL for *Cryptosporidium* (equivalent to 800 oocysts detected per test), and 3,300 PZ/mL for *E. histolytica* (equivalent to 66 PZs detected per test). Because the *TRI-COMBO PARASITE SCREEN* test detects soluble antigen in fecal specimens in addition to cysts, oocysts, and trophozoites, this LOD study represents an estimate of analytical sensitivity based on purified cysts, oocysts, and trophozoites.

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The Alere Logo and Alere are trademarks of the Alere group of companies.
The TECHLAB Logo and TECHLAB are trademarks of TECHLAB®, Inc., under license.
© 2014 TECHLAB®, Inc. All rights reserved.
Alere North America, LLC.