TECHLAB® *E. HISTOLYTICA* QUIK CHEK™

TECHLAB® E. HISTOLYTICA QUIK CHEK™

ENGLISH p. 3

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TECHLAB® *E. HISTOLYTICA* QUIK CHEK™

INTENDED USE

The TECHLAB® *E. HISTOLYTICA* QUIK CHEK™ test is a rapid membrane enzyme immunoassay for the qualitative detection of adhesin from *Entamoeba histolytica* in a single use cassette. It is intended for use with human fecal specimens from patients with diarrhea or dysentery as an aid in the diagnosis of *E. histolytica* gastrointestinal infection. Test results should be considered in conjunction with patient history.

EXPLANATION

Entamoeba histolytica and *Entamoeba dispar* are intestinal parasites that infect approximately half a billion people worldwide each year (1). It is necessary to distinguish between the two species because *E. histolytica* is pathogenic, causing intestinal amebiasis (e.g. diarrhea, dysentery, colitis) and extra-intestinal amebiasis (e.g. liver abscesses). *E. dispar* is not associated with symptomatic disease and inaccurate diagnosis may result in unnecessary treatment. The most common method used to diagnose amebiasis has been wet mount microscopy, which suffers from poor sensitivity and specificity. Trophozoites and cysts are not easily identified in a single fecal specimen and it is difficult to visually distinguish between *E. dispar* and *E. histolytica* when they are observed. Detection of *Entamoeba* spp. by immunoassay provides an alternative method of diagnosis with greater sensitivity (2). Immunoassays specific for *E. histolytica*, such as the *E. HISTOLYTICA* QUIK CHEK test, provide the added benefit of only identifying *E. histolytica* infections. Large numbers of specimens can be tested rapidly and objectively, and the procedure is less labor-intensive than most other methods of diagnosis.

PRINCIPLE OF THE TEST

The *E. HISTOLYTICA* QUIK CHEK test utilizes antibodies specific for *E. histolytica*. The *Membrane Device* contains a *Reaction Window* with two vertical lines of immobilized antibodies. The test line ("T") contains monoclonal antibodies specific for *E. histolytica* adhesin. The control line ("C") contains antibodies to horseradish peroxidase (HRP). The *Conjugate* consists of antibodies to *E. histolytica* coupled to horseradish peroxidase. To perform the test, the sample is added to a tube containing a mixture of *Diluent* and *Conjugate*. The diluted sample-conjugate mixture is added to the *Sample Well* and the device is allowed to incubate at room temperature for 15 minutes. During the incubation, any *E. histolytica* adhesin in the sample binds to the antibody-peroxidase conjugate. The antigen-antibody-peroxidase complexes migrate through a filter pad to a membrane where they are captured by the immobilized anti-adhesin antibodies in the line. The *Reaction Window* is subsequently washed with *Wash Buffer*, followed by the addition of *Substrate*. After a 10 minute incubation period, the *Reaction Window* is examined visually for the appearance of vertical blue lines on the "C" and "T" sides of the *Reaction Window*. A blue line on the "T" side of the *Reaction Window* indicates a positive result. A positive "C" reaction, indicated by a vertical blue line on the "C" side of the *Reaction Window*, monitors/ confirms that the sample and reagents were added correctly, the reagents were active at the time of performing the assay, and that the sample migrated properly through the *Membrane Device*. It also confirms the reactivity of the other reagents associated with the assay.

MATERIALS PROVIDED

MEM DEV

Membrane Devices – Each pouch contains 1 device

CONJ ENZ

Conjugate (2 mL) – Antibody specific for *E. histolytica* coupled to horseradish peroxidase in a buffered protein solution

DIL SPE

Diluent (16 mL) – Buffered protein solution with black graduated dropper assembly

CONTROL +

Positive Control (1 mL) – *E. histolytica* antigen in a buffered protein solution

SUBS REAG

Substrate (3.5 mL) – Solution containing tetramethylbenzidine

WASH REAG

Wash Buffer (12 mL) – Buffered solution with white graduated dropper assembly

Disposable plastic pipettes (50) – Graduated at 25 µL, 100 µL, 200 µL, 300 µL, 400 µL, and 500 µL

IVD

In Vitro Diagnostic Medical Device

MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

Small test tubes (e.g., plastic microcentrifuge tubes)

Wooden Applicator Sticks

Disposable gloves

Pipettor and tips

Timer

SHELF LIFE AND STORAGE

The expiration date of the kit is printed on the kit box label. Store the kit between 2°C and 8°C. Return the kit to the refrigerator as soon as possible after use.

PRECAUTIONS

- Upon arrival, inspect the kit to ensure that components are not frozen or warm to the touch due to improper shipping conditions, and that there are no signs of leakage.
- Bring all components to room temperature before use to ensure proper kit reactivity.
- The *Substrate* reagent should be colorless. If the *Substrate* reagent changes to a dark blue/violet color, discard and call Technical Services for a replacement.
- Reagents from different kits should not be mixed or interchanged. Do not use a kit past the expiration date.
- Caps, tips and dropper assemblies are color-coded; do NOT mix or interchange!
- Use fecal specimens within 72 hours of collection to obtain optimal results. Specimens that are frozen may lose reactivity due to freezing and thawing.
- The test has been optimized for sensitivity and specificity. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test. Do not deviate from the specified procedure.
- Fecal specimens and used membrane devices may contain potentially infectious agents and should be handled at "Biosafety Level 2" as recommended in the CDC/NIH Manual "Biosafety in Microbiological and Biomedical Laboratories." Wear disposable gloves when performing the test.
- For *in vitro* diagnostic use.

COLLECTION, HANDLING AND STORAGE OF FECAL SPECIMENS

Acceptable Sample Types
Fresh Fecal Specimens
Frozen Fecal Specimens (frozen undiluted)

Do Not Use
Fecal specimens in Formalin-based fixative (e.g. sodium acetate formalin, 10% formalin)
Fecal specimens in alcohol-based fixative (e.g. polyvinyl alcohol)
Fecal specimens in transport media (e.g. Cary Blair)

- Use standard in-house collection and handling procedures for fecal specimens. Collect fecal specimens in clean, leak-proof containers.
- Store fresh fecal specimens between 2°C and 8°C up to 72 hours. Freeze specimens and store at ≤ -10°C if testing cannot be performed within 72 hours of collection. Test specimens that are less than 24 hours old, whenever possible.
- Repeated freeze/thaw cycles should be avoided. If using frozen specimens, thaw at room temperature.
- Do not store fecal specimens in the *Diluent*.

TEST PROCEDURE

- Be attentive to the total assay time when testing more than one fecal specimen.
- Bring all reagents and devices to room temperature before use. Remove the reagents from the foam insert to reduce the time needed to warm to room temperature
- Set up and label one small test tube for each specimen and optional external control.**
- Add 500µL of *Diluent* to each tube using the black graduated dropper assembly.**

Sample Type	Volume of <i>Diluent</i>
Fresh Fecal Specimens	500 µL
Frozen Fecal Specimens (frozen undiluted)	500 µL

- Add one drop of *Conjugate* (red capped bottle) to each tube.** Hold the dropper bottle vertically to ensure proper drop size. The *Diluent* and *Conjugate* should be added to all tubes prior to adding the specimens.

- Obtain one disposable plastic transfer pipette (supplied with the kit) for each sample.

Graduated Transfer Pipette:

- For Liquid/Semi-Solid Specimens** - Mix specimen thoroughly. Using a transfer pipet, add 25 µL of specimen to the *Diluent/Conjugate* mixture in the tube.

For Formed/Solid specimens - Mix specimen thoroughly using a wooden applicator stick and transfer a small portion (approximately 2 mm diameter, the equivalent of 25 µL) of the specimen into the *Diluent/Conjugate* mixture. Emulsify the specimen using the applicator stick.

NOTE: Transferring too little sample, or failure to mix and completely suspend the sample in the *Diluent/Conjugate* mixture, may result in a false-negative test result. The addition of too much sample may cause invalid results due to restricted flow.


- Optional External Controls:**

External Positive Control - add one drop of *Positive Control* (gray-capped bottle) to the appropriate test tube.

External Negative Control - add 25 µL *Diluent* to the appropriate test tube.

- For all test and control samples, close the tubes and mix thoroughly using a vortex mixer or by inverting the tube several times. Samples or controls diluted in the *Diluent/Conjugate* mixture may be incubated at room temperature up to 2 hours prior to addition to the *Membrane Device*.
- Open one room temperature *Membrane Device* pouch for each diluted specimen and external control (as necessary). Label each device appropriately and orient it on a flat surface so the "E HISTO QUIK CHEK" print is at the bottom of the device, and the small *Sample Well* is located in the top right corner of the device.


Membrane Device

- Make sure that each diluted sample is thoroughly mixed (See Step 9) before adding to the *Membrane Device*. **Using a new transfer pipette, transfer 500 µL (topmost graduation) from each tube into the *Sample Well* (smaller hole in the top right corner of the device) of a *Membrane Device*.** When adding the sample into the *Sample Well*, make sure that the tip of the transfer pipette is inside the *Sample Well* hole and angled towards the *Reaction Window*. Expel the diluted sample onto the wicking pad inside the *Membrane Device*.

12. **Incubate the device at room temperature for 15 minutes** – the sample will wick through the device and a wet area will spread across the *Reaction Window*. The 15-minute incubation step begins after the last diluted sample-conjugate mixture has been transferred to the last *Membrane Device*.
NOTE FOR SAMPLES THAT FAIL TO MIGRATE:
Occasionally, a diluted sample fails to migrate properly and the Reaction Window does not fully wet. If the Reaction Window does not appear to be completely wet within 5 minutes of adding the sample to the Sample Well, then add 100 µL (4 drops) of Diluent to the Sample Well and wait an additional 5 minutes (for a total of 20 minutes). Continue with the next step of the Test Procedure.
13. **After the incubation, add 300 µL of Wash Buffer to the central Reaction Window using the graduated white dropper assembly.** Allow the Wash Buffer to be absorbed completely.
14. **Add 2 drops of Substrate (white-capped bottle) to the central Reaction Window.**
15. Incubate 10 minutes at room temperature. Read visually and record results after the incubation.

INTERPRETATION OF RESULTS

NOTE: The External Controls are intended to monitor for gross systematic errors only. If the expected results are not observed with the positive and negative Controls, this could indicate that the controls themselves are not working correctly.

1. Interpretation of the test is most reliable when the device is read immediately at the end of the reaction, in a well-lit area, and from directly over the device at a normal working distance.
2. **Positive Result:** Two vertical blue lines are visible, the control line on the "C" (left) side of the *Reaction Window* and the test line on the "T" (right) side of the *Reaction Window*. The lines may appear faint to dark in intensity - any line on the "T" side is considered positive. Do not interpret membrane discoloration as a positive result. A positive result indicates the presence of *E. histolytica* antigen, and that there is a properly reactive positive control line.
3. **Negative Result:** A single vertical blue line is visible on the "C" (left) side of the *Reaction Window* and no test line is visible on the "T" side of the *Reaction Window*. A negative result indicates that *E. histolytica* antigen is either absent in the sample or is below the detection limit of the test, and that there is a properly reactive positive control line.
4. **Invalid Result:** A single line is visible on the "T" side of the *Reaction Window*, or no lines are visible in the *Reaction Window*. The test is invalid if a control line is not present at the completion of the test reaction.

QUALITY CONTROL

The validity of the test results using the *E. HISTOLYTICA QUIK CHEK* test is dependent upon the proper reaction of the internal and external controls.

Internal: A vertical blue control line must be visible on the "C" side of the *Reaction Window* on every *Membrane Device* that is tested. This confirms that the sample and reagents were added correctly and reacted properly in the assay. A clear background in the result area is considered an internal negative control. It may appear white to light blue and any developed lines will be clearly visible.

External: The reactivity of the *E. HISTOLYTICA QUIK CHEK* test should be verified on receipt using the *Positive Control* and negative control (*Diluent*). The *Positive Control* confirms the reactivity of the other reagents associated with the assay, and is not intended to ensure precision at the analytical assay cut-off. Additional tests can be performed with the controls to meet the requirements of local, state and/or federal regulations and/or accrediting organizations.

LIMITATIONS

1. A negative test result does not preclude the possibility of the presence of *E. histolytica* adhesin in the specimen which may occur if the level of antigen is below the detection limit of the test.
2. The *E. HISTOLYTICA QUIK CHEK* test is qualitative. The intensity of the color should not be interpreted quantitatively.

EXPECTED VALUES

Normal healthy individuals should not be infected with *E. histolytica* and should test negative in the *E. HISTOLYTICA QUIK CHEK* test. A positive test result in the *E. HISTOLYTICA QUIK CHEK* test indicates that the person is shedding detectable amounts of *E. histolytica* antigen. The incidence of *E. histolytica* infection varies significantly between populations and geographic regions. It is estimated that *Entamoeba histolytica* infects about 50 million people around the world (2). Roughly 90% of these persons remain asymptomatic, whereas about 10% develop clinical symptoms ranging from gastrointestinal disease to liver abscesses. High risk groups include persons who have traveled abroad, immigrants, immunocompromised persons, migrant workers, and active male homosexuals (2, 3). Nonpathogenic strains (*E. dispar*) are predominant among male homosexuals (4). The disease often is transmitted by asymptomatic carriers of *E. histolytica*.

PERFORMANCE CHARACTERISTICS

The performance of the *E. HISTOLYTICA QUIK CHEK* test was compared to the *E. HISTOLYTICA II* test and Remel™ ProSpec™ *Entamoeba histolytica* test with discrepant samples resolved by PCR. Fecal samples included 449 frozen samples. Age information was available for 350 patients. Of the 350 patients, 87% were < 18 years. Sex information was available for 424 patients, and 46% were male and 54% were female. Table 1 shows a summary of the clinical performance of the *E. HISTOLYTICA QUIK CHEK* test at all 3 sites.

Table 1. Summary of clinical performance comparing the *E. HISTOLYTICA* QUIK CHEK test to the *E. HISTOLYTICA* II test, and to the Remel™ ProSpec™ *Entamoeba histolytica* test.

n = 449	<i>E. HISTOLYTICA</i> II Positive	<i>E. HISTOLYTICA</i> II Negative
<i>E. HISTOLYTICA</i> QUIK CHEK Positive	94	3
<i>E. HISTOLYTICA</i> QUIK CHEK Negative	3	349

	95% Confidence Interval	
Percent Positive Agreement	96.9%	90.6% - 99.2%
Percent Negative Agreement	99.1%	97.3% - 99.8%
Overall Percent Agreement	98.7%	98.5% - 98.9%
*Resolved by PCR – Overall Percent Agreement	99.6%	99.6% - 99.6%

n = 449	ProSpec™ <i>E. histolytica</i> Positive	ProSpec™ <i>E. histolytica</i> Negative
<i>E. HISTOLYTICA</i> QUIK CHEK Positive	87	10
<i>E. HISTOLYTICA</i> QUIK CHEK Negative	24	328

	95% Confidence Interval	
Percent Positive Agreement	78.4%	69.4% - 85.4%
Percent Negative Agreement	97.0%	94.5% - 98.5%
Overall Percent Agreement	92.4%	91.0% - 93.6%
*Resolved by PCR – Overall Percent Agreement	99.3%	99.2% - 99.4%

*The 40 discrepant specimens were resolved by purification of DNA from the fecal specimens and polymerase chain reaction (PCR) amplification for *E. histolytica*, *E. dispar*, and *E. moshkovskii*.

For the 6 discrepant specimens between *E. HISTOLYTICA* II and *E. HISTOLYTICA* QUIK CHEK: The 3 specimens positive by *E. HISTOLYTICA* II and negative by *E. HISTOLYTICA* QUIK CHEK included 2 *E. histolytica* PCR-positive and 1 *E. histolytica* PCR-negative. The 3 specimens negative by *E. HISTOLYTICA* II and positive by *E. HISTOLYTICA* QUIK CHEK were *E. histolytica* PCR-positive.

For the 34 discrepant specimens between the Remel™ ProSpec™ *Entamoeba histolytica* test and the *E. HISTOLYTICA* QUIK CHEK: The 24 specimens positive by Remel™ ProSpec™ and negative by *E. HISTOLYTICA* QUIK CHEK included 3 *E. histolytica* PCR-positives, 13 *E. dispar* PCR-positives, and 8 *Entamoeba* PCR-negatives. The 10 specimens negative by Remel™ ProSpec™ and positive by *E. HISTOLYTICA* QUIK CHEK were *E. histolytica* PCR-positive.

REPRODUCIBILITY

The reproducibility of the *E. HISTOLYTICA* QUIK CHEK test was determined using 12 fecal specimens that were coded to prevent their identification during testing. Testing was performed at 2 independent laboratories and on-site at TECHLAB®, Inc. The samples were tested, twice a day over a 5-day period by multiple technicians at each site using 2 different kit lots. A positive and negative control was run with each panel of the masked samples. The results from each laboratory were submitted to TECHLAB®, Inc. and compared with in-house results. The results were consistent among the different locations, and exhibited a correlation of 100%. The samples produced the expected results 100% of the time.

CROSS REACTIVITY

The *E. HISTOLYTICA* QUIK CHEK test was evaluated for cross-reactivity with the bacterial and viral strains listed below. None of the strains were shown to interfere with the performance of *E. HISTOLYTICA* QUIK CHEK test.

<i>Aeromonas hydrophila</i>	<i>Clostridium difficile</i>	<i>Shigella dysenteriae</i>
<i>Bacillus cereus</i>	<i>Enterococcus faecalis</i>	<i>Shigella flexneri</i>
<i>Bacillus subtilis</i>	<i>Escherichia coli</i>	<i>Shigella sonnei</i>
<i>Bacteroides fragilis</i>	<i>Escherichia coli</i> 0157:H7	<i>Staphylococcus aureus</i>
<i>Campylobacter coli</i>	<i>Escherichia coli</i> ETEC	<i>Staphylococcus aureus</i> (Cowan's)
<i>Campylobacter fetus</i>	<i>Escherichia coli</i> EPEC	<i>Staphylococcus epidermidis</i>
<i>Campylobacter jejuni</i>	<i>Escherichia coli</i> EIEC	<i>Vibrio parahaemolyticus</i>
<i>Candida albicans</i>	<i>Klebsiella pneumonia</i>	<i>Yersinia enterocolitica</i>
<i>Clostridium bifermentans</i>	<i>Salmonella typhimurium</i>	

Adenovirus, 2, 5, 40, 41	Human Adenovirus 1, 3	Human Enterovirus 69, 70, 71
Coxsackievirus B5	Human Coronavirus	Human parechovirus 1
Echovirus 11, 18, 22, 33	Human Coxsackievirus B2, B3, B4	[Echovirus 22]
Enterovirus 68, 69	Human Echovirus 9	Human Rotavirus

Additionally, the *E. HISTOLYTICA* QUIK CHEK test was run on fecal specimens documented to be positive for other parasites by microscopy. No cross-reactivity was seen with the following organisms.

<i>Entamoeba coli</i>	<i>Chilomastix mesnili</i>	Hookworm eggs
<i>Entamoeba dispar</i>	<i>Cyclospora</i> spp.	<i>Iodamoeba bütschlii</i>
<i>Entamoeba moshkovskii</i>	<i>Endolimax nana</i>	<i>Trichomonas hominis</i>
<i>Ascaris lumbricoides</i> eggs	<i>Giardia</i> spp.	<i>Trichuris trichiura</i> eggs
<i>Blastocystis hominis</i>	<i>Hymenolepis nana</i>	

INTERFERING SUBSTANCES (U.S. FORMULATION)

The following substances had no effect on positive or negative test results analyzed at the concentrations indicated: Hog gastric mucin (3.5% w/v), Human blood (40% v/v), Hydrocortisone (40% v/v), Barium sulfate (5% w/v), Imodium® (5% v/v), Kaopectate® (5% v/v), Pepto-Bismol® (5% v/v), Maalox® Advanced (5% v/v) Steric Acid, (40% w/v), Palmitic Acid (40% w/v), Metronidazole (0.25% w/v), Phenylephrine (40% w/v), Naproxen Sodium (40% w/v), Nonoxynol-9 (40% v/v), Vancomycin (0.25% w/v), Priolsec OTC® (5 µg/mL), TUMS (50 µg/mL), Tagamet® (5 µg/mL), Ciprofloxacin (0.25% w/v), Mylanta® (4.2 mg/mL), and Human Urine (5% v/v).

PRECISION – INTRA-ASSAY

For the determination of intra-assay performance, twelve human fecal specimens were analyzed by the *E. HISTOLYTICA QUIK CHEK* test. Of these twelve specimens, six were positive for *E. histolytica* of varying levels (low, moderate, and high) and six were negative for *E. histolytica*. Each specimen was assayed five times in the same test run, using two different kit lots. A positive and negative control was run with each panel. All positive samples remained positive and all negative samples remained negative.

PRECISION – INTER-ASSAY

For the determination of inter-assay performance, twelve human fecal specimens were analyzed by the *E. HISTOLYTICA QUIK CHEK* test. Of these twelve specimens, six were positive for *E. histolytica* of varying levels (low, moderate, and high) and six were negative for *E. histolytica*. The samples were tested, twice a day by multiple technicians over a 12-day period using 2 different kit lots. A positive and negative control was run on each day. All positive samples remained positive and all negative samples remained negative.

ANALYTICAL SENSITIVITY

The Limit of Detection (LoD) for the *E. HISTOLYTICA QUIK CHEK* test was established at a concentration of 0.2 ng/mL for *E. histolytica* adhesin.

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© 2014 TECHLAB®, Inc. All rights reserved.
The Alere Logo and Alere are trademarks of the Alere group of companies.
The *E. HISTOLYTICA QUIK CHEK*, TECHLAB Logo and TECHLAB are trademarks of TECHLAB, Inc., under license.
All trademarks referenced are trademarks of their respective owners.
Alere North America, LLC