



What is Amebiasis?




Intestinal invasion by *E. histolytica* trophozoites




Importance of *E. histolytica* diagnostics

- Intestinal colonization by the protozoan parasite *Entamoeba histolytica* is an important cause of diarrheal disease globally, particularly in the developing world.
- Diagnostics targeting protozoan parasites are important tools for surveillance and determination of appropriate therapeutic intervention in clinical settings.
- False positive results lead to inappropriate treatments with added risk to the patient and wasted cost for the health care system.
- In the developing world, limited resources and transitory populations make sensitive and rapid diagnosis of infectious diseases essential for early and effective intervention.

Multi-site evaluation of the *E. HISTOLYTICA* QUIK CHEK (TECHLAB®)



E. HISTOLYTICA QUIK CHEK (TECHLAB®)



RAPID: Results in less than 30 minutes

SIMPLE: Minimal training required

Fresh or frozen samples

CONTROL: Internal positive control in every cassette

E. HISTOLYTICA QUIK CHEK Evaluation Workflow

Fecal Specimens from South Africa and Bangladesh Suspected to Contain *Entamoeba* spp.

ProSpecT <i>E. histolytica</i> Microplate Assay Remel Lenexa, KS	<i>E. HISTOLYTICA</i> QUIK CHEK TECHLAB®, Inc. Blacksburg, VA	<i>E. HISTOLYTICA</i> II ELISA TECHLAB®, Inc. Blacksburg, VA
--	---	--

Comparison of results from 3 assays

PCR for *E. histolytica*, *dispar*, and *moshkovskii* on discrepant samples

Results

449 Total Fecal Specimens were tested

Remel ProSpecT <i>E. histolytica</i> Assay 337 Negative 112 Positive	<i>E. HISTOLYTICA</i> QUIK CHEK 352 Negative 97 Positive	<i>E. HISTOLYTICA</i> II ELISA 349 Negative 100 Positive
--	--	--

Contact Info: Hans Verkerke
 William Petri
 TechLab, Inc.
 hv7a@virginia.edu
 wap3g@cms.mail.virginia.edu
 techlab@techlab.com

Results, Continued

Discrepant Results (N=38)						
PCR for <i>E. histolytica</i> , <i>E. dispar</i> , and <i>E. moshkovskii</i> for discrepant samples						
ProSpecT ELISA discrepant		<i>EH</i> QUIK CHEK discrepant		<i>EH</i> II ELISA discrepant		
31 Samples		3 Samples		4 Samples		
22 +	9 -	0 +	3 -	1 +	3 -	
PCR Result						
EH	0	9	0	3	0	3
ED	13	0	0	0	0	0
EM	0	0	0	0	0	0

PCR for *E. histolytica*, *E. dispar*, and *E. moshkovskii* was performed on 38 discrepant results.

<i>QUIK CHEK</i>	<i>E. HISTOLYTICA</i> II ELISA PCR		ProSpecT	<i>E. HISTOLYTICA</i> II ELISA or PCR	
	Positive	Negative		Positive	Negative
Positive	97	0	Positive	91	21
Negative	3	349	Negative	9	328
	PPA 97%	NPA 100%		PPA 91%	NPA 94%

Results Summary:

- 85 samples tested positive by all three assays
- 327 tested negative by all three assays
- 38 discrepant samples were identified and tested by PCR for *E. histolytica* and other *Entamoeba* spp.
- The *E. HISTOLYTICA* QUIK CHEK exhibited a positive percent agreement (PPA) of 97% and a negative percent agreement (NPA) of 100%
- The ProSpecT™ ELISA exhibited a PPA of 91% and NPA of 94%

Conclusions:

- Our results indicate that the *E. HISTOLYTICA* QUIK CHEK outperforms the ProSpecT™ ELISA in sensitive and specific detection of *E. histolytica* in clinical isolates from the developing world.
- Our analysis of discrepant results suggests that the use of specific antibodies targeting the *E. histolytica* adhesin in the QUIK CHEK enables *Entamoeba* species discrimination—reducing the frequency of false positive results.
- This assay represents a strong candidate for rapid diagnosis of *E. histolytica* infection.

Funding and Conflict of Interest

Funding by NIH grants R42A1058450, R43A1094733 to WAP and JFH as well as AI 43596 to WAP, who receives licensing fees from TECHLAB, Inc., which are donated in full to the Am. Society of Trop. Med. and Hyg. Additional funding provided by a Fulbright Research Award to HV