TECHLAB® GIARDIA/CRYPTOSPORIDIUM CHEK®

A Monoclonal ELISA for Detecting *Giardia* and
Cryptosporidium Antigen in Fecal Specimens
Catalog No. 30401/T30401SKG (96 Tests)

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GIARDIA/CRYPTOSPORIDIUM CHEK®

INTENDED USE

The *GIARDIA/CRYPTOSPORIDIUM CHEK*® test is an enzyme immunoassay for the qualitative detection of *Giardia* cyst and *Cryptosporidium* oocyst antigen in human fecal specimens. It is indicated for use as an aid in the diagnosis of patients with diarrhea suspected of *Giardia* and/or *Cryptosporidium* gastrointestinal infections.
FOR *IN VITRO* DIAGNOSTIC USE.

EXPLANATION

Giardia spp. is a binucleated flagellated protozoan parasite which exists in two forms: a noninfectious, pear-shaped trophozoite (9 to 20 µm) inhabiting the small intestine and the highly infectious cyst form which is elliptical in shape and ranges in size from 8 to 12 µm (1). Survival outside its host varies greatly between the two forms: the trophozoite, which is extremely labile, lasts only a matter of hours outside the body, while the cyst form may survive for several days in an external environment (1). The parasite is responsible for infections due to water contamination and travelers have been found to contract giardiasis from endemic areas (2-5). Transmission also occurs by direct contact especially by asymptomatic carriers and by food contamination (6,7). High-risk categories include young children, immunocompromised patients and those without previous exposure (8). More recently, giardiasis has become a common sexually transmitted disease (9). Animal fecal contamination especially of water is another route of transmission in humans (4,10,11).

Clinical manifestations of giardiasis range from asymptomatic carriage with the passing of cysts to chronic debilitating diarrhea, weight loss and malabsorption (8,12,13). Diagnosis by microscopy has been the most common method used for giardiasis. However, this procedure requires extensive experience and the presence of intact cysts in the feces. An alternative test is the ELISA. This procedure is fairly simple to perform and exhibits increased sensitivity compared to microscopic examination. With microscopy, fecal specimens are examined using various procedures (14) and accuracy of results is dependent on the skill of the technician. The success rate of fecal examination varies between 50 and 70% (15, 16) and several specimens are usually necessary to establish a diagnosis. Giardiasis also may be present in the absence of whole organisms and can be confused with other illnesses such as Crohn's disease and ulcerative colitis (17). When infection is present but parasites are not detected, sampling and testing of duodenal fluid may detect trophozoites, but this method is invasive and expensive (14,16). Detection of the organism and antigens by ELISA provides an alternative method of establishing a diagnosis and is sensitive and specific (18-20). Large numbers of specimens can be tested rapidly and objectively, and the procedure is less labor intensive than most microscopy methods.

Cryptosporidium spp. is a protozoan parasite of vertebrates previously thought to cause diarrhea only in animals (23). In 1976, the first human infection was reported (24). Since then it has been found to be associated with diarrheal illness in most parts of the world and is a frequent cause of travellers' diarrhea (23,25). The disease is transmitted by the thick-walled oocyst form, 2-6 µm in diameter, which is remarkably resistant to common disinfectants and routine chlorination of drinking water. Person-to-person transmission, especially among children, is common (26). *Cryptosporidium* has little or no host specificity and animals such as rodents, cattle and domestic pets serve as a reservoir for zoonotic transmission to humans (23,27). This occurs either by direct contact or by contamination of water supplies with fecal matter (23,28-30). Cryptosporidiosis is a serious opportunistic infection in patients with acquired immunodeficiency syndrome (AIDS) and is potentially a sexually transmitted disease (28,31).

Clinical manifestations of cryptosporidiosis include cholera-like diarrhea, abdominal pain, nausea, vomiting and weight loss. In normal persons the infection is usually self-limiting and short term. In AIDS and immunosuppressed patients, cryptosporidiosis can

result in prolonged and life-threatening illness due to excessive fluid loss. In these patients the infection may also spread to the respiratory and biliary tracts (28).

BIOLOGICAL PRINCIPLES

The *GIARDIA/CRYPTOSPORIDIUM CHEK*[®] test uses monoclonal and polyclonal antibodies to cell-surface antigens of *Giardia* and an oocyst antigen of *Cryptosporidium*. The *Microassay Plate* in the kit contains immobilized monoclonal antibodies against the antigens and the *Conjugate* consists of polyclonal antibodies against the antigens. In the assay, an aliquot of a diluted fecal specimen is transferred to a microassay well. The immobilized monoclonal antibodies bind the *Giardia* and *Cryptosporidium* antigens if the antigens are present. Upon addition, *Conjugate* then binds to the antigen/antibody complex. Any unbound materials are removed during the washing steps. Following the addition of substrate, a color is detected due to the enzyme-antibody-antigen complexes that formed in the presence of *Giardia* and/or *Cryptosporidium* antigens and *Conjugate*.

MATERIALS PROVIDED

CONJ	ENZ	Conjugate , 7 mL (Rabbit polyclonal antibody to a cell-surface antigen of <i>Giardia</i> and rabbit polyclonal antibody to a cell-surface antigen of <i>Cryptosporidium</i> in a buffered protein solution containing 0.02% thimerosal)
DIL	SPE	Diluent , 50 mL (Buffered protein solution containing 0.02% thimerosal). The <i>Diluent</i> is also to be used as the negative control solution (see TEST PROCEDURE).
H ₂ SO ₄	0.6N	Stop Solution , 7 mL (0.6N sulfuric acid). Caution: Avoid contact with skin. Flush with water immediately if contact occurs.
CONTROL	+	Giardia Positive Control , 3.5 mL (<i>Giardia</i> antigen in a buffered protein solution containing 0.02% thimerosal)
CONTROL	+	Cryptosporidium Positive Control , 3.5 mL (Heat-inactivated bovine fecal material containing <i>Cryptosporidium</i> antigen in a buffered protein solution with 0.02% thimerosal)
SUBS	REAG	Substrate , 14 mL (solution containing tetramethylbenzidine and peroxide)
WASHBUF	20X	Wash Buffer Concentrate , 50 mL (20X concentrate containing phosphate-buffered saline, detergent and 0.2% thimerosal)
MA	PLT	Microassay Plate , 12 strips, each consisting of 8 wells coated with monoclonal antibodies to <i>Giardia</i> cell-surface antigen and monoclonal antibodies to <i>Cryptosporidium</i> oocyst antigen (stored with desiccant)
		2 plastic adhesive sheets
		100 graduated disposable pipettes

MATERIALS REQUIRED BUT NOT PROVIDED

Squirt bottle for wash reagent	Vortex mixer
950 mL distilled water for diluting wash reagent	Discard container
Absorbent paper	Applicator sticks
Small test tubes (e.g., microcentrifuge tubes)	
ELISA reader capable of reading at 450 nm or 450/620 nm	

PRECAUTIONS

1. Reagents from different kits should not be mixed. Do not use a kit past the assigned expiration date.
2. Reagents should be brought to room temperature before use.
3. Caps and tips are color-coded; do not mix.
4. Each component in the kit should be inspected for any signs of leakage. Upon arrival, the kit should be inspected to ensure that components are not frozen or warm to the touch due to improper shipping conditions.
5. When handling assay wells, avoid scratching the bottom of the wells because this may result in elevated absorbance readings.
6. Hold dropper bottles vertically when dispensing to ensure proper drop size.
7. Microassay wells should be handled and disposed of as potential biohazards after use. Wear gloves when doing the test.

- 4
8. Reagents contain thimerosal as a preservative and should be handled with normal laboratory caution.
 9. The *Stop Solution* contains 0.6 N sulfuric acid. Flush with water immediately if contact with skin or eyes occurs.
 10. Unused microassay wells must be placed back inside the resealable pouch with the desiccant to protect them from moisture.
 11. Perform the washing procedure as directed to avoid high background reactions.
 12. The *Substrate* is light sensitive and should be protected from direct sunlight or UV sources. If the *Substrate* is exposed to light and develops a color, it must be discarded.
 13. Optimal results are obtained by following the specified test procedure. The concentrations, incubation conditions, and processing specifications have been optimized for sensitivity and specificity. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test.
 14. Microbial contamination of reagents may decrease the accuracy of the assay. Avoid microbial contamination of reagents by using sterile disposable pipettes if removing aliquots from reagent bottles.

SHELF LIFE AND STORAGE

The expiration date of the kit is given on the box label. Expiration dates for each component are listed on individual labels. The kit should be stored between 2° and 8°C and should be returned to the refrigerator as soon as possible after use.

SPECIMEN HANDLING

1. Standard collection and handling procedures used in-house for fecal specimens for culture are appropriate. No modification of collection methods used for standard microscopic O&P examinations is needed. Fecal specimens may be used as unpreserved or frozen, or in preservation media of 10% buffered formalin or Sodium Acetate Formalin (SAF).
2. Unpreserved specimens should be kept between 2° and 8°C and tested within 24 hours of collection. Specimens that cannot be tested within this time should be stored at -20°C or less until tested.
3. Preserved specimens may be kept at room temperature and tested within 18 months of collection.
4. Concentration steps are not needed (or recommended) for fecal specimens.
5. Make sure that specimens are thoroughly mixed (vortexed) prior to performing the test. This includes complete mixing of the specimen prior to transfer to the *Diluent* and/or microassay well.
6. All dilutions must be made with *Diluent*.

REAGENT PREPARATION

1. The contents of the kit should be brought to room temperature before use.
2. Prepare 1X *Wash Solution*. The *Wash Buffer Concentrate* is supplied as a 20X concentrate (a precipitate may be noticed). It should be diluted to a total volume of 1 liter by adding 50 mL of the concentrate to 950 mL of distilled water. The 1X *Wash Solution* can be stored between 2° and 8°C.
3. All reagents, with the exception of the *Wash Buffer Concentrate*, are supplied in ready-to-use bottles. Reagents can be dispensed directly from the dropper bottles or decanted for use with multichannel pipettes. If excess reagent has been decanted, the excess should be discarded. Do not pour back into the bottle. The *Substrate* should be stored in and used from the light-protected bottle in which it is supplied. If an aliquot is removed from the original bottle for any reason, do not return unused *Substrate* to the original bottle.

SPECIMEN PREPARATION

1. Fresh/Frozen Fecals: Frozen fecal specimens should be thawed. Add 400 µL of *Diluent* to a microcentrifuge tube (one tube per sample), then add 100 µL (2nd

graduation mark on pipette) of fecal sample to the tube and mix well. If the specimen cannot be pipetted, use an applicator stick to transfer approximately 0.1 gram of feces. This is about the size of a small pea (about 4 mm in diameter).

2. Preserved Fecals: Mix (vortex) contents of container thoroughly before transferring specimen. No further processing or dilution is necessary.

TEST PROCEDURE

1. Prepare three control wells to serve as positive and negative controls. To the first positive control well, add 1 drop (50 μ L) of *Giardia Positive Control* (black cap), to the second positive control well add 1 drop (50 μ L) of *Cryptosporidium Positive Control* (white cap), and to the negative control well add 2 drops (100 μ L) of *Diluent*.
2. If fresh or frozen samples (i.e. unpreserved) are to be used, make a dilution as stated in SPECIMEN PREPARATION above.
3. Transfer 100 μ L of *Diluent* to each test well of the *Microassay Plate*. Using plastic pipettes, transfer 1 drop (50 μ L, 1st graduation mark on pipette) of sample (preserved or diluted as above) to each test well already containing *Diluent* and gently tap the wells to mix. Seal with a plastic adhesive sheet and incubate for 1 hour at room temperature.
4. Shake out contents of assay wells into a discard pan. Wash each well using the 1X *Wash Solution* in a squirt bottle with a fine-tipped nozzle, directing the *Wash Solution* to the bottom of the well with force. Fill the wells, shake the wash solution out of the well into a discard pan and slap the plate hard onto a dry paper towel. Repeat the washing step 4X (for a total of 5 washes). If any fecal material remains in the wells, wash the plate until it appears clean. (Go to www.techlab.com for additional information on ELISA well washing)
5. After washing, completely remove any residual liquid in the wells by striking the plate onto a dry paper towel until no liquid comes out. Dispose of paper towels and specimen containers properly.
6. Add 1 drop (50 μ L) of *Conjugate* (red cap) to each well and gently tap to mix. Seal with a plastic adhesive sheet. Incubate the wells for 30 minutes at room temperature.
7. Repeat the washing procedure (Steps 4 and 5).
8. Add 2 drops (100 μ L) of *Substrate* (blue cap) to each well. Gently tap the wells to mix. Incubate the wells at room temperature for 10 minutes.
9. Add 1 drop (50 μ L) of *Stop Solution* (yellow cap) to each well. Gently tap the wells to mix and wait 2 minutes before reading. The addition of the *Stop Solution* converts the blue color to a yellow color, which may be quantitated by measuring the absorbance at 450 nm on a microplate ELISA reader. The instrument should be blanked against air. Wipe the underside of each well to remove moisture before measuring the absorbance. If a dual reader is used, blank against air at 620 nm and read at 450 nm. Visual readings should be noted. Read within ten minutes after adding *Stop Solution*.

QUALITY CONTROL

1. A *Giardia Positive Control*, *Cryptosporidium Positive Control*, and a negative control must be run with each series of test specimens.
2. Each positive control well should be an easily visible yellow color and should give an absorbance of 0.500 or higher. Any well that gives a positive reading without visible color should be repositioned, wiped on the underside of the well to remove moisture, and read again.
3. Negative control wells should be colorless or they may have a faint yellow color but they must have an absorbance value of <0.150 when read on single wavelength (OD_{450nm}) or <0.090 when read on dual wavelength ($OD_{450/620nm}$).
4. Test results are not valid unless the performance characteristics of the positive and negative controls are met. If these results are not observed, call Technical Services.

6

5. Test results along with control absorbance values should be recorded and reported according to in-house procedures and should be stored according to in-house procedures for future reference.

INTERPRETATION OF RESULTS

Absorbance

<u>450 nm</u>	<u>450/620 nm</u>	<u>Visual Color</u>	<u>Interpretation</u>
< 0.150	< 0.090	Clear to slight yellow	Negative - Antigen is not detected. Below the detectable limits of the assay.
≥ 0.150	≥ 0.090	Pale yellow to strong yellow	Positive - Specimen contains <i>Giardia</i> antigen and/or <i>Cryptosporidium</i> antigen

Visual Interpretation

Negative: Any sample that is colorless or resembles the negative control well in intensity of color. Antigen is not detected. Below the detectable limits of the assay.

Positive: Any sample that is obviously more yellow than the negative control well. Specimen contains *Giardia* antigen and/or *Cryptosporidium* antigen.

NOTE: The negative control, as well as some test wells, may show some slight yellow color. A sample well must be obviously more yellow than the negative control well to be called a positive result.

Spectrophotometric Interpretation

1. Determine the absorbance value of the negative control.
The negative control reading should be < 0.150 OD₄₅₀ or < 0.090 OD_{450/620}. If not, the test is invalid and should be repeated, paying attention to the wash procedure.
2. The reading for the *Positive Control* should be ≥ 0.500.
3. Test results
Negative: < 0.150 (absorbance at 450 nm) or < 0.090 (absorbance at 450/620 nm).
Antigen is not detected. Below the detectable limits of the assay.
Positive: ≥ 0.150 (absorbance at 450 nm) or ≥ 0.090 (absorbance at 450/620 nm).
Specimen contains *Giardia* antigen and/or *Cryptosporidium* antigen.

LIMITATIONS OF THE PROCEDURE

1. The *GIARDIA/CRYPTOSPORIDIUM CHEK*® test detects the presence of *Giardia* and *Cryptosporidium* antigen in fecal specimens.
2. The test results should be interpreted by a physician in consideration of other laboratory results and clinical history.
3. Concentrated fecal specimens should not be tested in the *GIARDIA/CRYPTOSPORIDIUM CHEK*® test and will not give accurate results.
4. The predictive value of a positive result decreases when testing in a low prevalence population.
5. The *GIARDIA/CRYPTOSPORIDIUM CHEK*® test is for the qualitative detection of *Giardia* and *Cryptosporidium* antigen in fecal specimens. It has not been evaluated for quantitative determinations of organism load, and the magnitude of the absorbance value is not intended to correlate with organism load.

EXPECTED VALUES

Normal healthy individuals should not be infected with *Giardia* or *Cryptosporidium* and should test negative in the *GIARDIA/CRYPTOSPORIDIUM CHEK*® test. A positive test result in the *GIARDIA/CRYPTOSPORIDIUM CHEK*® test indicates that the person is shedding detectable amounts of *Giardia* or *Cryptosporidium* antigen. The incidence of *Giardia* and *Cryptosporidium* infection varies significantly between populations and geographic regions. Children in daycare settings have exhibited higher rates of infection with *Giardia* than the normal population (21). In addition, homosexual men have shown higher rates of infection (22). In general, laboratory-confirmed incidence of cryptosporidiosis in developed countries ranges from 1 to 2% overall with a higher incidence in children (32).

PERFORMANCE CHARACTERISTICS

An independent diagnostics laboratory evaluated the *GIARDIA/CRYPTOSPORIDIUM CHEK*® test for its ability to detect *Giardia* and *Cryptosporidium* antigens in preserved fecal specimens. Of the 217 specimens analyzed by microscopy, 74 were positive for *Giardia* and 44 were positive for *Cryptosporidium*. Three specimens were positive for both *Giardia* and *Cryptosporidium*. Each specimen was subsequently tested using the *GIARDIA/CRYPTOSPORIDIUM CHEK*® test. Results are summarized in Tables 1, 2 and 3.

Table 1. Comparison of the *GIARDIA/CRYPTOSPORIDIUM CHEK*® test with microscopy for detection of *Giardia* and/or *Cryptosporidium*.

Total Samples (n = 217)		Microscopy	
		Positive	Negative
<i>GIARDIA/CRYPTOSPORIDIUM CHEK</i> ®	Positive	121	0
	Negative	3	93
		95% CI	
Sensitivity	97.6%	92.6% - 99.4%	
Specificity	100%	95.1% - 100%	
Predictive Positive Value	100%	96.2% - 100%	
Predictive Negative Value	96.9%	90.5% - 99.2%	
Correlation	98.6%	98.2% - 98.9%	

Table 2. Comparison of the *GIARDIA/CRYPTOSPORIDIUM CHEK*® test with microscopy for the detection of *Giardia* only. All samples in this data set have been shown to be negative for *Cryptosporidium* by microscopy.

Total Samples (n = 167)		Microscopy - <i>Giardia</i>	
		Positive	Negative
<i>GIARDIA/CRYPTOSPORIDIUM CHEK</i> ®	Positive	74	0
	Negative	0	93
		95% CI	
Sensitivity	100%	93.9% - 100%	
Specificity	100%	95.1% - 100%	
Predictive Positive Value	100%	93.9% - 100%	
Predictive Negative Value	100%	95.1% - 100%	
Correlation	100%	100% - 100%	

Table 3. Comparison of the *GIARDIA/CRYPTOSPORIDIUM CHEK*® test with microscopy for the detection of *Cryptosporidium* only. All samples in this data set have been shown to be negative for *Giardia* by microscopy.

Total Samples (n = 140)		Microscopy - <i>Cryptosporidium</i>	
		Positive	Negative
<i>GIARDIA/CRYPTOSPORIDIUM CHEK</i> ®	Positive	44	0
	Negative	3	93
		95% CI	
Sensitivity	93.6%	81.4% - 98.3%	
Specificity	100%	95.1% - 100%	
Predictive Positive Value	100%	90.0% - 100%	
Predictive Negative Value	96.9%	90.5% - 99.2%	
Correlation	97.9%	97.1% - 98.4%	

CROSS REACTIVITY

An independent diagnostics laboratory evaluated the *GIARDIA/CRYPTOSPORIDIUM CHEK*[®] test using fecal specimens found to be positive for a variety of intestinal pathogens. No cross reactivity was observed with fecal specimens that contained any of the pathogens listed below. The number of specimens tested with each organism is shown in parentheses.

<i>Ascaris lumbricoides</i> eggs (26)	<i>Entamoeba coli</i> (17)
<i>Blastocystis hominis</i> (31)	<i>Entamoeba hartmanni</i> (4)
<i>Chilomastix mesnili</i> (2)	<i>Enterobius vermicularis</i> eggs (8)
<i>Cyclospora cayetanensis</i> (1)	<i>Hymenolepis nana</i> eggs (4)
<i>Dientamoeba fragilis</i> (10)	<i>Iodamoeba bütschlii</i> (4)
<i>Diphyllobothrium latum</i> eggs (1)	<i>Strongyloides stercoralis</i> larvae (2)
<i>Endolimax nana</i> (36)	<i>Taenia</i> spp. eggs (2)
<i>Entamoeba histolytica/dispar</i> (9)	<i>Trichuris trichiura</i> eggs (20)

The *GIARDIA/CRYPTOSPORIDIUM CHEK*[®] test was evaluated for crossreactivity with the bacterial and viral strains listed below. None of the strains were shown to crossreact with the *GIARDIA/CRYPTOSPORIDIUM CHEK*[®] test.

<i>Escherichia coli</i>	<i>Escherichia coli</i> ETEC (enterotoxic)	
<i>Escherichia coli</i> 0157H7	<i>Escherichia coli</i> EPEC (enteropathogenic)	
<i>Yersinia enterocolitica</i>	<i>Escherichia coli</i> EIEC (enteroinvasive)	
<i>Aeromonas hydrophila</i>	<i>Salmonella typhimurium</i>	
<i>Shigella dysenteriae</i>	<i>Shigella flexneri</i>	
<i>Salmonella typhimurium</i>	<i>Campylobacter coli</i>	
<i>Campylobacter fetus</i>	<i>Clostridium difficile</i>	
<i>Vibrio parahaemolyticus</i>	<i>Staphylococcus aureus</i> (Cowan's)	
<i>Staphylococcus aureus</i>	<i>Staphylococcus epidermidis</i>	
<i>Klebsiella pneumoniae</i>	<i>Enterococcus faecalis</i>	
<i>Clostridium bifermentans</i>	<i>Bacteroides fragilis</i>	
<i>Bacillus subtilis</i>	<i>Bacillus cereus</i>	
Adenovirus type 1	Adenovirus type 2	Adenovirus type 3
Adenovirus type 5	Adenovirus type 40	Adenovirus type 41
Human coronavirus	Coxsackievirus B2	Coxsackievirus B3
Coxsackievirus B4	Coxsackievirus B5	Echovirus 9
Echovirus 11	Echovirus 18	Echovirus 22
Echovirus 33	Enterovirus type 68	Enterovirus type 69
Enterovirus type 70	Enterovirus type 71	

INTERFERING SUBSTANCES

The following substances had no effect on positive or negative test results analyzed at the concentrations indicated: mucin (3.5% w/v), human blood (40% w/v), Imodium[®] (5% w/v), Kaopectate[®] (5 mg/mL), Pepto-Bismol[®] (5% w/v), fecal fat (stearic acid - 40% w/v), Metronidazole (0.25% w/v), Vancomycin (0.25% w/v).

SPECIMEN PRESERVATIVE COMPATIBILITY

The *GIARDIA/CRYPTOSPORIDIUM CHEK*[®] test was examined for compatibility with preserved fecal specimens during an in-house evaluation. No decrease in sensitivity or specificity was observed with specimens preserved in 10% buffered formalin or Sodium Acetate Formalin (SAF) when compared to fresh and frozen specimens.

PRECISION – INTRA-ASSAY

The intra-assay % Coefficient of Variation (CV) was determined by analyzing six positive fecal specimens (three positive for *Giardia* and three positive for *Cryptosporidium*) and two negative fecal specimens. Each specimen was assayed in twelve test wells. The % CV of positive specimens ranged from 1.5 to 12.2, with an average of 6.3. The % CV of negative specimens ranged from 3.2 to 3.5, with an average of 3.4.

PRECISION - INTER-ASSAY

For the determination of inter-assay performance, six positive fecal specimens (three positive for *Giardia* and three positive for *Cryptosporidium*) and two negative fecal specimens were assayed a total of five times over a five-day period using the *GIARDIA/CRYPTOSPORIDIUM CHEK*® kit. The % Coefficient of Variation (CV) of positive specimens ranged from 0.3 to 18.1, with an average % CV of 7.4. The % CV of negative specimens ranged from 4.1 to 7.7, with an average of 5.9.

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REFERENCES

1. Erlandsen, L. S. and E. A. Meyer (Ed.). 1984. *Giardia* and giardiasis: biology, pathogenesis and epidemiology. Plenum Press, New York.
2. World Health Organization Scientific Working Group. 1980. Parasite-related diarrhoeas. Bull W.H.O. 58:819-830.
3. Petersen, H. 1972. Giardiasis (lambliasis). Scand. J. Gastroenterol. 7 (Suppl 14):1-44.
4. Craun, G. F. 1979. Waterborne outbreaks of giardiasis. In Jakubowski, W., J. Hoff (Ed.), Waterborne transmission of giardiasis. U.S. Environmental Protection Agency, Washington, D.C., pp. 127-149.
5. Brodsky, R. E., H. C. Spencer, and M. G. Schultz. 1974. Giardiasis in American travelers to the Soviet Union. J. Infect. Dis. 130:319-23.
6. White, K. E., C. W. Hedberg, L. M. Edmons, D. B. W. Jones, M. T. Osterholm, and K. L. MacDonald. 1989. An outbreak of giardiasis in a nursing home with evidence for multiple modes of transmission. J. Infect. Dis. 160:298-304.
7. Pickering, L. K., W. E. Woodward, H. L. DuPont, and P. Sullivan. 1984. Occurrence of *Giardia lamblia* in day care centres. J. Ped. 104:522-526.
8. Stevens, D. P., and A. A. Mahmoud. 1980. Giardiasis: the rediscovery of an ancient pathogen. Curr. Clin. Top. in Infect. Dis. 1:195-207.
9. Phillips, S. C., D. Mildvan, and D. C. Williams. 1981. Sexual transmission of enteric protozoa and helminths in a venereal disease clinic population. New Eng. J. Med. 305:603-606.
10. Friend, D. S. 1966. The fine structure of *Giardia muris*. J. Cell Biol. 29:317-332.
11. Faubent, G. M. 1988. Evidence that giardiasis is a zoonosis. Parasitol. Today 4:66-89.
12. Raizman, R. E. 1976. Giardiasis: an overview for the clinician. Digest. Dis. 21:70-74.
13. Kay, R., G. L. Barnes, and R. R. W. Townley. 1977. *Giardia lamblia* infestation in 154 children. Aust. Paed. J. 13:98-104.
14. Sun, T. 1980. The diagnosis of giardiasis. Am. J. Surg. Path. 4:265-271.
15. Burke, J. A. 1977. The clinical and laboratory diagnosis of giardiasis. Crit. Rev. Clin. Lab. Sci. 7:373-391.
16. Kamath, K. R., and R. Murugasu. 1974. A comparative study of four methods for detecting *Giardia lamblia* in children with diarrhoeal disease and malabsorption. Gastroenterology 66:16-21.
17. Allison, M. C., E. L. Green, D. N. Bhattacharya, A. Smith, and R. E. Ponder. 1988. A microscopic and immunodiagnostic search for giardiasis in patients with gastrointestinal disorders. Scan. J. Gastroenterol. 23:209-212.
18. Nash, T. E., D. A. Herrington, and M. M. Levine. 1987. Usefulness of an enzyme-linked immunosorbent assay for detection of *Giardia* antigen in faeces. J. Clin. Microbiol. 27:1169-1174.
19. Janoff, E. N., J. C. Croft, L. K. Pickering, T. Novotny, M. J. Blaser, C. V. Knisley, and L. B. Reller. 1989. Diagnosis of *Giardia lamblia* infections by detection of parasite specific antigens. J. Clin. Microbiol. 27:431-435.
20. Stibbs, H. 1989. Monoclonal antibody-based enzyme immunoassay for *Giardia lamblia* antigen in human stool. J. Clin. Microbiol. 27:2582-2588.
21. Novotny, T. E., R. S. Hopkins, P. Shillam, and E. N. Janoff. 1990. Prevalence of *Giardia lamblia* and risk factors for infection among children attending day-care facilities in Denver. Public Health Rep. 105(1): 72-5.
22. William, D.C. 1981. Enteric Diseases. Cutis. 27(3):278-81, 283-5.
23. Fayer, R. and L. P. Ungar. 1986. *Cryptosporidium* spp and Cryptosporidiosis. Micro. Rev. 50:458-483.
24. Meisel, J. L., D. R. Perera, C. Meligro, and C. E. Rubin. 1976. Overwhelming watery diarrhea associated with *Cryptosporidium* in an immunosuppressed patient. Gastroenterology 70:1156-60.
25. Sterling, C. R. 1986. *Cryptosporidium* as a causative agent of travellers' diarrhoea. J. Infect. Dis. 153:380-1.
26. Alpert, G., L. M. Bell, C. E. Kirkpatrick, L. D. Budnick, J. M. Campos, H. M. Friedman, and S. A. Plotkin. 1984. Cryptosporidiosis in a day-care centre. New Eng. J. Med. 311:860-1.
27. Pitlik, S. D., V. Fainstein, D. Garza, R. Bolivar, A. Rios, R. L. Hopfer, and P. A. Mansell. 1983. Human cryptosporidiosis: spectrum of disease (1983). Arch. Intern. Med. 143:2269-74.

28. Current, W. L. 1989. Cryptosporidiosis. In: New Strategies in Parasitology.(Ed. K. P. W. J. McAdam) Churchill Livingstone pp 257-73.
29. Hayes, E. B., T. D. Matter, T. R. O'Brien, T. W. McKinley, G. S. Logsdon et al. 1989. Large community outbreak of Cryptosporidiosis due to contamination of a filtered public water supply. New Eng. J. Med. 320:1372-76.
30. Badenoch, J. 1990. *Cryptosporidium* in water supplies. London H.M.S.O. pp 37-45.
31. Angus, K. W. 1990. Cryptosporidiosis and AIDS. Clinical Gastroenterol. 4:425-41.
32. Fayer, R. 1997. *Cryptosporidium* and Cryptosporidiosis. CRC Press, New York.

FOR INFORMATIONAL USE
ONLY

Technical Support

Advice Line

Further information can be obtained from your distributor, or by contacting Alere

Technical Support on:

US	+ 1 877 866 9335	TS.SCR@alere.com
Africa, Russia, CIS	+972 8 9429 683	ARCISproductsupport@alere.com
Asia Pacific	+61 7 3363 7711	APproductsupport@alere.com
Canada	+1 800 818 8335	CANproductsupport@alere.com
Europe & Middle East	+44 161 483 9032	EMEproductsupport@alere.com
Latin America	+57 01800 094 9393	LAPproductsupport@alere.com

The Alere Logo and Alere are trademarks of the Alere group of companies.

GIARDIA/CRYPTOSPORIDIUM CHEK and *TECHLAB* are trademarks of *TECHLAB*®, Inc. under license.