## REFERENCES

1. Lyerly, D. M., H. C. Krivan, and T. D. Wilkins. 1988. *Clostridium difficile*: its disease and toxins. Clin. Microbiol. Rev. 1:1-18.
2. Lyerly, D. M., K. E. Saum, D. K. MacDonald, and T. D. Wilkins. 1985. Effects of *Clostridium difficile* toxins given intragastrically to animals. Infect. Immun. 47: 349-352.
3. Borriello, S. P., F. E. Barclay, A. R. Welch, J. M. Ketley, T. J. Mitchell, J. Stephen, and G. E. Griffin. 1985. Host and microbial determinants of the spectrum of *Clostridium difficile* mediated gastrointestinal disorders. Microecol. Ther. 15:231-236.
4. Lyerly, D. M., N. M. Sullivan, and T. D. Wilkins. 1983. Enzyme-linked immunosorbent assay for *Clostridium difficile* toxin A. J. Clin. Microbiol. 17:72-78.
5. Laughon, B. E., R. P. Viscidi, S. L. Gdovin, R. H. Yolken, and J. G. Bartlett. 1984. Enzyme immunoassays for detection of *Clostridium difficile* toxins A and B in fecal specimens. J. Infect. Dis. 149: 781-788.
6. Lyerly, D. M., L. A. Barroso, and T. D. Wilkins. 1992. Characterization of a toxin A-/toxin B+ isolate of *Clostridium difficile*. Infect. Immun. 60: 4633-4639.
7. Dove, C. H., S.-Z. Wang, S. B. Price, C. J. Phelps, D. M. Lyerly, T. D. Wilkins, and J. L. Johnson. 1990. Molecular characterization of the *Clostridium difficile* toxin A gene. Infect. Immun. 58: 480-488.

## Technical Support

Further information can be obtained from contacting TECHLAB® Technical Support:

US + 1 800 TECHLAB  
 Phone (540) 953-1664  
 FAX (540) 953-1665

TOX A/B QUIK CHEK, the TECHLAB Logo and TECHLAB are trademarks of TECHLAB, Inc.

All other trademarks referenced are trademarks of their respective owners.

© 2018 TECHLAB, Inc. All rights reserved.

RMS #92-033-02

Issued: 06/2018

# TOX A/B QUIK CHEK®

A rapid test for the detection of  
*C. difficile* toxins A and B in fecal specimens  
 Catalog No. T5033 (25 Tests)  
 IVD In Vitro Diagnostic Medical Device  
 For Canadian Users: For Laboratory Use Only

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
## TOX A/B QUIK CHEK®

### INTENDED USE

The *TOX A/B QUIK CHEK*® test is a rapid immunoassay for detecting *Clostridium difficile* toxins A and B in fecal specimens from persons suspected of having *C. difficile* disease. The test is to be used as an aid in the diagnosis of *C. difficile* disease and results should be considered in conjunction with the patient history.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician.

### EXPLANATION

After treatment with antibiotics, many patients develop gastrointestinal problems ranging from mild diarrhea to severe pseudomembranous colitis. Many cases of the milder forms of gastrointestinal illness and most cases of pseudomembranous colitis are caused by toxigenic strains of *Clostridium difficile* (1). This organism is an opportunistic anaerobic bacterium that grows in the intestine once the normal flora has been altered by the antibiotic. Toxigenic strains of *C. difficile* carry the genes encoding the toxins while non-toxigenic strains do not carry the toxin genes. The disease results from the toxins that the toxigenic organism produces. The clinical symptoms associated with the disease are believed to be primarily due to toxin A, which is a tissue-damaging enterotoxin (2,3). *C. difficile* also produces a second toxin, designated toxin B. Toxin B, which has been referred to as the cytotoxin of the organism, is the toxin detected by the tissue culture assay currently used by many laboratories. Toxigenic *C. difficile* strains produce both toxins or only toxin B (4-7).

### PRINCIPLE OF THE TEST

The *TOX A/B QUIK CHEK*® uses antibodies specific for toxins A and B of *C. difficile*. The device contains a *Reaction Window* with two vertical lines of immobilized antibodies (Fig. 1a). The test line ("T") contains antibodies against *C. difficile* toxins A and B. The control line ("C") contains anti-IgG antibodies. The *Conjugate* consists of antibodies to toxins A and B coupled to horseradish peroxidase. To perform the test, the sample is added to a tube containing a mixture of *Diluent* and *Conjugate*. The diluted sample-conjugate mixture is added to the *Sample Well* and the device is allowed to incubate at room temperature for 15 minutes. During the incubation, any toxin A and toxin B in the sample bind to anti-toxin antibody-peroxidase conjugate. The toxin-antibody complexes migrate through a filter pad to a membrane where they are captured by the immobilized anti-toxin antibodies in the line. The *Reaction Window* is subsequently washed with *Wash Buffer*, and the test is developed with the addition of *Substrate*. After a 10 minute incubation period, the "T" reaction is examined visually for the appearance of a vertical blue line on the "T" side of the *Reaction Window*. A blue line indicates a positive test. A positive "C" reaction, indicated by a vertical blue line on the "C" side of the *Reaction Window*, confirms that the test is working properly and the results are valid.

### MATERIALS PROVIDED

MEM	DEV	<b>Membrane Devices</b> –25 pouches, each containing 1 device
DIL	SPE	<b>Diluent (14 mL)</b> – Buffered protein solution containing 0.02% thimerosal with graduated dropper assembly*
WASH	REAG	<b>Wash Buffer (10 mL)</b> – A buffered solution containing 0.02% thimerosal with graduated dropper assembly*
SUBS	REAG	<b>Substrate (3.5 mL)</b> – Solution containing tetramethylbenzidine
CONJ	ENZ	<b>Conjugate (2 mL)</b> – Mouse monoclonal antibody specific for toxin A coupled to horseradish peroxidase and goat polyclonal antibody specific for toxin B coupled to horseradish peroxidase in a buffered protein solution containing 0.02% thimerosal*
CONTROL	+	<b>Positive Control (1 mL)</b> – Antigen in a buffered protein solution
<b>Disposable plastic transfer pipettes</b> – 50 (graduated at 25 µL, 400 µL and 500 µL)		

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
\*contains mercury

#### MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

*Small test tubes (e.g., plastic Eppendorf tubes)*  
*Timer*  
*Disposable gloves for handling fecal samples*

*Applicator sticks*  
*Vortex mixer*  
*Pipettor and tips*

#### SHELF LIFE AND STORAGE

The expiration date of the kit is given on the label. Expiration dates for each component are listed on the individual labels. The kit should be stored between 2° and 8°C.

#### PRECAUTIONS

- Rx Only - Prescription Only
- Reagents from different kits should not be mixed. Do not use a kit past the expiration date.
- Bring all components to ROOM TEMPERATURE BEFORE USE!
- Caps, tips and dropper assemblies are color-coded; do NOT mix or interchange!
- Do not freeze the reagents. The kit should be stored between 2°C and 8°C.
- The pouch containing the *Membrane Device* should be at room temperature before opening, and opened just before use. Keep the membrane devices dry before use.
- Use fecal specimens within 72 hours of collection to obtain optimal results. Specimens that are frozen may lose activity due to freezing and thawing.
- Specimens that have been preserved in 10% Formalin, merthiolate formalin, sodium acetate formalin or polyvinyl alcohol cannot be used.
- Specimens in transport media such as Cary Blair and C&S can be used as specified in the specimen preparation protocol.
- Hold reagent bottles vertically to dispense reagents to ensure consistent drop size.
- Specimens and membrane devices should be handled and disposed of as potential biohazards after use. Wear disposable gloves when doing the test.
- Membrane devices cannot be reused.
- The test has been optimized for sensitivity and specificity. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test. Do not deviate from the specified procedure.
- Microbial contamination of reagents may decrease the accuracy of the assay. Avoid microbial contamination of reagents by using sterile disposable pipettes if removing aliquots from reagent bottles.
- Be attentive to the total assay time when testing more than one fecal specimen. Add *Diluent* first, then add the *Conjugate* to each tube of *Diluent*. Then add specimen to the tube of *Diluent/Conjugate*. Thoroughly mix all of the diluted specimens, and then transfer to the *Membrane Device*. The 15-minute incubation step begins after the last diluted sample-conjugate mixture has been transferred to the final *Membrane Device*.
- If the *Substrate* reagent changes to a dark blue/violet color call technical services for replacement.
- Reagents contain thimerosal as a preservative. Handle reagents according to existing regulations for laboratory safety and good laboratory practice. Safety Data Sheets for this product are available upon request, contact technical support.
- Follow your national, regional, and local ordinances accordingly for waste disposal regulations. Do not place in trash, dispose of as hazardous waste.

#### COLLECTION AND HANDLING OF FECAL SPECIMENS


- Standard collection and handling procedures used in-house for fecal specimens are appropriate. Specimens should be stored between 2° and 8°C; test specimens that are less than 24 hours old, whenever possible.

2. Store specimens frozen ( $\leq -10^{\circ}\text{C}$ ) if the test cannot be performed within 72 hours of collection, but note that freezing and thawing of the specimen may result in loss of activity due to degradation of the toxins.
3. Make sure that specimens are thoroughly mixed PRIOR to performing the assay.
4. Storing fecal specimens in the *Diluent* is NOT recommended.
5. Do not allow the fecal specimens to remain in the *Diluent/Conjugate* for > 24 hours.

#### SPECIMEN PREPARATION

1. Bring all reagents and the required number of devices to room temperature before use.
2. Set up and label one small test tube for each specimen, and optional external controls as necessary.
3. Add 500  $\mu\text{L}$  *Diluent* to each tube for fecal specimens using the graduated black dropper assembly (or equivalent). For specimens in transport media such as Cary Blair or C&S, add 425  $\mu\text{L}$  of *Diluent* to the tube.
4. Add one drop of *Conjugate* (red capped bottle) to each tube.
5. Obtain one disposable plastic transfer pipette (supplied with the kit) for each sample – the pipettes have raised graduations at 25  $\mu\text{L}$ , 400  $\mu\text{L}$  and 500  $\mu\text{L}$ .

#### Graduated Transfer Pipette:

6. Mix all specimens thoroughly regardless of consistency- it is essential that the specimens be evenly suspended before transferring.

**Liquid/Semi-solid specimens** – pipette 25  $\mu\text{L}$  of specimen with a transfer pipette (graduated at 25  $\mu\text{L}$ , 400  $\mu\text{L}$  and 500  $\mu\text{L}$ ) and dispense into the *Diluent/Conjugate* mixture. Use the same transfer pipette to mix the diluted specimen.

**Formed/Solid specimens** – Care must be taken to add the correct amount of formed feces to the sample mixture. Mix the specimen thoroughly using a wooden applicator stick and transfer a small portion (approximately 2 mm diameter, the equivalent of 25  $\mu\text{L}$ ) of the specimen into the *Diluent/Conjugate* mixture. Emulsify the specimen using the applicator stick.

**Fecal specimens in Cary Blair or C&S transport media** - pipette 100  $\mu\text{L}$  of sample into the *Diluent/Conjugate* mixture.

7. **Optional External Control Samples:**

**External Positive Control** - add one drop of *Positive Control* (gray-capped bottle) to the appropriate test tube.

**External Negative Control** - add 25  $\mu\text{L}$  *Diluent* to the appropriate test tube.

**NOTE:** Transferring too little specimen, or failure to mix and completely suspend the specimen in the *Diluent* mixture, may result in a false-negative test result. The addition of too much fecal specimen may cause invalid results due to restricted sample flow.

#### TEST PROCEDURE

1. Obtain one *Membrane Device* per specimen, and one device per optional external positive or negative control as necessary. The foil bags containing the devices should be brought to room temperature before opening. Label each device appropriately and orient it on a flat surface so the letter “C” on the device is on the left, the letter “T” is on the right, and the small *Sample Well* is located in the top right corner of the device (Fig. 1a).
2. Close each tube of diluted specimen and mix thoroughly. Proper mixing can be achieved by vortexing or inverting the tube. Once a patient sample or *Positive Control* has been diluted in the *Diluent/Conjugate* mixture, it may be incubated at room temperature for any period of time up to 24 hours prior to addition to the *Membrane Device*.

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- Using a transfer pipette (graduated at 25 µL, 400 µL and 500 µL), transfer 400 µL of the diluted sample-conjugate mixture into the **Sample Well** (smaller hole in the top right corner of the device) of a *Membrane Device*, making certain to expel the liquid sample onto the wicking pad inside of the *Membrane Device*.
- Incubate the device at room temperature for 15 minutes – the sample will wick through the device and a wet area will spread across the *Reaction Window* (larger hole in the middle of the device).

### NOTE FOR SAMPLES THAT FAIL TO MIGRATE:


*Occasionally, a diluted fecal specimen cannot be tested because it clogs the membrane and the Reaction Window does not wet properly. If the diluted fecal specimen fails to migrate properly within 5 minutes of adding the sample to the Sample Well (i.e. the membrane in the Reaction Window does not appear to be completely wet), then add 100 µL (4 drops) of Diluent to the Sample Well and wait an additional 5 minutes (for a total of 20 minutes).*

- After the incubation, add 300 µL of *Wash Buffer* to the **Reaction Window** using the graduated white dropper assembly (or equivalent). Allow the *Wash Buffer* to flow through the *Reaction Window* membrane and be absorbed completely.
- Add 2 drops of *Substrate* (white-capped bottle) to the **Reaction Window**. Read and record results visually after 10 minutes.

### INTERPRETATION OF RESULTS

- Interpretation of the test is most reliable when the device is read immediately at the end of the reaction period. Read the device at a normal working distance in a well-lit area. View with a line of vision directly over the device.
- Observe device for the appearance of a blue line on the "C" side of the *Reaction Window* representing the internal positive control line. Observe device for the appearance of a blue line on the "T" side of the *Reaction Window* representing the test line. The lines may appear faint to dark in intensity.
- Positive Result:** A positive result may be interpreted at any time between the addition of *Substrate* and the 10-minute read time. Two blue lines are visible, the control line ("C") and the test line ("T"). The lines may appear faint to dark in intensity. The appearance of a blue line on the "T" side along with a blue control line is interpreted as a positive result. An obvious partial line is interpreted as a positive result. Do not interpret membrane discoloration or shadow as a positive result. A positive result indicates the presence of *C. difficile* toxin.
- Negative Result:** A test cannot be interpreted as negative or invalid until 10 minutes following the addition of *Substrate*. A single blue line is visible on the control ("C") side of the *Reaction Window* and no test line is visible on the "T" side of the *Reaction Window* (Fig. 1b). A negative result indicates *C. difficile* toxin is either absent in the specimen or is below the detection limit of the test.
- Invalid Result:** A single line is visible on the test ("T") side of the *Reaction Window*, or no lines are visible in the *Reaction Window* (Fig. 1c, 1d). The test result is invalid if a control line is not present at the completion of the reaction period.

**FIGURE 1: TOX A/B QUIK CHEK® INTERPRETATION OF RESULTS**

## QUALITY CONTROL

**Internal:** A blue control line must be visible on the “C” side of the *Reaction Window* on every *Membrane Device* that is tested. The appearance of the blue control line confirms that the sample and reagents were added correctly, that the reagents were active at the time of performing the assay, and that the sample migrated properly through the *Membrane Device*. A clear background in the result area is considered an internal negative control. If the test has been performed correctly and reagents are working properly, the background will be clear to give a discernible result.

**External:** The reactivity of the *TOX A/B QUIK CHEK*® test should be verified on receipt using the *Positive Control* and negative control (*Diluent*). The *Positive Control* is supplied with the kit (gray-capped bottle). The *Positive Control* confirms the reactivity of the other reagents associated with the assay, and is not intended to ensure precision at the analytical assay cut-off. *Diluent* is used for the negative control. Additional tests can be performed with the controls to meet the requirements of local, state and/or federal regulations and/or accrediting organizations.

## LIMITATIONS

1. The *TOX A/B QUIK CHEK*® test is used to detect *C. difficile* toxin(s) in fecal specimens. The test confirms the presence of toxin in feces and this information should be taken under consideration by the physician in light of the clinical history and physical examination of the patient. The *TOX A/B QUIK CHEK*® test will detect levels of toxin A at  $\geq 0.63$  ng/mL and toxin B at  $\geq 1.25$  ng/mL.
2. Fecal specimens are extremely complex. Optimal results with the *TOX A/B QUIK CHEK*® test are obtained with specimens that are less than 24 hours old. Most undiluted specimens can be stored between 2° and 8°C for 72 hours before significant degradation of the toxin is noted. If specimens are not assayed within this time period, they may be frozen and thawed. However, repeated freezing and thawing may result in loss in the immunoreactivity of toxins A and B.
3. Some specimens may give weak reactions. This may be due to a number of factors such as the presence of low levels of toxin, the presence of binding substances, or inactivating enzymes in the feces. *Under these conditions, a fresh specimen should be tested.* Additional tests that may be used in conjunction with the *TOX A/B QUIK CHEK*® test include culture with toxigenic testing or tissue culture cytotoxicity assay for the detection of *C. difficile* or its toxin(s).
4. Fecal specimens preserved in 10% Formalin, merthiolate formalin, sodium acetate formalin, or polyvinyl alcohol cannot be used.
5. The *TOX A/B QUIK CHEK*® test is qualitative. The intensity of the color should not be interpreted quantitatively.
6. Some isolates of *C. sordellii* may react in the *TOX A/B QUIK CHEK*® test due to the production of immunologically related toxins (1).
7. Colonization rates of up to 50% have been reported in infants. A high rate has also been reported in cystic fibrosis patients (1,3).

## EXPECTED VALUES

The reported incidence of *C. difficile*-associated disease in patients with antibiotic-associated diarrhea is 10 to 20%. In our studies, the incidence ranged from 10% to 22%. The prevalence of a positive *TOX A/B QUIK CHEK*® test will vary from location to location and hospitals may experience rates lower or higher than those observed at the sites used in this evaluation. *Clostridium difficile* disease is primarily a nosocomial disease of elderly patients, and hospitals that have higher numbers of elderly patients may experience higher rates. It is important to consider any test results in conjunction with clinical symptoms because some healthy adults and large numbers of healthy infants (up to 50%) will be positive for *C. difficile* toxin.

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#### PERFORMANCE CHARACTERISTICS

The *TOX A/B QUIK CHEK*® test was compared with the tissue culture test at three U.S. hospitals and in-house at TECHLAB, Inc. Specimens included in the evaluation were submitted to the clinical laboratory for routine testing. The tissue culture test was done according to the in-house procedure. The table below shows a summary of the clinical performance of the *TOX A/B QUIK CHEK*® test. The test exhibited a sensitivity and specificity of 90.2% and 99.7%, respectively. The predictive positive and negative values were 98.6% and 97.9%, respectively, and the correlation was 98.0%.

TABLE 1. Correlation of the *TOX A/B QUIK CHEK*® test with tissue culture.

N = 842	Tissue Culture positive	Tissue Culture negative
<i>TOX A/B QUIK CHEK</i> ® positive	138	2
<i>TOX A/B QUIK CHEK</i> ® negative	15	687

		95% CI
<b>Sensitivity</b>	90.2%	84.1 - 94.2
<b>Specificity</b>	99.7%	98.8 - 99.9
<b>Predictive Positive Value</b>	98.6%	94.4 - 99.8
<b>Predictive Negative Value</b>	97.9%	96.4 - 98.7
<b>Correlation</b>	98.0%	97.8 - 98.2

Of the 2 tissue culture-negative/*TOX A/B QUIK CHEK*®-positive samples, 1 was negative in a commercial toxin A+B ELISA. Of the 15 specimens that were tissue culture-positive/*TOX A/B QUIK CHEK*®-negative, 12 were negative in commercial toxin A+B ELISAs. There were 9 specimens that were unreadable. All of the specimens were negative by PCR analysis for the genes of toxin A (*tcdA*) and toxin B (*tcdB*).

A total of 51 fecal specimens diluted in Cary Blair and 32 fecal specimens diluted in C&S Transport Media were tested in the *TOX A/B QUIK CHEK*® test and the results were compared to those obtained by routine testing. The test exhibited an agreement of 97.6% for the detection of *C. difficile* toxins in specimens prepared in Transport Media.

#### ANALYTICAL SENSITIVITY

The test was consistently positive at a concentration of 0.63 ng/mL for toxin A and 1.25 ng/mL for toxin B.

#### REPRODUCIBILITY

The reproducibility of the *TOX A/B QUIK CHEK*® test was determined using known positive (n=6) and negative (n=2) fecal specimens that were coded and sorted to prevent their identification during testing. Testing was performed on-site at 3 laboratories, which tested the samples on 3 days. The samples produced the expected results 100% of the time.

#### CROSS REACTIVITY

Fecal specimens inoculated with the following microorganisms to a final concentration of approximately 10<sup>8</sup> or higher organisms per mL did not react in the *TOX A/B QUIK CHEK*®:

## TOX A/B QUIK CHEK® - FRANÇAIS

**Bactéries:** *Aeromonas hydrophila*, *Bacillus cereus*, *Bacillus subtilis*, *Bacteroides fragilis*, *Campylobacter coli*, *Campylobacter fetus*, *Campylobacter jejuni*, *Candida albicans*, *Clostridium bifermentans*, *Clostridium butyricum*, *Clostridium perfringens*, *Clostridium septicum*, *Clostridium sordellii* (nontoxigenic), *Clostridium sporogenes*, *Enterococcus faecalis*, *Escherichia coli* EIEC, *Escherichia coli*, *Escherichia coli* 0157 H7, *Escherichia coli* ETEC, *Klebsiella pneumoniae*, *Peptostreptococcus anaerobius*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Salmonella typhimurium*, *Shigella dysenteriae*, *Shigella flexneri*, *Shigella sonnei*, *Staphylococcus aureus*, *Staphylococcus aureus* (Cowans), *Staphylococcus epidermidis*, *Vibrio parahaemolyticus*, *Yersinia enterocolitica*

**Viruses:** Adenovirus types 1,2,3,5,40,41, Human coronavirus, Coxsackievirus B2,B3,B4,B5, Echovirus 9,11,18,22,33, Enterovirus type 68,69,70,71.

The only non-*C. difficile* organism to react with the TOX A/B QUIK CHEK® was *C. sordellii* VPI 9048. This strain produces toxins HT and LT, which are homologous to toxins A and B, respectively.

### INTERFERING SUBSTANCES

The following substances had no effect on test results when present in feces in the concentrations indicated: mucin (3.5% w/v), human blood (40% v/v), barium sulfate (5% w/v), Imodium® (5% w/v), Kaopectate® (5 mg/mL), Pepto-Bismol® (5% w/v), steric/palmitic acid (40% w/v), Metronidazole (0.25% w/v), Vancomycin (0.25% w/v).

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