FOR INFORMATIONAL USE
ONLY

TECHLAB® *E. HISTOLYTICA* QUIK CHEK™

TECHLAB® E. HISTOLYTICA QUIK CHEK™

ENGLISH p. 3

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TECHLAB® *E. HISTOLYTICA* QUIK CHEK™

INTENDED USE

The TECHLAB® *E. HISTOLYTICA* QUIK CHEK™ test is a rapid membrane enzyme immunoassay for the qualitative detection of adhesin from *Entamoeba histolytica* in a single use cassette. It is intended for use with human fecal specimens from patients with diarrhea or dysentery as an aid in the diagnosis of *E. histolytica* gastrointestinal infection. Test results should be considered in conjunction with patient history.

FOR IN VITRO DIAGNOSTIC USE

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician

EXPLANATION

Entamoeba histolytica and *Entamoeba dispar* are intestinal parasites that infect approximately half a billion people worldwide each year (1). It is necessary to distinguish between the two species because *E. histolytica* is pathogenic, causing intestinal amebiasis (e.g. diarrhea, dysentery, colitis) and extra-intestinal amebiasis (e.g. liver abscess). *E. dispar* is not associated with symptomatic disease and inaccurate diagnosis may result in unnecessary treatment. The most common method used to diagnose amebiasis has been wet mount microscopy, which suffers from poor sensitivity and specificity.

Trophozoites and cysts are not easily identified in a single fecal specimen and it is difficult to visually distinguish between *E. dispar* and *E. histolytica* when they are observed. Detection of *Entamoeba* spp. by immunoassay provides an alternative method of diagnosis with greater sensitivity (2). Immunoassays specific for *E. histolytica*, such as the *E. HISTOLYTICA* QUIK CHEK™ test, provide the added benefit of only identifying *E. histolytica* infections. Large numbers of specimens can be tested rapidly and objectively, and the procedure is less labor-intensive than most other methods of diagnosis.

PRINCIPLE OF THE TEST

The *E. HISTOLYTICA* QUIK CHEK™ test utilizes antibodies specific for *E. histolytica*. The *Membrane Device* contains a *Reaction Window* with two vertical lines of immobilized antibodies. The test line ("T") contains monoclonal antibodies specific for *E. histolytica* adhesin. The control line ("C") contains antibodies to horseradish peroxidase (HRP). The *Conjugate* consists of antibodies to *E. histolytica* coupled to horseradish peroxidase. To perform the test, the sample is added to a tube containing a mixture of *Diluent* and *Conjugate*. The diluted sample-conjugate mixture is added to the *Sample Well* and the device is allowed to incubate at room temperature for 15 minutes. During the incubation, any *E. histolytica* adhesin in the sample binds to the antibody-peroxidase conjugate. The antigen-antibody-peroxidase complexes migrate through a filter pad to a membrane where they are captured by the immobilized anti-adhesin antibodies in the line. The *Reaction Window* is subsequently washed with *Wash Buffer*, followed by the addition of *Substrate*. After a 10 minute incubation period, the *Reaction Window* is examined visually for the appearance of vertical blue lines on the "C" and "T" sides of the *Reaction Window*. A blue line on the "T" side of the *Reaction Window* indicates a positive result. A positive "C" reaction, indicated by a vertical blue line on the "C" side of the *Reaction Window*, monitors/ confirms that the sample and reagents were added correctly, the reagents were active at the time of performing the assay, and that the sample migrated properly through the *Membrane Device*. It also confirms the reactivity of the other reagents associated with the assay.

MATERIALS PROVIDED

MEM DEV

Membrane Devices – Each pouch contains 1 device

CONJ ENZ

Conjugate (2 mL) – Antibody specific for *E. histolytica* coupled to horseradish peroxidase in a buffered protein solution (contains 0.05% ProClin® 300)*

DIL SPE

Diluent (16 mL) – Buffered protein solution with grey graduated dropper assembly (contains 0.05% ProClin® 300)*

CONTROL +

Positive Control (1 mL) – *E. histolytica* antigen in a buffered protein solution (contains 0.05% ProClin® 300)*

SUBS REAG

Substrate (3.5 mL) – Solution containing tetramethylbenzidine

WASH REAG

Wash Buffer (12 mL) – Buffered solution with white graduated dropper assembly (contains 0.05% ProClin® 300)*

*(contains 0.05% ProClin® 300)

Signal Word: Warning

H317: May cause an allergic skin reaction

P261, P272, P280, P302, P352, P333, P313, P321, P362, P364, P501

Disposable plastic pipettes (50) – Graduated at 25 µL, 100 µL, 200 µL, 300 µL, 400 µL, and 500 µL

MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- Small test tubes (e.g., plastic microcentrifuge tubes)
- Disposable gloves
- Pipettor and tips
- Wooden Applicator Sticks
- Timer

SHELF LIFE AND STORAGE

The expiration date of the kit is printed on the kit box label. Store the kit between 2°C and 8°C. Return the kit to the refrigerator as soon as possible after use.

PRECAUTIONS

1. Rx Only – Prescription Only
2. Upon arrival, inspect the kit to ensure that components are not frozen or warm to the touch due to improper shipping conditions, and that there are no signs of leakage.
3. Bring all components to room temperature before use to ensure proper kit reactivity.
4. The *Substrate* reagent should be colorless. If the *Substrate* reagent changes to a dark blue/violet color, discard and call Technical Services for a replacement.
5. Reagents from different kits should not be mixed or interchanged. Do not use a kit past the expiration date.
6. Caps, tips and dropper assemblies are color-coded; do NOT mix or interchange!
7. Use fecal specimens according to recommendations in the table below to obtain optimal results. Specimens that are frozen may lose reactivity due to freezing and thawing. Raw fecal specimens that are stored frozen may be thawed up to 5 times. Fecal specimens in transport media that are stored frozen may be thawed one time. When storing specimens, avoid temperature extremes and keep out of direct sunlight.

- The test has been optimized for sensitivity and specificity. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test. Do not deviate from the specified procedure.
- Fecal specimens and used membrane devices may contain potentially infectious agents and should be handled at "Biosafety Level 2" as recommended in the CDC/NIH Manual "Biosafety in Microbiological and Biomedical Laboratories."
- Wear disposable gloves when doing the test.
- The *Conjugate*, *Diluent*, *Positive Control* and *Wash Buffer* reagents contain 0.05% ProClin® 300 as a preservative. Although the concentration is low, ProClin® 300 is known to be harmful. If skin irritation or rash occurs, get medical advice/attention. Take off contaminated clothing and wash it before reuse. Handle reagents according to existing regulations for laboratory safety and good laboratory practice. Safety Data Sheets for this product are available upon request, contact technical support.
- Follow your national, regional, and local ordinances accordingly for waste disposal regulations.

COLLECTION, HANDLING AND STORAGE OF FECAL SPECIMENS

Acceptable Sample Types
Fresh Fecal Specimens
Frozen Fecal Specimens (frozen undiluted)
Fecal specimens in transport media (e.g. Cary Blair, C&S)

Do Not Use
Fecal specimens in Formalin-based fixative (e.g. Sodium Acetate Formalin, 10% formalin)
Fecal specimens in alcohol-based fixative (e.g. Polyvinyl Alcohol)


Sample Storage Temperature	Acceptable Length of Storage	Comments
Room Temperature (18°C - 25°C)	24 hours	Fresh samples that will be tested within 24 hours may remain at room temperature (18°C - 25°C). If testing is not planned in <24 hours, refrigerate (2°C - 8°C) as soon as possible following collection.
Refrigerated (2°C - 8°C)	1 week	
Frozen ≤ -10°C	7 months	Freeze specimens and store at ≤ -10°C if testing cannot be performed within 1 week of collection. Thaw at room temperature. Freezing and thawing multiple times may result in loss of specimen activity due to antigen degradation.

- Use standard in-house collection and handling procedures for fecal specimens. Collect fecal specimens in clean, leak-proof containers.
- Do not store fecal specimens in the *Diluent*.

TEST PROCEDURE


- Be attentive to the total assay time when testing more than one fecal specimen.
- Bring all reagents and devices to room temperature before use. Remove the reagents from the foam insert to reduce the time needed to warm to room temperature.
- Set up and label one small test tube for each specimen and optional external control.
- Using the gray graduated dropper assembly, add 500 µL of *Diluent* (2nd graduation from the tip) to each tube for fresh and frozen specimens, and external controls. For specimens in transport media such as Cary Blair or C&S, add 400 µL (1st graduation from the tip) of *Diluent* to the tube.

Sample Type	Volume of Diluent
Fresh Fecal Specimens	500 µL (2 nd graduation from tip)
Frozen Fecal Specimens (frozen undiluted)	500 µL (2 nd graduation from tip)
Specimens in transport media (Cary Blair, C&S)	400 µL (1 st graduation from tip)
External Controls (positive and negative)	500 µL (2 nd graduation from tip)

- Add one drop of *Conjugate* (red capped bottle) to each tube. Hold the dropper bottle vertically to ensure proper drop size. The *Diluent* and *Conjugate* should be added to all tubes prior to adding the specimens.
- Obtain one disposable plastic transfer pipette (supplied with the kit) for each sample.

Graduated Transfer Pipette:

- For Liquid/Semi-Solid Specimens** - Mix specimen thoroughly. Using a transfer pipette, add 25 µL of specimen to the *Diluent/Conjugate* mixture in the tube.

For Formed/Solid specimens - Mix specimen thoroughly using a wooden applicator stick and transfer a small portion (approximately 2 mm diameter, the equivalent of 25 µL) of the specimen into the *Diluent/Conjugate* mixture. Emulsify the specimen using the applicator stick.

Fecal specimens in Cary Blair or C&S transport media - pipette 100 µL (2 drops from transfer pipette) of sample into the *Diluent/Conjugate* mixture.

NOTE: Transferring too little sample, or failure to mix and completely suspend the sample in the *Diluent/Conjugate* mixture, may result in a false-negative test result. The addition of too much sample may cause invalid results due to restricted flow.

8. **Optional External Controls:**


Option control cassettes may be run concurrently with patient samples.

External Positive Control - add one drop of *Positive Control* (gray-capped bottle) to the appropriate test tube.

External Negative Control - add 25 µL *Diluent* to the appropriate test tube.

9. For all test and control samples, close the tubes and mix thoroughly using a vortex mixer or by inverting the tube several times. Samples or controls diluted in the *Diluent/Conjugate* mixture may be incubated at room temperature up to 2 hours prior to addition to the *Membrane Device*.
10. Open one room temperature *Membrane Device* pouch for each diluted specimen and external control (as necessary). Label each device appropriately and orient it on a flat surface so the "E HISTO QUIK CHEK" print is at the bottom of the device, and the small *Sample Well* is located in the top right corner of the device.

Membrane Device

11. Make sure that each diluted sample is thoroughly mixed (See Step 9) before adding to the *Membrane Device*. **Using a new transfer pipette, transfer 500 µL (topmost graduation) from each tube into the *Sample Well* (smaller hole in the top right corner of the device) of a *Membrane Device*.** When adding the sample into the *Sample Well*, make sure that the tip of the transfer pipette is inside the *Sample Well* hole and angled towards the *Reaction Window* and onto the wicking pad.
12. **Incubate the device at room temperature for 15 minutes** – the sample will wick through the device and a wet area will spread across the *Reaction Window*. The 15-minute incubation step begins after the last diluted sample-conjugate mixture has been transferred to the last *Membrane Device*.
NOTE FOR SAMPLES THAT FAIL TO MIGRATE: Occasionally, a diluted sample fails to migrate properly and the *Reaction Window* does not fully wet. If the *Reaction Window* does not appear to be completely wet within 5 minutes of adding the sample to the *Sample Well*, then add 100 µL (4 drops) of *Diluent* to the *Sample Well* and wait an additional 5 minutes (for a total of 20 minutes). Continue with the next step of the Test Procedure.
13. **After the incubation, add 300 µL of *Wash Buffer* to the central *Reaction Window* using the graduated white dropper assembly.** Allow the *Wash Buffer* to be absorbed completely.
14. **Add 2 drops of *Substrate* (white-capped bottle) to the central *Reaction Window*.**
15. Incubate 10 minutes at room temperature. Read visually and record results after 10 minutes.

INTERPRETATION OF RESULTS

Positive Result

Negative Result

Invalid Result

Invalid Result

1. Interpretation of the test is most reliable when the device is read at the end of the 10 minute reaction period. Read the device at a normal working distance in a well-lit area. View with a line of vision directly over the device.
2. Observe device for the appearance of a blue line on the "C" side of the *Reaction Window* representing the internal positive control line. Observe device for the appearance of a blue line on the "T" side of the *Reaction Window* representing the test line. The lines may appear faint to dark in intensity.
3. **Positive Result:** A positive result may be interpreted at any time between the addition of *Substrate* and the 10-minute read time. Two blue lines are visible, the control line ("C") and the test line ("T"). The lines may appear faint to dark in intensity. The appearance of a blue line on the "T" side along with a blue control line is interpreted as a positive result. An obvious partial line is interpreted as a positive result. Do not interpret membrane discoloration or shadow as a positive result. A positive result indicates the presence of *E. histolytica*.
4. **Negative Result:** A test cannot be interpreted as negative or invalid until 10 minutes following the addition of *Substrate*. A single blue line is visible on the control ("C") side of the *Reaction Window* and no test line is visible on the "T" side of the *Reaction Window*. A negative result indicates *E. histolytica* is either absent in the specimen or is below the detection limit of the test.
5. **Invalid Result:** A single line is visible on the test ("T") side of the *Reaction Window*, or no lines are visible in the *Reaction Window*. The test result is invalid if a control line is not present at the completion of the reaction period.

QUALITY CONTROL

The validity of the test results using the *E. HISTOLYTICA QUIK CHEK™* test is dependent upon the proper reaction of the internal and external controls.

Internal: A vertical blue control line must be visible on the "C" side of the *Reaction Window* on every *Membrane Device* that is tested. This confirms that the sample and reagents were added correctly and reacted properly in the assay. A clear background in the result area is considered an internal negative control. It may appear white to light blue and any developed lines will be clearly visible.

External: The reactivity of the *E. HISTOLYTICA QUIK CHEK™* test should be verified on receipt using the *Positive Control* and negative control (*Diluent*). The *Positive Control* confirms the reactivity of the other reagents associated with the assay, and is not intended to ensure precision at the analytical assay cut-off. Additional tests can be performed with the controls to meet the requirements of local, state and/or federal regulations and/or accrediting organizations.

LIMITATIONS

1. A negative test result does not preclude the possibility of the presence of *E. histolytica* adhesin in the specimen which may occur if the level of antigen is below the detection limit of the test.
2. The *E. HISTOLYTICA QUIK CHEK™* test is qualitative. The intensity of the color should not be interpreted quantitatively.
3. Due to the small number of positive specimens collected during the prospective clinical study, performance characteristics for *E. histolytica* were also established with retrospective clinical specimens.
4. Transferring too little sample, or failure to mix and completely suspend the sample in the *Diluent/Conjugate* mixture, may result in a false-negative test result. The addition of too much sample may cause invalid results due to restricted flow.

EXPECTED VALUES

Normal healthy individuals should not be infected with *E. histolytica* and should test negative in the *E. HISTOLYTICA QUIK CHEK™* test. A positive test result in the *E. HISTOLYTICA QUIK CHEK™* test indicates that the person is shedding detectable amounts of *E. histolytica* antigen. The incidence of *E. histolytica* infection varies significantly between populations and geographic regions. It is estimated that *Entamoeba histolytica* infects about 50 million people around the world (2). Roughly 90% of these persons remain asymptomatic, whereas about 10% develop clinical symptoms ranging from gastrointestinal disease to liver abscesses. High risk groups include persons who have traveled abroad, immigrants, immunocompromised persons, migrant workers, and active male homosexuals (2, 3). Nonpathogenic strains (*E. dispar*) are predominant among male homosexuals (4). The disease often is transmitted by asymptomatic carriers of *E. histolytica*.

PERFORMANCE CHARACTERISTICS

The performance of the *E. HISTOLYTICA QUIK CHEK™* test was compared to a Composite Reference Method (CRM), which included molecular detection of *Entamoeba histolytica*. A total of 851 fecal samples were evaluated and included 96 retrospective samples. Age information was available for 851 patients. Of the 851 patients, 18.9% were ≤ 20 years. Sex information was available for 851 patients, and 42.7% were male and 57.3% were female. Table 1 and 2 show a summary of the clinical performance of the *E. HISTOLYTICA QUIK CHEK™* test. Table 1 represents the results for the prospective samples, of the 755 prospectively collected samples, 100 of the prospectively collected samples were diluted in Protocol™ Cary-Blair and Para-Pak® C&S. All of the samples diluted in transport media and tested were negative. The prospective testing did show that the *E. HISTOLYTICA QUIK CHEK™* test exhibited a sensitivity of 40.0%, with a specificity of 100%, a predictive positive value of 100%, and a predictive negative value of 99.6% with the CRM. Table 2 represents the results for the retrospective samples, which showed that the *E. HISTOLYTICA QUIK CHEK™* test exhibited a sensitivity and specificity of 100% with the CRM.

Table 1. Summary of prospective clinical performance comparing the *E. HISTOLYTICA QUIK CHEK™* test to the Composite Reference Method (CRM) Prospective Samples

N = 755	Composite Reference Method Positive	Composite Reference Method Negative
<i>E. HISTOLYTICA QUIK CHEK™</i> Positive	2	0
<i>E. HISTOLYTICA QUIK CHEK™</i> Negative	3	750

95% Confidence Limits		
Sensitivity	40.0%	7.3% - 83.0%
Specificity	100%	99.4% - 100%
Predictive Positive Value	100%	19.8% - 100%
Predictive Negative Value	99.6%	98.7% - 99.9%

All three of the false negative results were PCR positive and antigen negative. Additional antigen testing was done with a previously cleared FDA device.

Table 2. Summary of retrospective clinical performance comparing the *E. HISTOLYTICA QUIK CHEK™* test to the Composite Reference Method (CRM) Retrospective Samples

N = 96	Composite Reference Method Positive	Composite Reference Method Negative
<i>E. HISTOLYTICA QUIK CHEK™</i> Positive	30	0
<i>E. HISTOLYTICA QUIK CHEK™</i> Negative	0	66

95% Confidence Limits		
Sensitivity	100%	85.9% - 100%
Specificity	100%	93.1% - 100%

REPRODUCIBILITY

The reproducibility of the *E. HISTOLYTICA QUIK CHEK™* test was determined using 8 fecal specimens that were coded to prevent their identification during testing. Testing was performed at 2 independent laboratories and on-site at TECHLAB, Inc. The samples included 2 negative samples, 2 high negative samples, 2 low positive samples and 2 moderate positive samples. The samples were tested in triplicate twice a day over a 5-day period by multiple technicians at each site using 2 different kit lots. A positive and negative control was run with each panel of the masked samples. The results from each laboratory were submitted to TECHLAB, Inc. and compared with in-house results. The results were consistent among the different locations, and exhibited a correlation of 100%. The samples produced the expected results 100% of the time.

CROSS REACTIVITY

The *E. HISTOLYTICA QUIK CHEK™* test was evaluated for cross-reactivity with the bacterial and viral strains listed below. None of the strains were shown to interfere with the performance *E. HISTOLYTICA QUIK CHEK™* test.

<i>Aeromonas hydrophila</i>	<i>Clostridium difficile</i>	<i>Shigella dysenteriae</i>
<i>Bacillus cereus</i>	<i>Enterococcus faecalis</i>	<i>Shigella flexneri</i>
<i>Bacillus subtilis</i>	<i>Escherichia coli</i>	<i>Shigella sonnei</i>
<i>Bacteroides fragilis</i>	<i>Escherichia coli</i> O157:H7	<i>Staphylococcus aureus</i>
<i>Campylobacter coli</i>	<i>Escherichia coli</i> EIEC	<i>Staphylococcus aureus</i> (Cowan's)
<i>Campylobacter fetus</i>	<i>Escherichia coli</i> EPEC	<i>Staphylococcus epidermidis</i>
<i>Campylobacter jejuni</i>	<i>Escherichia coli</i> ETEC	<i>Vibrio parahaemolyticus</i>
<i>Candida albicans</i>	<i>Klebsiella pneumoniae</i>	<i>Yersinia enterocolitica</i>
<i>Clostridium bifermentans</i>	<i>Salmonella typhimurium</i>	
Adenovirus, 2, 5, 40, 41	Human Adenovirus 1, 3	Human Enterovirus 69, 70, 71
Coxsackievirus B5	Human Coronavirus	Human parechovirus 1
Echovirus 11, 18, 22, 33	Human Coxsackievirus B2, B3, B4	[Echovirus 22]
Enterovirus 68, 69	Human Echovirus 9	Human Rotavirus

Cross reactivity with Norovirus is unknown because it was not tested in analytical studies. However, Norovirus GI/GII was identified in 50 clinical specimens using an FDA cleared multiplex NAAT assay during clinical testing and no cross reactivity was found using the *E. HISTOLYTICA QUIK CHEK™* in those samples.

Additionally, the *E. HISTOLYTICA QUIK CHEK™* test was run on fecal specimens documented to be positive for other parasites by microscopy. The number in parenthesis is the number of clinical specimens in which each organism was identified. No cross-reactivity was seen with the following organisms.

<i>Ascaris lumbricoides</i> and with eggs (21)	<i>Entamoeba bangladeshi</i> (3)	<i>Giardia spp.</i> (45)
<i>Blastocystis hominis</i> (12)	<i>Entamoeba coli</i> (13)	<i>Iodamoeba bütschlii</i> (10)
<i>Cryptosporidium spp.</i> (30)	<i>Entamoeba moshkovskii</i> (3)	<i>Trichuris trichiura</i> eggs (11)

STRAIN SPECIFIC STUDY

The specificity of the *E. HISTOLYTICA QUIK CHEK™* test was also evaluated by examining the reactivity of pathogenic (*Entamoeba histolytica*) and non-pathogenic (*Entamoeba dispar*) zymodemes (strains) for reactivity by standard curve dilutions. *E. histolytica* results were positive from 244 to 30.5 pathogenic zymodemes (PZs)/mL and the *E. dispar* results were negative at all dilutions beginning at 2440 non-pathogenic zymodemes (NPZs)/mL. The *E. HISTOLYTICA QUIK CHEK™* test demonstrates proper reactivity with *Entamoeba histolytica* and does not cross-react with *Entamoeba dispar*.

Additionally, due to the similarity in morphology, 3 specimens identified by PCR as positive for *Entamoeba moshkovskii* and 3 positive for *Entamoeba bangladeshi* were evaluated using the *E. HISTOLYTICA QUIK CHEK™* test. These 6 specimens all tested negative in the *E. HISTOLYTICA QUIK CHEK™* test.

INTERFERING SUBSTANCES (U.S. FORMULATION)

The following substances had no effect on positive or negative test results analyzed at the concentrations indicated: Barium sulfate (5% w/v), Benzalkonium Chloride (1% w/v), Ciprofloxacin (0.25% w/v), Ethanol (1% w/v), Hog gastric mucin (3.5% w/v), Human blood (40% v/v), Hydrocortisone (1% w/v), Imodium® (5% v/v), Kaopectate® (5% v/v), Leukocytes (0.05% w/v), Maalox® Advanced (5% v/v), Mesalazine (10% w/v), Metronidazole (0.25% w/v), Mineral Oil (10% w/v), Mylanta® (4.2 mg/mL), Naproxen Sodium (5% w/v), Nonoxonyl-9 (40% w/v), Nystatin (1% w/v), Palmitic Acid/Fecal Fat (40% w/v), Pepto-Bismol® (5% v/v), Phenylephrine (1% w/v), Polyethylene glycol 3350 (10% w/v), Prilosec OTC® (5 µg/mL), Sennosides (1% w/v), Simethicone (10% w/v), Steric Acid/Fecal Fat (40% w/v), Tagamet® (5 µg/mL), TUMS (50 µg/mL), Human Urine (5% v/v), and Vancomycin (0.25% w/v).

PRECISION – INTRA-ASSAY

For the determination of intra-assay performance, twelve human fecal specimens were analyzed by the *E. HISTOLYTICA QUIK CHEK™* test. Of these twelve specimens, six were positive for *E. histolytica* of varying levels (low, moderate, and high) and six were negative for *E. histolytica*. Each specimen was assayed five times in the same test run, using two different kit lots. A positive and negative control was run with each panel. All positive samples remained positive and all negative samples remained negative. The overall correlation between the results was 100%.

PRECISION – INTER-ASSAY

For the determination of inter-assay performance, eight human fecal specimens were analyzed by the *E. HISTOLYTICA QUIK CHEK™* test. The samples included 2 negative samples, 2 high negative samples, 2 low positive samples and 2 moderate positive samples. The samples were tested, twice a day by multiple technicians over a 12-day period using 2 different kit lots. A positive and negative control was run on each day. All positive samples remained positive and all negative samples remained negative.

ANALYTICAL SENSITIVITY

The Limit of Detection (LoD) for the *E. HISTOLYTICA QUIK CHEK™* test was established at 320 pathogenic zymodemes (PZs)/mL for *E. histolytica* (equivalent to 15 PZs detected per test). For specimens in Protocol™ Cary Blair media, the LoD was established at 275 PZs/mL for *E. histolytica* (equivalent to 14 PZs detected per test). For specimens in Para-Pak® C&S media, the LoD was established at 245 PZs/mL for *E. histolytica* (equivalent to 12 PZs detected per test).

FRESH VERSUS FROZEN SAMPLES

The effect of long term frozen specimen storage on antigen stability was evaluated. For the analysis, a total of 15 fecal specimens were tested with the *E. HISTOLYTICA QUIK CHEK™* test. The fecal specimens consisted of 3 negative fecal samples, 3 *E. histolytica* high negative fecal samples, 3 *E. histolytica* low positive fecal samples, 3 *E. histolytica* moderate positive fecal samples, and 3 *E. histolytica* high positive fecal samples. Samples were prepared and stored < -10°C at 0, 1, and 4, 8, 12, 16, 20, 24 and 28 weeks. No conversion of positive-to-negative or negative-to-positive was observed in any of the samples at the specified time points.

PROZONE

To ensure that a high concentration of *E. histolytica* antigen does not interfere with a positive reaction in the *E. HISTOLYTICA QUIK CHEK™* test, high samples were prepared by spiking a negative fecal pool at a concentration possibly observed in clinical specimens. A total of 5 different dilutions of the antigen, up to and including the clinically observed high concentration, were prepared and tested in triplicate. The results demonstrated that there was no overall prozone affect, that elevated levels of antigen did not affect the detection of the antigen.

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TECHNICAL SUPPORT

Advice Line

Further information can be obtained from your distributor, or by contacting Alere Technical Support on:

US +	+1 800 257 9525	TS.SCR@alere.com
Africa, Russia, CIS	+972 8 9429 683	ARCISproductsupport@alere.com
Asia Pacific	+61 7 3363 7711	APproductsupport@alere.com
Canada	+1 800 818 8335	CANproductsupport@alere.com
Europe & Middle East	+44 161 483 9032	EMEproductsupport@alere.com
Latin America	+57 2 6618797	LAPproductsupport@alere.com

REFERENCES / BIBLIOGRAFÍA / QUELLENANGABE / RÉFÉRENCES / LITERATURA / REFERENCER / ΒΙΒΛΙΟΓΡΑΦΙΑ / REFERENCIÁK / BIBLIOGRAFIA / REFERENTIES / HENVISNINGER / REFERENSER / REFERANSLAR / REFERINTE / REFERÊNCIAS

1. Haque, R., K. Kress, S. Wood, T. Jackson, D. Lyerly, T. Wilkins, and W. A. Petri, Jr. 1993. Diagnosis of pathogenic *Entamoeba histolytica* infection using a stool ELISA based on monoclonal antibodies to the galactose-specific adhesin. J. Infect. Dis. 167:247-249.
2. Tanyuksel, M. and W.A. Petri, Jr. 2003. Laboratory Diagnosis of Amebiasis. Clin. Microbiol. Rev. 16(4):713-729.
3. Krogstad, D. J., H. C. Spencer, G. R. Healy, N. N. Gleason, D. J. Sexton, and C. A. Herron. 1978. Amebiasis: epidemiologic studies in the United States, 1971-1974. Ann. Intern. Med. 88:89-97.
4. Tannich, E., and G. D. Burchard. 1991. Differentiation of pathogenic from nonpathogenic *Entamoeba histolytica* by restriction fragment analysis of a single gene amplified *in vitro*. J. Clin. Microbiol. 29:250-255.

FOR INFORMATIONAL USE
ONLY

Developed and Manufactured by:

 **TECHLAB, Inc.**
2001 Kraft Drive
Blacksburg, VA 24060-6358, USA

Distributed by:

Alere North America, LLC
30 South Keller Road
Orlando, Florida 32810
TEL 1-877-441-7440
1-321-441-7200 OUTSIDE USA

Made in the USA

© 2017 TECHLAB, Inc. All rights reserved.

The Alere Logo and Alere are trademarks of the Alere group of companies.

E. HISTOLYTICA QUIK CHEK, the TECHLAB Logo and TECHLAB are trademarks of TECHLAB, Inc.

ProClin is a trademark of Rohm and Hass Company.

All trademarks referenced are trademarks of their respective owners.

RMS #91-409-01 Issued: 2017/10

EC REP Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands