CRYPTOSPORIDIUM II™

A Monoclonal ELISA for the Detection of
Cryptosporidium Oocyst Antigen In Fecal Specimens

Catalog No. PT5014 (96 Tests)

[IVD] *In Vitro* Diagnostic Medical Device

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CRYPTOSPORIDIUM II™

INTENDED USE

The *Cryptosporidium II™* test is an enzyme immunoassay for the qualitative detection of *Cryptosporidium* oocyst antigen in human fecal specimens. It is indicated for use as an aid in the diagnosis of patients with diarrhea suspected of *Cryptosporidium* gastrointestinal infection.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician.

EXPLANATION

Cryptosporidium spp. is a protozoan parasite of vertebrates previously thought to cause diarrhea only in animals (1). In 1976, the first human infection was reported (2). Since then it has been found to be associated with diarrheal illness in most parts of the world and is a frequent cause of travelers' diarrhea (1,3).

The disease is transmitted by the thick-walled oocyst form, 2-6 µm in diameter, which is remarkably resistant to common disinfectants and routine chlorination of drinking water. Person-to-person transmission, especially among children, is common (4). *Cryptosporidium* has little or no host specificity, and animals such as rodents, cattle and domestic pets serve as a reservoir for zoonotic transmission to humans (1,5). Such transmission occurs either by direct contact or by contamination of water supplies with fecal matter (1,6-8). Cryptosporidiosis is a serious opportunistic infection in patients with acquired immunodeficiency syndrome (AIDS) and is potentially a sexually transmitted disease (6,9). Clinical manifestations of cryptosporidiosis include cholera-like diarrhea, abdominal pain, nausea, vomiting and weight loss. In normal persons the infection is usually self-limiting and of short term. In AIDS and other immunosuppressed patients, cryptosporidiosis can result in prolonged and life-threatening illness due to excessive fluid loss. In these patients the infection may also spread to the respiratory and biliary tracts (6).

BIOLOGICAL PRINCIPLES

The *Cryptosporidium II™* test uses monoclonal and polyclonal antibodies to *Cryptosporidium* oocyst antigen. The *Microassay Plate* in the kit contains an immobilized monoclonal antibody against *Cryptosporidium* oocyst antigen, and the *Conjugate* consists of a polyclonal antibody against *Cryptosporidium* oocyst antigen. In the assay, an aliquot of a diluted fecal specimen is transferred to a microassay well. If *Cryptosporidium* oocyst antigen is present, it binds to the immobilized monoclonal antibody. Upon addition, *Conjugate* then binds to the antigen/antibody complex. Any unbound materials are removed during the washing steps. Following the addition of substrate, a color is detected due to the enzyme-antibody-antigen complexes that formed in the presence of *Cryptosporidium* oocyst antigen and *Conjugate*.

MATERIALS PROVIDED

CONJ [ENZ]

Conjugate, 7 mL (Rabbit polyclonal antibody to a *Cryptosporidium* oocyst antigen coupled to horseradish peroxidase in a protein buffered solution containing 0.02% thimerosal)*

DIL [SPE]

Diluent, 50 mL (Buffered protein solution containing 0.02% thimerosal)*. The **Diluent** is also to be used as the negative control solution (see TEST PROCEDURE).

H₂SO₄ 0.6N

Stop Solution, 7 mL (0.6N sulfuric acid). Caution: Avoid contact with skin. Flush with water immediately if contact occurs.

Signal Word: Danger

H314: Causes severe skin burns and eye damage

P260, P264, P280, P301, P330, P331, P303, P361, P353, P363, P304, P340, P310, P321, P305, P351, P338, P501

CONTROL +

Positive Control, 3.5 mL (Heat-inactivated bovine fecal material containing *Cryptosporidium* oocyst antigen in a protein buffered solution with 0.02% thimerosal)*

SUBS REAG **Substrate**, 14 mL (solution containing tetramethylbenzidine and peroxide)
WASHBUF 20X **Wash Buffer Concentrate**, 50 mL (20X concentrate containing phosphate-buffered saline, detergent and 0.2% thimerosal)*

Signal Word: Warning

H373: May cause damage to organs through prolonged or repeated exposure P260, P314, P501

MA PLT **Microassay Plate**, 12 strips, each consisting of 8 wells coated with monoclonal antibody to *Cryptosporidium* oocyst antigen (stored with desiccant)

2 plastic adhesive sheets

*contains mercury

100 graduated disposable pipettes

Vortex mixer
Discard container
Microcentrifuge tubes
Applicator sticks

MATERIALS REQUIRED BUT NOT PROVIDED

Squirt bottle for wash reagent

950 mL distilled water for diluting wash reagent

Absorbent paper

ELISA reader capable of reading at 450 nm or 450/620 nm

PRECAUTIONS

1. Rx Only - Prescription Only
2. Reagents from different kits should not be mixed. Do not use a kit past the assigned expiration date.
3. Reagents should be brought to room temperature before use.
4. Caps and tips are color-coded; do not mix.
5. Each component in the kit should be inspected for any signs of leakage. Upon arrival, the kit should be inspected to ensure that components are not frozen or warm to the touch due to improper shipping conditions.
6. When handling assay wells, avoid scratching the bottom of the wells because this may result in elevated absorbance readings.
7. Hold dropper bottles vertically when dispensing to ensure proper drop size.
8. Unused microassay wells must be immediately placed back inside the resealable pouch with the desiccant and sealed to protect from moisture.
9. Perform the washing procedure as directed to avoid high background reactions.
10. The *Substrate* is light sensitive and should be protected from direct sunlight or UV sources. If the *Substrate* is exposed to light and develops a color, it must be discarded.
11. Optimal results are obtained by following the specified test procedure. The concentrations, incubation conditions, and processing specifications have been optimized for sensitivity and specificity. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test.
12. Microbial contamination of reagents may decrease the accuracy of the assay. Avoid microbial contamination of reagents by using sterile disposable pipettes if removing aliquots from reagent bottles.
13. Handle specimens and used microassay wells as if capable of transmitting infectious agents. Wear disposable gloves when doing the test.
14. The 20X Wash Buffer Concentrate contains 0.2% Thimerosal as a preservative. Once diluted to normal use concentration this solution is classified as non-hazardous. The Stop Solution contains 0.6N sulfuric acid. Flush with water immediately if contact occurs. Take off contaminated clothing and wash it before reuse. Handle reagents according to existing regulations for laboratory safety and good laboratory practice. Safety Data Sheets for this product are available upon request, contact technical support.

15. Follow your national, regional, and local ordinances accordingly for waste disposal regulations. Do not place in trash, dispose of as hazardous waste.

SHELF LIFE AND STORAGE

The expiration date of this kit is given on the box label. Expiration dates for each component are listed on individual labels. The kit should be stored between 2° and 8°C and should be returned to the refrigerator as soon as possible after use.

SPECIMEN HANDLING

1. Standard collection and handling procedures used in-house for fecal specimens for culture are appropriate. No modification of collection methods used for standard microscopic O&P examinations is needed.
2. Fecal specimens may be used as unpreserved or frozen, or in preservation media of 10% buffered formalin, Sodium Acetate Formalin (SAF) or transport medium such as Cary Blair and C&S.
3. Unpreserved specimens should be kept between 2° and 8°C and tested within 24 hours of collection. Specimens that cannot be tested within this time should be stored at -20°C or less until tested.
4. Preserved specimens may be kept at room temperature and tested within 18 months of collection.
5. Concentration steps are not needed (or recommended) for fecal specimens.
6. Make sure that specimens are thoroughly mixed (vortexed) prior to performing the test. This includes complete mixing of the specimen prior to transfer to the *Diluent* and/or microassay well.
7. All dilutions of specimens must be made with *Diluent*.

REAGENT PREPARATION

1. The contents of the kit should be brought to room temperature before use.
2. Prepare 1X *Wash Solution*. The *Wash Buffer Concentrate* is supplied as a 20X concentrate (a precipitate may be noticed). It should be diluted to a total volume of 1 liter by adding 50 mL of the concentrate to 950 mL of distilled water. The 1X *Wash Solution* can be stored between 2° and 8°C.
3. All reagents, with the exception of the *Wash Buffer Concentrate*, are supplied in ready-to-use bottles. Reagents can be dispensed directly from the dropper bottles or decanted for use with multichannel pipettes. If excess reagent has been decanted, the excess should be discarded. Do not pour back into the bottle. The *Substrate* should be stored in and used from the light-protected bottle in which it is supplied. If an aliquot is removed from the original bottle for any reason, do not return unused *Substrate* to the original bottle.

SPECIMEN PREPARATION

1. Fresh/Frozen Fecals: Frozen fecal specimens should be thawed. Add 400 µL of *Diluent* to a microcentrifuge tube (one tube per sample), then add 100 µL (2nd graduation mark on pipette) of fecal sample to the tube and mix well. If the specimen cannot be pipetted, use an applicator stick to transfer approximately 0.1 gram of feces. This is about the size of a small pea (about 4 mm in diameter).
2. Preserved Fecals: Mix (vortex) contents of container thoroughly before transferring specimen. No further processing or dilution is necessary.

TEST PROCEDURE

1. Prepare two control wells to serve as positive and negative controls each time the test is performed. For the positive control, shake the *Cryptosporidium Positive Control* bottle (white cap) for several seconds, then add one drop (50µL) to the positive control well, and add 2 drops (100 µL) of *Diluent* to the negative control well.

2. If fresh or frozen samples (i.e. unpreserved) are to be used, make a dilution as stated in SPECIMEN PREPARATION above.
3. Transfer 100 µL of *Diluent* to each test well of the *Microassay Plate*. Using plastic pipettes, transfer 1 drop (50 µL, 1st graduation mark on pipette) of sample (preserved or diluted as above) to each test well already containing *Diluent* and gently tap the wells to mix. Seal with a plastic adhesive sheet and incubate for 1 hour at room temperature.
4. Shake out contents of assay wells into a discard pan. Wash each well using the 1X *Wash Solution* in a squirt bottle with a fine-tipped nozzle, directing the *Wash Solution* to the bottom of the well with force. Fill the wells, shake the wash solution out of the wells into a discard pan and slap the plate hard onto a dry paper towel. Repeat the washing step 4X (for a total of 5 washes). If any fecal material remains in the wells, wash the plate until it appears clean. (Go to www.techlab.com for additional information on ELISA well washing)
5. After washing, completely remove any residual liquid in the wells by striking the plate onto a dry paper towel until no liquid comes out. Dispose of paper towels and specimen containers properly.
6. Add 1 drop (50 µL) of *Conjugate* (red cap) to each well and gently tap to mix. Seal with a plastic adhesive sheet. Incubate the wells for 30 minutes at room temperature.
7. Repeat the washing procedure (Steps 4 and 5).
8. Add 2 drops (100 µL) of *Substrate* (blue cap) to each well. Gently tap the wells to mix. Incubate the wells at room temperature for 10 minutes.
9. Add 1 drop (50 µL) of *Stop Solution* (yellow cap) to each well. Gently tap the wells to mix and wait 2 minutes before reading. The addition of the *Stop Solution* converts the blue color to a yellow color, which may be quantitated by measuring the absorbance at 450 nm on a microplate ELISA reader. The instrument should be blanked against air. Wipe the underside of each well before measuring the absorbance. If a dual reader is used, blank against air at 620 nm and read at 450 nm. Visual readings should be noted. Read within ten minutes after adding *Stop Solution*.

QUALITY CONTROL

1. A positive and a negative control must be run with each series of test specimens. The positive control demonstrates that the assay is functioning properly for the detection of *Cryptosporidium* oocyst antigen in fecal specimens. The negative control demonstrates that the assay is reacting specifically.
2. Each positive control well should be an easily visible yellow color and should give an absorbance of 0.500 or higher. Any well that gives a positive reading without visible color should be repositioned, wiped on the underside of the well, and read again.
3. Negative control wells should be colorless or may have a faint yellow color but they must have an absorbance value of <0.150 when read on single wavelength (OD_{450nm}) or <0.090 when read on dual wavelength (OD_{450/620nm}).
4. Test results are not valid unless the performance characteristics of the positive and negative controls are met. If these results are not observed, call Technical Services.
5. Test results along with control absorbance values should be recorded and reported according to in-house procedures and should be stored according to in-house procedures for future reference.

INTERPRETATION OF RESULTS

Absorbance

			<u>Visual Color</u>		<u>Interpretation</u>
450 nm	450/620 nm	< 0.150	Clear to slight yellow		Negative – Antigen is not detected.
≥0.150	≥0.090		Pale yellow to strong yellow		Below the detectable limits of the assay. Positive - Specimen contains <i>Cryptosporidium</i> antigen

Visual Interpretation

Negative: Any sample that is colorless or resembles the negative control well in intensity of color. Antigen is not detected. Below the detectable limits of the assay.

Positive: Any sample that is obviously more yellow than the negative control well. Specimen contains *Cryptosporidium* antigen.

NOTE: The negative control, as well as some test wells, may show some slight yellow color. A sample well must be obviously more yellow than the negative control well to be called a positive result.

Spectrophotometric Interpretation

- Determine the absorbance value of the negative control.

The negative control reading should be $< 0.150 \text{ OD}_{450}$ or $< 0.090 \text{ OD}_{450/620}$. If not, the test is invalid and should be repeated, paying attention to the wash procedure.

- The reading for the *Positive Control* should be ≥ 0.500 .

- Test results

Negative: < 0.150 (absorbance at 450 nm) or < 0.090 (absorbance at 450/620 nm). Antigen is not detected. Below the detectable limits of the assay.

Positive: ≥ 0.150 (absorbance at 450 nm) or ≥ 0.090 (absorbance at 450/620 nm). Specimen contains *Cryptosporidium* antigen.

LIMITATIONS OF THE PROCEDURE

- The *CRYPTOSPORIDIUM II™* test detects the presence of *Cryptosporidium* oocyst antigen in fecal specimens.
- The test results should be interpreted by a physician in consideration of other laboratory results and clinical history.
- Concentrated fecal specimens should not be tested in the *CRYPTOSPORIDIUM II™* test and will not give accurate results.
- The *CRYPTOSPORIDIUM II™* test is intended for the qualitative detection of *Cryptosporidium* oocyst antigen in fecal specimens. It has not been evaluated for quantitative determinations of organism load, and the magnitude of the absorbance value is not intended to correlate with organism load.

PERFORMANCE CHARACTERISTICS

An independent diagnostics laboratory evaluated the *CRYPTOSPORIDIUM II™* test for its ability to detect *Cryptosporidium* antigens in preserved fecal specimens. Of the 185 specimens analyzed by microscopy 44 were positive for *Cryptosporidium*. Each specimen was subsequently tested using the *CRYPTOSPORIDIUM II™* test. Results are summarized in Table 1.

Table 1. Comparison of the *CRYPTOSPORIDIUM II™* test with microscopy for detection of *Cryptosporidium*.

Total Samples (n = 185)		Microscopy	
		Positive	Negative
<i>CRYPTOSPORIDIUM II™</i>	Positive	43	0
	Negative	1	141

		95% CI
Sensitivity	97.7%	86.5% - 99.9%
Specificity	100%	96.7% - 100%
Predictive Positive Value	100%	89.8% - 100%
Predictive Negative Value	99.3%	95.6% - 99.9%
Correlation	99.5%	99.3% - 99.6%

CROSS REACTIVITY

An independent diagnostics laboratory evaluated the *CRYPTOSPORIDIUM II™* test using fecal specimens found to be positive for a variety of intestinal pathogens. No cross reactivity was observed with fecal specimens that contained any of the pathogens listed below. The number of specimens tested with each organism is shown in parentheses.

<i>Ascaris lumbricoides</i> eggs (5)	<i>Blastocystis hominis</i> (35)	<i>Chilomastix mesnili</i> (1)
<i>Cyclospora cayetanensis</i> (2)	<i>Dientamoeba fragilis</i> (23)	<i>Diphyllobothrium latum</i> eggs (1)
<i>Endolimax nana</i> (37)	<i>Entamoeba coli</i> (22)	<i>Entamoeba hartmanni</i> (3)
<i>Entamoeba histolytica/E. dispar</i> (7)	Hookworm eggs (1)	<i>Enterobius vermicularis</i> eggs (4)
<i>Giardia lamblia</i> (40)	<i>Strongyloides stercoralis</i> larvae (3)	<i>Hymenolepis nana</i> eggs (4)
<i>Iodamoeba butschlii</i> (4)	<i>Trichuris trichiura</i> eggs (2)	
<i>Taenia</i> spp. eggs (1)		

The *CRYPTOSPORIDIUM II™* test was evaluated for crossreactivity with the bacterial and viral strains listed below. None of the strains were shown to crossreact with the *CRYPTOSPORIDIUM II™* test.

<i>Aeromonas hydrophila</i>	<i>Bacillus cereus</i>	<i>Bacillus subtilis</i>
<i>Bacteroides fragilis</i>	<i>Campylobacter coli</i>	<i>Campylobacter fetus</i>
<i>Clostridium bifermentans</i>	<i>Clostridium difficile</i>	<i>Enterococcus faecalis</i>
<i>Escherichia coli</i>	<i>Escherichia coli</i> 0157H7	
<i>Escherichia coli</i> EIEC (enteroinvasive)		
<i>Escherichia coli</i> EPEC (enteropathogenic)		
<i>Escherichia coli</i> ETEC (enterotoxic)		<i>Klebsiella pneumoniae</i>
<i>Salmonella typhimurium</i>	<i>Shigella dysenteriae</i>	<i>Shigella flexneri</i>
<i>Staphylococcus aureus</i>	<i>Staphylococcus aureus</i> (Cowan's)	
<i>Staphylococcus epidermidis</i>	<i>Vibrio parahaemolyticus</i>	<i>Yersinia enterocolitica</i>
Adenovirus types 1, 2, 3, 5, 40, 41		Coxsackievirus B2, B3, B4, B5
Echovirus 9, 11, 18, 22, 33		Enterovirus types 68, 69, 70, 71
Human coronavirus		

INTERFERING SUBSTANCES

The following substances had no effect on positive or negative test results analyzed at the concentrations indicated: mucin (3.5% w/v), human blood (40% w/v), Imodium® (5% w/v), Kaopectate® (5 mg/mL), Pepto-Bismol® (5% w/v), fecal fat (stearic acid 40% w/v), Metronidazole (0.25% w/v), Vancomycin (0.25% w/v).

SPECIMEN PRESERVATIVE COMPATIBILITY

The *CRYPTOSPORIDIUM II™* test was examined for compatibility with preserved fecal specimens during an in-house evaluation. No effect on sensitivity or specificity was observed with specimens preserved in 10% buffered formalin or Sodium Acetate Formalin (SAF) or with specimens stored in transport media such as Cary Blair and C&S.

PRECISION – INTRA-ASSAY

The intra-assay % Coefficient of Variation (CV) was determined for the *CRYPTOSPORIDIUM II™* test by analyzing four fecal specimens positive for *Cryptosporidium* and four negative fecal specimens. Each specimen was assayed in eight test wells. The average %CV for positive specimens was 4.91. The average %CV for negative specimens was 2.71.

PRECISION - INTER-ASSAY

The inter-assay % Coefficient of Variation (CV) was determined for the *CRYPTOSPORIDIUM II™* test by analyzing four fecal specimens positive for *Cryptosporidium* and for negative fecal specimens. Each specimen was assayed during five independent assay runs conducted over a four-day period using the *CRYPTOSPORIDIUM II™* kit. The average %CV for positive specimens was 18.14. The average %CV for negative specimens was 2.51.

EXPECTED VALUES

Normal healthy individuals should not be infected with *Cryptosporidium* and should test negative in the *CRYPTOSPORIDIUM II™* test. A positive test result in the *CRYPTOSPORIDIUM II™* test indicates that the person is shedding detectable amounts of *Cryptosporidium* oocyst antigen. The incidence of *Cryptosporidium* infection varies significantly between populations and geographic regions. In general, laboratory-confirmed incidence of cryptosporidiosis in developed countries ranges from 1 to 2% overall, with a higher incidence in children (10).

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Cryptosporidium II, the TECHLAB Logo and TECHLAB are trademarks of TECHLAB, Inc., under license.
All other trademarks referenced are trademarks of their respective owners.