# TECHLAB® C. DIFF CHEK™ - 60

An ELISA for the detection of *Clostridium difficile* Antigen  
Catalog No. TL5025 (96 Tests) and T5025B (960 Tests)

**IVD** In Vitro Diagnostic Medical Device

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**C. DIFF CHEK™ - 60****INTENDED USE**

The C. *DIFF CHEK™* - 60 test is an enzyme immunoassay for use as a screening test to detect *Clostridium difficile* antigen, glutamate dehydrogenase, in fecal specimens from persons suspected of having C. *difficile* disease. The test does not distinguish toxigenic from nontoxigenic strains of C. *difficile*. With the use of additional tests that detect C. *difficile* toxins, this test is to be used as an aid in the diagnosis of C. *difficile* disease. As with other C. *difficile* tests, results should be considered in conjunction with patient history.

**Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician**

**EXPLANATION**

After treatment with antibiotics, many patients develop gastrointestinal problems ranging from mild diarrhea to severe pseudomembranous colitis. Many cases of the milder forms of gastrointestinal illness and most cases of pseudomembranous colitis are caused by toxigenic strains of *Clostridium difficile* (1). This organism is an opportunistic anaerobic bacterium that grows in the intestine once the normal flora has been altered by the antibiotic. Toxigenic strains of C. *difficile* carry the genes encoding the toxins while non-toxigenic strains do not carry the toxin genes. The disease results from the toxins that the toxigenic organism produces. The clinical symptoms associated with the disease are believed to be primarily due to toxin A, which is a tissue-damaging enterotoxin (2,3). C. *difficile* also produces a second toxin, designated toxin B. Toxin B, which has been referred to as the cytotoxin of the organism, is the toxin detected by the tissue culture assay currently used by many laboratories. Toxigenic C. *difficile* strains produce both toxins or only toxin B (4,5,6,7). The glutamate dehydrogenase of C. *difficile* is a good antigen marker for the organism in feces because it is produced in high amounts by all strains, toxigenic or non-toxigenic (8-11). The antigen can be detected in fecal specimens by using the C. *DIFF CHEK™* - 60 test. A positive result in the test, which is highly specific for the glutamate dehydrogenase of C. *difficile*, confirms the presence of this organism in a fecal specimen; a negative result indicates the absence of the organism. A positive result should be followed by a toxin-specific test to confirm the presence of toxigenic C. *difficile*.

**PRINCIPLE OF THE TEST**

The C. *DIFF CHEK™* - 60 test uses antibodies specific for the glutamate dehydrogenase of C. *difficile*. The *Microassay Plate* supplied with the kit contains immobilized polyclonal antibody against the antigen. The *Conjugate* consists of a highly specific monoclonal antibody conjugated to horseradish peroxidase. In the assay, an aliquot of a fecal specimen is emulsified in the *Diluent* and the diluted specimen is transferred to the microwell containing the *Conjugate*. If the antigen is present in the specimen, it will bind to the *Conjugate* and to the immobilized polyclonal antibody during the incubation phase. Any unbound material is removed during the washing steps. Following the addition of *Substrate*, a color is detected due to the enzyme-antibody-antigen complexes that form in the presence of antigen.

**MATERIALS PROVIDED****DIL** **SPE**

**Diluent (40 mL each)** – buffered protein solution containing preservative\*  
H412 - Harmful to aquatic life with long lasting effects.  
P273, P501

**CONJ** **ENZ**

**Conjugate (7 mL each)** – antigen-specific mouse monoclonal antibody coupled to horseradish peroxidase, in a buffered protein solution (contains 0.05% ProClin® 300)

Signal Word: Warning

H317: May cause an allergic skin reaction

P261, P272, P280, P302, P352, P333, P313, P321, P362, P363, P364, P501

**SUBS** **REAG**

**Substrate (14 mL each)** – solution containing tetramethylbenzidine and peroxide

|                                         |                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                        |
|-----------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>CONTROL</b> +                        | <b>Positive Control (4.5 mL each)</b> – antigen in a buffered protein solution containing preservative                                                                                                                                                                                                                         |                                                                                                                                                                        |
| <b>WASHBUF</b> 20X                      | <b>Wash Buffer Concentrate (50 mL each)</b> – 20X Concentrate containing phosphate-buffered saline, detergent and 0.2% thimerosal*<br>Signal Word: Warning<br>H373: May cause damage to organs through prolonged or repeated exposure<br>H411: Toxic to aquatic life with long lasting effects<br>P260, P273, P314, P391, P501 | <br> |
| <b>H<sub>2</sub>SO<sub>4</sub></b> 0.6N | <b>Stop Solution (7 mL each)</b> – 0.6 N Sulfuric acid. CAUTION: Avoid contact with skin; flush with water immediately if contact occurs<br>Signal Word: Danger<br>H314: Causes severe skin burns and eye damage<br>P260, P264, P280, P301, P330, P331, P303, P361, P353, P363, P304, P340, P310, P321, P305, P351, P338, P501 |                                                                                       |
| <b>MA</b> <b>PLT</b>                    | <b>Microassay Plate</b> – 12 strips each, each consisting of 8 wells coated with antibody specific for antigen                                                                                                                                                                                                                 |                                                                                                                                                                        |

\*contains mercury

**IVD** FOR IN VITRO DIAGNOSTIC USE

## ACCESSORIES

Applicator Sticks (50/500)  
Plastic adhesive sheets (2/20)

Disposable plastic pipettes (100/1000)  
Wash Solution Label (1)

## MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

|                          |                        |                                         |
|--------------------------|------------------------|-----------------------------------------|
| <i>Wash bottle</i>       | <i>Timer</i>           | <i>Paper towels or absorbent sheets</i> |
| <i>Discard container</i> | <i>Distilled water</i> | <i>Vortex mixer</i>                     |

*Microplate Reader (Spectrophotometer)* – capable of reading dual wavelength at 450/620 nm or single wavelength at 450 nm (a dual wavelength plate reader is recommended; absorbances should be measured at 450 nm and Referenced at 620 nm)

*Refrigerator* - set between 2° and 8°C      *Incubator* - set at 37°C ± 2°C

## SHELF LIFE AND STORAGE

The expiration date of the kit is given on the label. Expiration dates for each component are listed on the individual labels. The kit containing the reagents with designated shelf life should be stored between 2° and 8°C and should be returned to the refrigerator as soon as possible after use.

## PRECAUTIONS

1. Rx Only – Prescription Only
2. For professional use only.
3. Reagents from different kits should not be mixed. Do not use a kit past the expiration date.
4. Each component in the kit should be inspected for any signs of leakage. Upon arrival, inspect the kit to ensure that components are not frozen or warm to the touch due to improper shipping conditions.
5. Reagents should be at room temperature before use.
6. Caps and tips are color-coded; do NOT mix!
7. When handling microwells, avoid scratching the bottom of the wells because this may result in elevated absorbance readings.
8. Hold dropper bottles vertically to ensure proper drop size.
9. Handle specimens and used microwells as if capable of transmitting infectious agents. Wear gloves when doing the test.
10. The 20X *Wash Buffer Concentrate* contains 0.2% Thimerosal as a preservative. Once diluted to normal use concentration this solution is classified as non-hazardous.

11. The *Stop Solution* contains 0.6 N sulfuric acid. Flush with water immediately if contact occurs.
12. *Unused microwells must be placed back inside of the resealable pouch with the desiccant to protect them from moisture.* Check the desiccant pack before using the microwells. The color indicator on the desiccant pack should be blue. If the color turns pink, the quality of the microwells may be compromised. Please do not use the microwells stored with pink desiccant.
13. Specimens that have been preserved in 10% Formalin, merthiolate formalin, sodium acetate formalin, or polyvinyl alcohol, or specimens that are in transport media such as Cary Blair or C&S cannot be used.
14. Perform the washing procedure as directed to avoid high background reactions.
15. Use fecal specimens within 72 hours of collection to obtain optimal results. Specimens that are frozen ( $\leq -10^{\circ}\text{C}$ ) may lose activity due to freezing and thawing.
16. The *Substrate* is light sensitive and should be protected from direct sunlight or UV sources. The appearance of the reagent should be clear and colorless.
17. The *Positive Control* confirms the reactivity of the other reagents associated with the assay, and is not intended to ensure precision at the analytical assay cut-off.
18. Optimal results are obtained by following the specified test procedure. The concentrations, incubation conditions, and processing specifications have been optimized for sensitivity and specificity. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test.
19. The *Conjugate* reagent contains 0.05% ProClin® 300 as a preservative. Although the concentration is low, ProClin® 300 is known to be harmful (skin sensitization may occur). If skin sensitization/irritation or rash occurs, get medical advice/attention.
20. Take off contaminated clothing and wash it before reuse. Handle reagents according to existing regulations for laboratory safety and good laboratory practice. Safety Data Sheets for this product are available upon request, contact technical support.
21. Follow your national, regional, and local ordinances accordingly for waste disposal regulations. Do not place in trash, dispose of as hazardous waste.

## COLLECTION AND HANDLING OF FECAL SPECIMENS

**NOTE:** Standard collection and handling procedures used in-house for fecal specimens are appropriate. Specimens should be stored between  $2^{\circ}$  and  $8^{\circ}\text{C}$ . Whenever possible, test specimens that are less than 24 hours old. Freeze the specimens ( $\leq -10^{\circ}\text{C}$ ) if the test cannot be performed within 72 hours of collection. Freezing and thawing of the specimen, especially multiple times, results in loss of activity due to degradation of the antigen. Fecal specimens that have been preserved in 10% Formalin, merthiolate formalin, sodium acetate formalin, or polyvinyl alcohol, or specimens that are in transport media such as Cary Blair or C&S cannot be used. Perform the washing procedure as directed to avoid high background reactions. Make sure that specimens are thoroughly mixed (vortexed) prior to performing the assay. This includes complete mixing of the specimen prior to transfer to *Diluent* as well as complete mixing of the diluted specimen prior to transfer to the microwell. Storing of fecal specimens in the *Diluent* is **not** recommended. Test immediately once the fecal specimen is diluted in *Diluent*. Disposable pipettes included in the accessories are graduated at 50, 100, 200, and 300  $\mu\text{L}$ .

1. Set up one dilution tube for each specimen to be tested. Add 200  $\mu\text{L}$  *Diluent* to each tube. *Label the tube directly on the side.*
2. For **formed fecal specimens**, use an applicator stick to transfer the fecal specimen to the tube. Transfer an amount equal to about 3 mm in diameter with the applicator stick to the *Diluent*. For **liquid fecal specimens**, use a plastic pipette to transfer 50  $\mu\text{L}$  specimen to the tube. Make sure the liquid specimens are evenly suspended (vortexed) before transferring.
3. Vortex the tubes for 10 seconds before transferring the diluted specimen to the microwell. This ensures thorough mixing of the specimen.
4. Automated equipment or semi-automated washing equipment may be used with specimens that have been centrifuged ( $5000 \times g$  for 10 minutes) to remove particulate

matter. Double the amount of fecal specimen and the amount of *Diluent* to ensure sufficient supernatant for testing.

## PRELIMINARY PREPARATIONS

1. **IMPORTANT:** All reagents must be at room temperature prior to use in the assay.
2. **Prepare 1X Wash Solution.** The *Wash Buffer Concentrate* is supplied as a 20X concentrate (a precipitate may be noticed). It should be diluted to a total volume of 1 liter by adding 50 mL of the concentrate to 950 mL of distilled water. Label the bottle. Store any unused 1X *Wash Solution* between 2° and 8°C.
3. **Assay Strip Preparation.** Each strip contains 8 wells coated with polyclonal antibodies specific for *C. difficile* antigen. Each specimen or control will employ one of these coated wells. Determine the number of wells to be used. Avoid contact with the base of the wells. Microwells not used must be returned to the resealable bag with desiccant and carefully resealed in the bag.

## TEST PROCEDURE

1. Add 1 drop (50 µL) of *Conjugate* (red cap) to each well. Be sure to hold each bottle vertically when adding the drops. Use 1 well for each fecal specimen, 1 well for the positive control and 1 well for the negative control. Identification marks may be written directly on side of well.
2. Transfer 100 µL (2 drops using a transfer pipette from the accessory kit) of diluted specimen to the microwell. Add 1 to 2 drops (50 to 100 µL) of the *Positive Control* (black cap) to the positive control well and 1 to 2 drops (50 to 100 µL) of the *Diluent* (negative control) to the negative control well. Gently tap to mix.  
**Note:** If using semi-automated or automated washing equipment, specimens must be centrifuged (5000 x g for 10 minutes) to remove particulate matter.
3. Cut the adhesive plastic sheet to the size necessary to cover the wells. Cover the wells and incubate them at 37°C ± 2°C for 50 minutes.
4. Shake out the contents of the microwells into a discard pan.
5. Wash each well using the 1X *Wash Solution* in a squirt bottle with a fine-tipped nozzle, directing the *Wash Solution* to the bottom of the well with force. Fill the wells, then shake the *Wash Solution* out of the well into a discard pan. Slap the inverted plate on a dry paper towel.
6. Repeat step 5 **four times** using a dry paper towel each time. If any particulate matter is seen in the wells, continue washing until all the particulate matter is removed. **Note:** If using semi-automated or automated washing equipment, add 350 µL of 1X *Wash Solution* to each well. Wash for a total of **7** times. If particulate matter is seen in the well, continue washing until all particulate matter is removed.
7. After washing, completely remove any residual liquid in the wells by striking the plate once again onto a dry paper towel until no liquid comes out. *Dispose of paper towels & specimen containers properly.*
8. Add 2 drops (100 µL) of *Substrate* (blue cap) to each well. Gently tap the wells to mix the substrate. Incubate the wells at room temperature for 10 minutes. Gently tap the wells at 5 minutes.
9. Add 1 drop (50 µL) of *Stop Solution* (yellow cap) to each well. The addition of the *Stop Solution* converts the blue color to a yellow color that may be quantitated by measuring the optical density at 450 nm on a microassay plate reader. The instrument should be blanked against air. If a dual wavelength reader is used, blank against air at 620 nm and read at 450 nm. Wipe the underside of each well before measuring the optical density. If a microassay plate reader is unavailable, the test may be read visually in good light against a white background. Read within two to ten minutes after adding *Stop Solution*.

## OPTIONAL TEST PROCEDURE/RAPID FORMAT

Perform the regular test procedure according to the instructions provided above replacing the 50 minute incubation at 37°C ± 2°C with 20 minutes at 37°C using the *Stat Fax 2200 Incubator/Shaker* or equivalent incubator/shaker. If the *Stat Fax 2200 Incubator/Shaker* is

6

used, set the shaker to speed 7 and the temperature to 37°C. If other shakers are used, a speed of 1500 rpm is recommended. Shaking should not cause spillage. If spillage occurs, reduce the speed accordingly (12).

### QUALITY CONTROL

1. A positive and negative control must be run with each series of test specimens.
2. Positive and negative controls must fall in their respective ranges or the test is not valid.
  - a) Positive Control must be a visible yellow color.  
If read on a microplate reader, the OD at 450 or 450/620 nm must be  $\geq 0.500$ .
  - b) Negative Control must be visually clear.  
If read on a microplate reader, the OD at 450 nm must be  $< 0.120$ .  
If read at 450/620 nm the absorbance must be  $< 0.080$ .
3. Wells that are clear visually but give absorbance  $> 0.120$  should be wiped on the underside and re-measured.
4. Visual readings must be taken in good light against a white background.
5. A sample that yields a weak positive result (i.e.,  $< 0.200$  O.D.) and is adjacent to a strong positive should be re-tested to assure carryover did not occur.

### INTERPRETATION OF RESULTS

1. **Visual Reading**  
Negative = Colorless  
Positive = Any yellow color
2. **Spectrophotometric Single Wavelength at 450 nm**  
Negative = OD  $< 0.120$   
Positive = OD  $\geq 0.120$
3. **Spectrophotometric Dual Wavelength at 450/620 nm**  
Negative = OD  $< 0.080$   
Positive = OD  $\geq 0.080$

A positive test result indicates that *C. difficile* antigen is present in the fecal specimen.

A negative test result indicates that *C. difficile* antigen is not present in the fecal specimen.

### LIMITATIONS OF THE *C. DIFF CHEK*™ - 60 TEST

1. The *C. DIFF CHEK*™ - 60 test is used to detect *C. difficile* antigen in fecal specimens. The test confirms the presence of antigen in feces and this information should be taken under consideration by the physician in light of the clinical history of the patient and the results of toxin detecting tests. A positive result with the *C. DIFF CHEK*™ - 60 test DOES NOT confirm the existence of toxigenic strain of *C. difficile*. If a positive result is seen using the *C. DIFF CHEK*™ - 60 test, it is suggested that the *C. DIFFICILE TOX A/B II*™ test (30397) be used to confirm the presence of toxin. Alternatively, tissue culture cytotoxicity assay may be performed to verify the presence of the toxins.
2. Fecal specimens represent an extremely complex clinical specimen. Optimal results with the *C. DIFF CHEK*™ - 60 test are obtained with specimens that are less than 24 hours old. Most undiluted specimens can be stored between 2° and 8°C for 72 hours before significant degradation of the antigen is noted. If specimens are not assayed within this time period, they maybe frozen and thawed. However, repeated freezing and thawing may result in loss in the immunoreactivity of the antigen.
3. Some specimens may give weak reactions. This may be due to a number of factors such as the presence of low levels of antigen or the presence of binding substances or inactivating enzymes in the feces. *Under these conditions, the specimen should be re-tested or a fresh specimen should be tested.* Additional tests that may be used in conjunction with the *C. DIFF CHEK*™ - 60 test include isolation of the organism on selective media, toxin-specific ELISAs, or tissue culture cytotoxicity assays for the detection of *C. difficile* toxin. The isolation of the organism does not confirm that the organism is toxigenic. Additional testing of the isolate using a toxin-specific ELISA or tissue culture assay must confirm this.
4. Fecal specimens preserved in 10% Formalin, merthiolate formalin, sodium acetate formalin, or polyvinyl alcohol, or in transport media such as Cary Blair, or C&S CANNOT be used.

## EXPECTED VALUES

Fecal specimens used in this study were submitted for routine laboratory testing for the presence of *C. difficile* and its toxins. Fecal specimens of all consistencies (solid, semisolid and liquid) were included in these studies. The prevalence of a positive *C. DIFF CHEK™* - 60 test at three independent study sites was: Site #2 (N=100), 20.0%; Site #3 (N=100), 33.0%; Site #4 (N=275), 14.5%.

## PERFORMANCE CHARACTERISTICS

### Clinical evaluation

The *C. DIFF CHEK™* - 60 test was compared to the bacterial culture test at three hospitals/reference laboratories and in-house at TECHLAB<sup>®</sup>. Specimens included in the evaluation were submitted to the clinical laboratory for routine testing. The presumptive bacterial culture test was performed according to the in-house procedure (13). The result is shown in Table 1 below.

Table 1. Summary of clinical performance comparing *C. DIFF CHEK™* - 60 test to presumptive bacterial culture.

| Summary<br>(n=642)                        | Presumptive Bacterial<br>Culture<br><b>positive</b> | Presumptive Bacterial<br>Culture<br><b>negative</b> |  |
|-------------------------------------------|-----------------------------------------------------|-----------------------------------------------------|--|
| <i>C. DIFF CHEK™</i> - 60 <b>positive</b> | 74                                                  | 47                                                  |  |
| <i>C. DIFF CHEK™</i> - 60 <b>negative</b> | 31                                                  | 490                                                 |  |

|                           | Estimated Value | 95% Confidence Interval |             |
|---------------------------|-----------------|-------------------------|-------------|
|                           |                 | Lower Limit             | Upper Limit |
| Sensitivity               | 70.5%           | 60.7%                   | 78.8%       |
| Specificity               | 91.2%           | 88.5%                   | 93.4%       |
| Negative Predictive Value | 91.2%           | 91.6%                   | 95.9%       |
| Correlation               | 91.2%           | 85.9%                   | 89.6%       |

Discrepant samples were resolved by purification of DNA from fecal specimens and polymerase chain reaction amplification of the *C. difficile* GDH gene. Thirty-eight of the 47 apparent false positive samples were PCR positive for the *C. difficile* GDH gene, and were considered as true positives. Nine remained false positives. Thirty of the 31 apparent false negative specimens were PCR negative for the GDH gene, and were considered true negatives. One remained a false negative. The results are presented in Table 2 below.

Table 2. Summary of clinical performance comparing *C. DIFF CHEK™* - 60 test to resolved bacterial culture.

| Summary<br>(n=642)                        | Resolved Bacterial<br>Culture<br><b>positive</b> | Resolved Bacterial<br>Culture<br><b>negative</b> |  |
|-------------------------------------------|--------------------------------------------------|--------------------------------------------------|--|
| <i>C. DIFF CHEK™</i> - 60 <b>positive</b> | 112                                              | 9                                                |  |
| <i>C. DIFF CHEK™</i> - 60 <b>negative</b> | 1                                                | 520                                              |  |

|             | Estimated Value | 95% Confidence Interval |             |
|-------------|-----------------|-------------------------|-------------|
|             |                 | Lower Limit             | Upper Limit |
| Correlation | 98.4%           | 98.2%                   | 98.6%       |

### Centrifugation

A total of 43 fecal specimens were evaluated in the *C. DIFF CHEK™* - 60 test to determine the effect of centrifugation on the test results. For the analysis, specimens were

diluted as described in the Package Insert, and the appropriate volume of each specimen was transferred into a microassay well for assay. The remaining diluted specimens were centrifuged (5,000 x g) for 1 minute and an appropriate volume of each supernatant fluid was transferred to a microassay well for assay. The results from the two methods were compared. The results demonstrated a correlation of 100% between noncentrifuged and centrifuged specimens.

### SPECIFICITY

Various organisms were examined for cross-reactivity in the *C. DIFF CHEK™* - 60 test. For the analysis, broth cultures mixed with *Diluent* were evaluated. Broth cultures at log phase containing  $>10^6$  CFU/mL were used. Organisms that did not react under any of the conditions are listed in Table 3 below. The only organisms that reacted were toxigenic and nontoxigenic *C. difficile*.

Table 3. Organisms that do not react in the *C. DIFF CHEK™* - 60 test

#### Bacterial and Yeast Isolates

|                                 |                                                    |                                     |
|---------------------------------|----------------------------------------------------|-------------------------------------|
| <i>Aeromonas hydrophila</i>     | <i>Clostridium histolyticum</i>                    | <i>Pseudomonas aeruginosa</i>       |
| <i>Bacillus cereus</i>          | <i>Clostridium novyi</i> types A, B, and C         | <i>Salmonella typhimurium</i>       |
| <i>Bacillus subtilis</i>        | <i>Clostridium perfringens</i> types A, B, D and E | <i>Shigella dysenteriae</i>         |
| <i>Bacteroides fragilis</i>     | <i>Clostridium septicum</i>                        | <i>Shigella flexneri</i>            |
| <i>Candida albicans</i>         | <i>Clostridium sordellii</i>                       | <i>Shigella sonnei</i>              |
| <i>Campylobacter coli</i>       | <i>Clostridium sporogenes</i>                      | <i>Staphylococcus aureus</i>        |
| <i>Campylobacter fetus</i>      | <i>Enterococcus faecalis</i>                       | <i>Staphylococcus epidermidis</i>   |
| <i>Campylobacter jejuni</i>     | <i>Escherichia coli</i>                            | <i>Staphylococcus saprophyticus</i> |
| <i>Clostridium bifermentans</i> | <i>Klebsiella pneumoniae</i>                       | <i>Veillonella parvula</i>          |
| <i>Clostridium butyricum</i>    | <i>Peptostreptococcus anaerobius</i>               | <i>Vibrio parahaemolyticus</i>      |
| <i>Clostridium innocuum</i>     | <i>Proteus vulgaris</i>                            | <i>Yersinia enterocolitica</i>      |
| <i>Clostridium haemolyticum</i> |                                                    |                                     |

#### Viruses

|                                     |                                   |                              |
|-------------------------------------|-----------------------------------|------------------------------|
| Adenovirus types 1, 2, 3, 5, 41, 4D | Echovirus types 9, 11, 18, 22, 33 | Enterovirus types 70, 71, 68 |
| Coxsackie virus B2, B3, B4, B5      | Enterovirus Toluca 1              | Human coronavirus 229E       |

### EFFECT OF FECAL CONSISTENCY

The reaction of fecal samples of varying consistencies in the *C. DIFF CHEK™* - 60 test and in resolved bacterial culture test is shown in Table 4 below. A total of 629 fecal samples of known consistency were included. The percentages of positive reactions using either the antigen test or the bacterial culture were similar in all three types of fecal specimens (liquid, semi-solid, and solid). All of the specimens were submitted for *C. difficile* testing. The basis of the submission was the clinical history of the patient and not the consistency of the specimen. The result shows the *C. DIFF CHEK™* - 60 test performed similarly to bacterial culture assay when testing fecals of different consistencies.

Table 4. Reaction of fecals of varying consistencies in the *C. DIFF CHEK™* - 60 test

| # of Specimens<br>(n=629)            | Liquid Fecals<br>(n=256) | Semi-solid<br>Fecals<br>(n=247) | Solid Fecals<br>(n=126) |
|--------------------------------------|--------------------------|---------------------------------|-------------------------|
| <i>C. DIFF CHEK™</i> - 60, positive  | 39 (15.2%)               | 49 (19.8%)                      | 29 (23.0%)              |
| Resolved Bacterial Culture, positive | 41 (16.0%)               | 49 (19.8%)                      | 25 (19.7%)              |
| <i>C. DIFF CHEK™</i> - 60, negative  | 217 (84.8%)              | 198 (80.2%)                     | 97 (77.0%)              |
| Resolved Bacterial Culture, negative | 215 (84.0%)              | 198 (80.2%)                     | 101 (80.3%)             |

### REPRODUCIBILITY AND PRECISION

Seven positive and five negative fecal specimens were sent to each of three independent laboratories for analysis using the *C. DIFF CHEK™* - 60 test. All specimens

were kept frozen at  $\leq -10^{\circ}\text{C}$  until the assay was performed. The results from each laboratory were compared with in-house results and found to be identical. The positive specimens were confirmed to be positive and the negative specimens were confirmed to be negative at each site.

The intra-assay Coefficient of Variation (CV) of the *C. DIFF CHEK™ - 60* test was determined by analyzing 48 positive and 48 negative control reactions along with 8 negative fecal specimens. The positive and negative reactions were tested using 1 or 2 drops of *Positive Control* or *Diluent*. Each fecal specimen was tested in 12 wells. The intra-assay CV was 2.0% to 2.4% with the positive control reactions, 9.3% to 9.4% for the negative control reactions, and 8.7% with the negative fecal specimens. The inter-assay CV was determined using 7 positive and 5 negative fecal specimens tested with three different lots of test kits. The CV ranged from 0.3% to 37.8%, with an average of 8.3%.

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© 2021 TECHLAB, Inc. All rights reserved.
C. DIFF CHEK, C. DIFFICILE TOX A/B II, the TECHLAB Logo, and TECHLAB are trademarks of TECHLAB, Inc. Proclin is a trademark of Rohm and Haas Company.