# ***Clostridium Perfringens* Enterotoxin Test**

An ELISA for detecting  
*Clostridium perfringens* enterotoxin in fecal specimens  
Catalog No. T5006 (96 Tests)

**IVD** *In Vitro* Diagnostic Medical Device

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## **C. perfringens Enterotoxin Test**

### **INTENDED USE**

The *Clostridium perfringens* Enterotoxin Test is an enzyme immunoassay for the rapid detection of enterotoxin produced by *Clostridium perfringens*. It can be used to detect enterotoxin in fecal specimens from persons suspected of having *C. perfringens* disease. The test is to be used as an aid in the diagnosis of *C. perfringens* enterotoxin-associated disease and results should be considered in conjunction with the patient history.

### **EXPLANATION**

Intestinal disease caused by *C. perfringens* enterotoxin is the result of food poisoning following the consumption of contaminated, precooked food which has been inadequately stored or reheated. The disease was first reported at the turn of the last century and has become increasingly well characterized since the mid-1940s. Diarrhea and abdominal cramps are typical features, whereas vomiting and death are uncommon. Symptoms appear from 12 to 24 hours after the suspect meal is eaten and last for about 24 hours. During food poisoning, the organism is consumed, sporulation takes place in the gut, enterotoxin is released, and disease follows. Enterotoxin is the product of the *cpe* gene and the two bacterial phenotypes are designated Ent+ and Ent-. Clearly, since most Ent- strains tested can sporulate, enterotoxin production is not a prerequisite for spore formation. However, although some enterotoxin is formed during exponential growth, the transcriptional factors activating the sporulation genes also up-regulate the *cpe* gene and maximum levels of enterotoxin occur during sporulation. Food poisoning as well as the molecular genetics and pathogenesis of *C. perfringens* enterotoxin have been well reviewed by Rood and Cole (1) and Stringer (2).

Even in the absence of typical food poisoning, *C. perfringens* enterotoxin can be produced *in vivo*, resulting in diarrhea. Following disruption of the normal intestinal microflora by agents such as antibiotics, clostridia can colonize the vacated ecological niches within the bowel. (*Clostridium difficile* and *Clostridium spiroforme* both infect susceptible hosts in this way.) In one study of *C. perfringens* diarrheas (3), 60% were antibiotic-associated. Just as the symptoms of foodborne botulism differ from those seen in infants with *in vivo* intoxication, the clinical features of *C. perfringens* food poisoning and antibiotic-associated disease also differ. In the latter, diarrhea lasts on average 10 days, fecal blood and mucus are common and relapses may occur.

### **PRINCIPLE OF THE TEST**

The *Clostridium perfringens* Enterotoxin Test uses two antibodies to *C. perfringens* enterotoxin. The microassay wells supplied with the kit contain immobilized polyclonal antibody against *C. perfringens* enterotoxin. The *Conjugate* consists of polyclonal antibody conjugated to horseradish peroxidase. In the assay, an aliquot of a fecal specimen is emulsified in the *Diluent* and the diluted specimen is then transferred to the microassay well. If enterotoxin is present in the specimen, it will bind to the immobilized polyclonal antibody and conjugate during the 2-hour incubation phase. Any unbound material is removed during the washing steps. Following the addition of *Substrate*, a color is detected due to the enzyme-antibody-antigen complexes that form in the presence of enterotoxin.

### **MATERIALS PROVIDED**

**[DIL | SPE]** ***Diluent***, 40 mL (buffered protein solution + 0.02% thimerosal\*). The *Diluent* is also to be used as the negative control solution (see TEST PROCEDURE).

H411: Toxic to aquatic life with long lasting effects  
P273, P501

**[CONJ | ENZ]** ***Conjugate***, 7 mL (polyclonal antibody specific for enterotoxin coupled to horseradish peroxidase - in a buffered protein solution + 0.02% thimerosal\*).

**[SUBS | REAG]** ***Substrate***, 14 mL (buffered solution containing tetramethylbenzidine and peroxide).

**[CONTROL | +]** ***Positive Control Reagent***, 3.5 mL (enterotoxin in a buffered protein solution containing 0.02% thimerosal\*).

H412 - Harmful to aquatic life with long lasting effects.  
P273, P501

**WASHBUF 20X** **Wash Buffer Concentrate**, 50 mL (20X concentrate containing phosphate-buffered saline, detergent and 0.2% thimerosal\*).

Signal Word: Warning

H373: May cause damage to organs through prolonged or repeated exposure

H411: Toxic to aquatic life with long lasting effects

P260, P273, P314, P391, P501

**H<sub>2</sub>SO<sub>4</sub> 0.6N** **Stop Solution**, 7 mL (0.6 N sulfuric acid). CAUTION: Avoid contact with skin. Flush with water immediately if contact occurs.

Signal Word: Danger

H314: Causes severe skin burns and eye damage

P260, P264, P280, P301, P330, P331, P303, P361, P353, P363, P304, P340,

P310, P321, P305, P351, P338, P501

**MA** **PLT** **12 Assay Well Strips**, each consisting of 8 wells coated with polyclonal antibody specific for enterotoxin (stored with desiccant).

\*contains mercury

## ACCESSORIES

100 disposable plastic pipettes

2 plastic adhesive sheets

50 applicator sticks

## MATERIALS NOT PROVIDED

Squirt bottle for *Wash Solution*

950 mL distilled water for diluting *Wash Solution*

Incubator set at 37°C ± 2°C

ELISA reader capable of reading dual wavelength at 450/620 nm or 450/550 nm, or single wavelength at 450 nm (a dual wavelength plate reader is recommended; absorbances should be measure at 450 nm and referenced at 620 nm or 550 nm)

Discard container and absorbent paper (e.g., paper towels)

Refrigerator for storage of kit and specimens (set between 2° and 8°C)

Tubes for dilutions

Vortex mixer

## PRECAUTIONS

1. For *in vitro* diagnostic use. For professional use only.
2. Reagents from different kits should not be mixed. Do not use kit past the expiration date.
3. Reagents should be at room temperature before use.
4. Caps and tips are color coded; do not mix!
5. When handling assay wells, avoid scratching the bottom of the wells because this may result in elevated absorbance readings.
6. Hold dropper bottles vertically to ensure proper drop size.
7. Unused microassay wells must be placed back inside of the resealable bag with the desiccant to protect them from moisture.
8. Perform the washing procedure as directed to avoid high background reactions.
9. The *Substrate* is light sensitive and should be protected from direct sunlight or UV sources.
10. Handle specimens and used microassay wells as if capable of transmitting infectious agents. Wear disposable gloves when doing the test.
11. The *20X Wash Buffer Concentrate* contains 0.2% Thimerosal as a preservative. Once diluted to normal use concentration this solution is classified as non-hazardous. The *Stop Solution* contains 0.6N sulfuric acid. Flush with water immediately if contact occurs. Take off contaminated clothing and wash it before reuse. Handle reagents according to existing regulations for laboratory safety and good laboratory practice. Safety Data Sheets for this product are available upon request, contact technical support.
12. Follow your national, regional, and local ordinances accordingly for waste disposal regulations. Do not place in trash, dispose of as hazardous waste.

## COLLECTION AND HANDLING OF FECAL SPECIMENS

**NOTE:** Standard collection and handling procedures used in-house for fecal specimens are appropriate. Specimens should be transported and diluted in kit *Diluent* as soon as possible. Diluted specimens should be stored between 2° and 8°C. Whenever possible, test samples which are less than 24 hours old. Refrigeration of pre-diluted samples is an effective way to reduce loss of ELISA reactivity for up to 7 days; this does not work for all fecal specimens (4). Freezing and thawing of the specimen, especially multiple times, may result in loss of activity due to degradation of the toxin. Make sure that specimens are thoroughly mixed prior to performing the assay. This includes complete mixing (vortexing) of the specimen prior to transfer to *Diluent* as well as complete mixing of the diluted specimen prior to transfer to the assay well. The *Diluent* has been formulated to stabilize the toxin in fecal specimens and minimize degradation. Disposable pipettes, which are included, are graduated at 50, 100, 200, and 300 µL.

1. Set up one dilution tube for each specimen to be tested. Add 200 µL *Diluent* to each tube. *Label the tube directly on the side.*
2. For **liquid fecal specimens**, use a plastic pipette to transfer 50 µL of specimen to the tube. Make sure the liquid specimens are evenly suspended before transferring. For **formed fecal specimens (solid specimens from persons without diarrhea are normally an inappropriate sample)**, use an applicator stick to transfer the fecal specimen to the tube. Transfer an amount equal to about 3 mm in diameter with the applicator stick to the *Diluent*.
3. Vortex the tubes for 10 seconds and store between 2° and 8°C until the ELISA is performed. Vortex again before transferring diluted specimen to microassay well. This ensures thorough mixing of the specimen.

## PRELIMINARY PREPARATIONS

1. **IMPORTANT** - All reagents must be at room temperature prior to use in the assay.
2. Prepare the 1X *Wash Solution*. The *Wash Buffer Concentrate* is supplied as a 20X concentrate (a precipitate may be noticed). It should be diluted to a total volume of 1 liter by adding 50 mL of the concentrate to 950 mL of distilled water. Label the bottle. Store any unused 1X *Wash Solution* between 2° and 8°C.
3. Assay Strip Preparation. Each strip contains 8 wells coated with affinity purified polyclonal antibody specific for *C. perfringens* enterotoxin. Each specimen or control will employ one of these coated wells. Determine the number of wells to be used. Assay wells not used must be returned to the plastic bag and carefully resealed with desiccant.

## TEST PROCEDURE

1. Add one drop (50 µL) of *Conjugate* (red cap) to each well. Be sure to hold each bottle vertically when adding the drops. Use 1 well for each fecal specimen, 1 well for the *Positive Control*, and 1 well for the Negative Control (i.e., *Diluent*). Identification marks may be written directly on the side of the well.
2. Transfer 100 µL (2 drops using a disposable pipette provided) of diluted specimen to the assay well. Add 1 drop (50 µL) of the *Positive Control Reagent* (black cap) to the positive control well and 50 µL (1 drop using a disposable pipette provided) of the negative control (i.e., *Diluent*) to the negative control well.
3. Cut the adhesive plastic sheet to the size necessary to cover the wells. Cover the wells and incubate them at 37°C ± 2°C for 2 hours.
4. Shake out the contents of the assay wells into a discard pan.
5. Wash each well using the 1X *Wash Solution* in squirt bottle with a fine-tipped nozzle, directing the *Wash Solution* to the bottom of the well with force. Fill the wells, then shake the *Wash Solution* out of the well into a discard pan. Slap the inverted plate on a dry paper towel and repeat steps #4 and #5 **four times** using a dry paper towel each time. If any particulate matter is seen in the wells, continue washing until all the particulate matter is removed.

6. After washing, completely remove any residual liquid in the wells by striking the plate once again onto a dry paper towel until no liquid comes out. *Dispose of paper towels and specimen containers properly.*
7. Add 2 drops (100  $\mu\text{L}$ ) of *Substrate* (blue cap) to each well. Incubate the wells at room temperature for 15 minutes.
8. Add 1 drop (50  $\mu\text{L}$ ) of *Stop Solution* (yellow cap) to each well. Gently tap the wells and wait 2 minutes before reading. The addition of the *Stop Solution* converts the blue color to a yellow color which may be quantitated by measuring the optical density at 450 nm on a microplate ELISA reader (reference 620 or 550 nm on dual wavelength reader). The instrument should be blanked against air. Wipe the underside of each well before measuring the optical density. If an ELISA reader is unavailable, the test may be read visually in good light against a white background. If a dual reader is used, blank against air at 620 or 550 nm and read at 450 nm. Read within ten minutes after adding *Stop Solution*.

### QUALITY CONTROL

1. A positive and a negative control must be run with each series of test specimens.
2. Positive and negative controls must fall in their respective ranges or the test is not valid. If these results are out of range, please call Technical Services.
  - a) Positive Control must be a visible yellow color.  
If read on a spectrophotometer, the  $\text{OD}_{450}$  or  $\text{OD}_{450/620}$  or  $\text{OD}_{450/550}$  must be  $\geq 0.500$ .
  - b) Negative Control must be visually clear.  
If read on a spectrophotometer, the  $\text{OD}_{450}$  must be  $< 0.120$ ,  
or the  $\text{OD}_{450/620}$  or  $\text{OD}_{450/550}$  must be  $< 0.080$ .
3. Wells that are clear visually but that give an absorbance  $\geq 0.120$  should be wiped on the underside and the absorbance should be measured again.
4. Visual readings must be taken in good light against a white background.
5. A sample that yields a weak positive result (i.e.,  $< 0.200$ ) and is adjacent to a strong positive should be repeated to assure carryover did not occur.

### INTERPRETATION OF RESULTS

1. **Visual Reading**  
Negative = Colorless  
Positive = Any yellow color
2. **Spectrophotometric Single Wavelength at 450 nm**  
Negative =  $\text{OD}_{450} < 0.120$   
Positive =  $\text{OD}_{450} \geq 0.120$
3. **Spectrophotometric Dual Wavelength at 450/620 nm or 450/550 nm**  
Negative =  $\text{OD}_{450/620}$  or  $\text{OD}_{450/550} < 0.080$   
Positive =  $\text{OD}_{450/620}$  or  $\text{OD}_{450/550} \geq 0.080$

A positive test result indicates that *C. perfringens* enterotoxin is present in the fecal specimen.

### SHELF LIFE AND STORAGE

Upon arrival, the kit should be inspected to ensure that the components are not frozen or warm-to-the-touch due to improper shipping conditions. Each component in the kit should be inspected for any signs of leakage. The expiration date of the kit is given on the label. Expiration dates for each component are listed on the individual labels. The kit containing the reagents with designated shelf life should be stored between 2° and 8°C and should be returned to the refrigerator as soon as possible after use.

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© 2021 TECHLAB, Inc. All rights reserved.
The TECHLAB Logo, and TECHLAB are trademarks of TECHLAB, Inc.