# ***E. HISTOLYTICA II***<sup>TM</sup>

A 2nd generation Monoclonal ELISA for detecting  
*E. histolytica* adhesin in fecal specimens

Catalog No. 30404 (96 Tests)

U.S. Patent #5,272,058

**IVD** *In Vitro* Diagnostic Medical Device

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## **E. HISTOLYTICA II™**

### **INTENDED USE**

The *E. HISTOLYTICA II™* test is an enzyme immunoassay for the rapid detection of the adhesin of *E. histolytica* in human fecal specimens. It is indicated for use with fecal specimens from patients with diarrhea or dysentery to determine the presence of *E. histolytica* gastrointestinal infection. The test can be used for fecal specimens submitted for routine clinical testing from adults or children. Conventional microscopy is not a prerequisite for use of the test.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician.

### **SUMMARY AND EXPLANATION**

*Entamoeba histolytica* (formerly referred to as pathogenic *E. histolytica*) and *E. dispar* (formerly referred to as nonpathogenic *E. histolytica*) are intestinal parasites that infect approximately a half billion people worldwide annually (1,2). Only malaria and schistosomiasis are believed to be more prevalent parasitic causes of morbidity and mortality. Of the huge number of persons infected, most are infected with *E. dispar*, which has not been associated with disease. Infection with *E. dispar* rather than *E. histolytica* is believed to explain, at least in part, the low rate of disease considering the high rate of infection. Even so, it is estimated that approximately 10% of the half billion people infected each year are infected with *E. histolytica*. These individuals become symptomatic and develop colitis and liver abscess, resulting in a mortality rate estimated between 40,000 and 120,000 persons annually. It is important to distinguish between these two species because *E. dispar* is not associated with colitis or liver abscess. In addition, inaccurate diagnosis may result in unwarranted and unnecessary treatment with harsh drugs with side effects.

Patients infected with pathogenic *E. histolytica* may show a wide range of conditions. Some are completely asymptomatic (1-7). These persons shed millions of cysts daily and represent a potential reservoir for dissemination. Some patients may show mild diarrhea that develops into bloody diarrhea with abdominal cramps, eventually resulting in fulminant colitis. Because of the tissue damage that can occur, perforation of the intestine may result and the amoeba may disseminate to other parts of the body. Approximately 10% of persons with invasive amebiasis develop liver abscess. The life cycle of *E. histolytica* is similar to that observed with other amoeba. The organism exists either as a trophozoite or as a cyst. Humans serve as the primary reservoir, with the organism being spread through ingested food and contaminated water, or by venereal transmission. The cyst is very stable. Once it enters the intestine, it begins to divide into trophozoites that bind to the intestinal mucosa via a galactose or N-acetyl-D-galactosamine-binding lectin referred to as the galactose adhesin. Once the trophozoites have attached, they release tissue-damaging enzymes and proteins that lyse the mucosal cell. The galactose adhesins of *E. histolytica* and *E. dispar* cross-react serologically but contain distinct epitopes. The adhesin is antigenically conserved but monoclonal antibodies can be used to distinguish the adhesin from *E. histolytica* and *E. dispar*.

### **PRINCIPLE OF THE PROCEDURE**

The *E. HISTOLYTICA II™* test uses antibodies to the adhesin. The microassay wells contain immobilized polyclonal antibody that binds adhesin of *E. histolytica/dispar*. The *Conjugate* is a monoclonal antibody-peroxidase conjugate specific for *E. histolytica* adhesin. In the assay, an aliquot of a fecal specimen is emulsified in *Diluent* and the diluted specimen is transferred to a microassay well. If adhesin is present in the specimen, it binds to the conjugate and immobilized polyclonal antibody during the incubation phase. Any unbound material is removed during the washing steps. Following the addition of substrate, a color develops due to the enzyme-antibody-antigen complexes that form in the presence of adhesin.

**MATERIALS PROVIDED**

<b>DIL</b>   <b>SPE</b>	<b>Diluent</b> , 40 mL (buffered protein solution with 0.02% thimerosal)*. The <i>Diluent</i> also is to be used as the negative control solution (see TEST PROCEDURE).	
<b>CONJ</b>   <b>ENZ</b>	<b>Conjugate</b> , 7.0 mL (mouse monoclonal antibody specific for adhesin from <i>E. histolytica</i> ; coupled to horseradish peroxidase and in a buffered protein solution with 0.02% thimerosal)*.	
<b>SUBS</b>   <b>REAG</b>	<b>Substrate</b> , 14.0 mL (solution containing tetramethylbenzidine and peroxide).	
<b>CONTROL</b>   <b>+</b>	<b>Positive Control Reagent</b> , 3.5 mL (purified adhesin from <i>E. histolytica</i> in a buffered protein solution containing 0.02% thimerosal)*.	
<b>WASHBUF</b>   <b>20X</b>	<b>20X Wash Buffer Concentrate</b> , 50 mL (20X concentrate containing phosphate-buffered saline, detergent, and 0.2% thimerosal)*. Signal Word: Warning H373: May cause damage to organs through prolonged or repeated exposure H411: Toxic to aquatic life with long lasting effects P260, P273, P314, P391, P501	  
<b>H<sub>2</sub>SO<sub>4</sub></b>   <b>0.6N</b>	<b>Stop Solution</b> , 7.0 mL (0.6 N sulfuric acid). CAUTION: Avoid contact with skin. Flush with water immediately if contact occurs. Signal Word: Danger H314: Causes severe skin burns and eye damage P260, P264, P280, P301, P330, P331, P303, P361, P353, P363, P304, P340, P310, P321, P305, P351, P338, P501	
<b>MA</b>   <b>PLT</b>	<b>12 Assay Well Strips</b> , each consisting of 8 wells coated with polyclonal antibody specific for adhesin from <i>E. histolytica</i> (stored with desiccant).	

1 plastic resealable bag

Plastic adhesive sheet

\*contains mercury

**MATERIALS REQUIRED BUT NOT PROVIDED**

Squirt bottle for Wash Solution	Vortex mixer
Refrigerator for storage (2°-8°C)	Plastic pipettes
Plastic tubes for dilution of fecal specimens	Discard container/absorbent paper
1-L bottle for diluted Wash Solution	Disposable swabs
950 mL distilled or deionized water for diluting Wash Buffer	
ELISA reader that reads 450 nm or 450/620 nm	

**WARNINGS AND PRECAUTIONS**

1. Rx Only - Prescription Only.
2. For *in vitro* diagnostic use. For professional use only.
3. Reagents from different kits should not be mixed. Do not use a kit past the assigned expiration date.
4. Reagents should be at room temperature before use.
5. Caps and tips are color coded; do not mix!
6. When handling assay wells, avoid scratching the bottom of the wells because this may result in elevated absorbance readings.
7. Hold dropper bottles vertically to ensure proper drop size.
8. **Unused microassay wells must be replaced in the resealable pouch with the desiccant to protect them from moisture.**
9. Perform the washing procedure as directed to avoid high background reactions.
10. Use fecal specimens within 24 hours of collection for optimal results. Frozen specimens (-20°C or lower) may lose activity due to freezing and thawing.
11. Handle specimens and used microassay wells as if capable of transmitting infectious agents. Wear disposable gloves when doing the test.
12. The *20X Wash Buffer Concentrate* contains 0.2% Thimerosal as a preservative. Once diluted to normal use concentration this solution is classified as non-hazardous. The

*Stop Solution* contains 0.6N sulfuric acid. Flush with water immediately if contact occurs. Take off contaminated clothing and wash it before reuse. Handle reagents according to existing regulations for laboratory safety and good laboratory practice. Safety Data Sheets for this product are available upon request, contact technical support.

13. Follow your national, regional, and local ordinances accordingly for waste disposal regulations. Do not place in trash, dispose of as hazardous waste.

## STORAGE AND HANDLING

The expiration date of the kit is given on the label. Expiration dates for each component are listed on individual labels. The kit should be stored between 2° and 8°C and should be returned to the refrigerator as soon as possible after use.

## SPECIMEN COLLECTION AND HANDLING

**NOTE:** Standard collection and handling procedures used in-house for fecal specimens are appropriate. Do not use specimens that have been collected or stored in 10% formalin, Sodium Acetate Formalin, or Polyvinyl Alcohol fixatives. Specimens should be transported as soon as possible and stored between 2° and 8°C. If possible, test fecal specimens which are less than 24 hours old. Store specimens at -20°C, or lower, if the test cannot be performed within 48 hours. Freezing and thawing of the specimen, especially multiple times, may result in loss of activity due to degradation of the adhesin. **Make sure that specimens are thoroughly mixed (vortexed) prior to performing the assay. This includes complete mixing of the specimen prior to transfer to *Diluent* as well as complete mixing of the diluted specimen prior to transfer to the microwell.** The *Diluent* has been formulated to stabilize the adhesin in fecal specimens and minimize degradation.

1. Set up one dilution tube for each specimen to be tested. Add 400 µL *Diluent* to each tube. *Label the tube directly on the side.*
2. **Thoroughly mix (vortex) the fecal specimen to ensure adequate sampling.**
3. For **formed samples**, use a swab to transfer the fecal specimen to the tube. Coat the swab completely before transferring the specimen. Mix the swab in the *Diluent* to remove as much sample as possible and squeeze the swab against the side of the tube to express any residual liquid. This procedure results in the transfer of approximately 0.15 to 0.20 g of specimen. If the swab is not used for formed specimens, transfer approximately 0.15 to 0.20 g of specimen (about the size of a small pea) to *Diluent*. For **liquid samples**, transfer 400 µL specimen to tube. Make sure the liquid specimens are evenly suspended (vortexed) before transferring.
4. Vortex the tubes for 10 seconds and store between 2° and 8°C until the ELISA is performed. Vortex again before transferring diluted specimen to microassay well. This ensures thorough mixing of the specimen.

## TEST PROCEDURE

1. Two control wells must be used each time the test is performed. These wells serve as positive and negative controls. One well is needed for each patient sample. Additional control wells are not needed to run a greater number of samples. Add 1 drop (50 µL) of *Conjugate* (red cap) to a positive control well, negative control well, and patient sample well. Hold the *Conjugate* bottle vertically when adding the drops. Identification marks may be written directly on side of well.
2. Add 1 drop (50 µL) of the *Positive Control Reagent* (black cap) to the positive control well and 100 µL of the negative control (i.e., *Diluent*) to the negative control well. Transfer 200 µL of diluted specimen to the test well. Cut the adhesive plastic sheet to the size necessary to cover the wells. Cover wells and incubate them at room temperature for 2 hours.
3. Shake out the contents of the assay wells into a discard pan. Wash each well using the diluted *Wash Solution* in a squirt bottle with a fine-tipped nozzle, directing the

*Wash Solution* to the bottom of the well with force. Fill the wells, then shake the wash solution out of the well into a discard pan. Slap the inverted plate on a dry paper towel and repeat the washing step **four times** using a dry paper towel each time. If any particulate matter is in the wells, continue washing until all the particulate matter is removed.

**NOTE:** The *Wash Solution* is supplied as a 20X concentrate (a precipitate may be noticed). It should be diluted to a total volume of 1000 mL by adding 50 mL of the concentrate to 950 mL of distilled water. Label the bottle. Any unused *Wash Solution* should be stored between 2° and 8°C.

4. After washing, completely remove any residual liquid in the wells by striking the plate once again onto a dry paper towel until no liquid comes out. *Dispose of paper towels & specimen containers properly.*
5. Add 2 drops (100 µL) of *Substrate* (blue cap) to each well. Gently tap the wells initially and again at 5 minutes to mix the substrate. Incubate the wells at room temperature for 10 minutes.
6. Add 1 drop (50 µL) of *Stop Solution* (yellow cap) to each well. Gently tap the wells and wait 2 minutes before reading. The addition of the *Stop Solution* converts the blue color to a yellow color, which may be quantitated by measuring the optical density at 450 nm on a microplate ELISA reader. The instrument should be blanked against air. Wipe the underside of each well before measuring the optical density. If a dual reader is used, blank against air at 620 nm and read at 450 nm. Read within ten minutes after adding *Stop Solution*.

## QUALITY CONTROL

A positive and a negative control must be run with each series of test specimens. Test results should be recorded within 10 minutes of adding the *Stop Solution*. Each positive control well should be an easily visible yellow color and should give an absorbance of 0.500 or higher after the absorbance value of the negative control has been subtracted (see **INTERPRETATION OF RESULTS**). Negative control wells must have an absorbance value of <0.150. Test results are not valid unless the performance characteristics of the positive and negative controls are met. If these results are not observed, call Technical Services.

Any wells that are clear visually but that give a positive reaction (see **INTERPRETATION OF RESULTS**) should be wiped on the underside and the absorbance should be measured again. Each component in the kit should be inspected for any signs of leakage. Upon arrival, the kit should be inspected to ensure that components are not frozen or warm-to-the-touch due to improper shipping conditions. Test results along with control absorbance values should be recorded and reported according to in-house procedures and should be stored according to in-house procedures for future reference.

## INTERPRETATION OF RESULTS

### Visual Interpretation

The negative control well should be colorless or have only a faint yellow color. The positive control well should give an obvious yellow color. If these results are not observed, call Technical Services. A test sample is considered positive if it has an obvious yellow color when compared to the negative control well. It may be less yellow or more yellow than the color observed in the positive control well. A test sample is considered negative if the reaction is colorless or less yellow than the negative control well.

### Spectrophotometric Interpretation

1. Set the microplate ELISA reader to read at 450 nm. If a dual wavelength reader is used, set the ELISA reader to read at 450 nm and reference at 620 nm.
2. Determine the absorbance value for the negative control. The negative control reading should be 0.150 or less. If not, the test is not valid and should be repeated, paying attention to the wash procedure.

3. Subtract the reading of the negative control well from the readings of the positive control well and test wells before interpreting the results (the microplate ELISA reader may be set to blank on the negative control well so this step is performed automatically).
4. The reading for the positive control should be 0.500 or higher after the negative control reading has been subtracted. If this value is not obtained, the test should be repeated.
5. A specimen is considered positive for adhesin if the reading is 0.050 or higher after the negative control reading has been subtracted. A specimen is considered negative for the adhesin if the reading is <0.050.
6. A positive test result indicates that *E. histolytica* adhesin is present in the fecal specimen. A negative result indicates that *E. histolytica* adhesin is not present in the fecal specimen.

### EXPECTED VALUES

It is estimated that *Entamoeba histolytica* infects about 500 million people around the world (7,8). Roughly 90% of these persons remain asymptomatic, whereas about 10% develop clinical symptoms ranging from gastrointestinal disease to liver abscesses. On average, about 3,500 cases of amebiasis per year are reported to the Centers for Disease Control in the United States (9). High risk groups include persons who have traveled abroad, immigrants, immunocompromised persons, migrant workers, and active male homosexuals (8). Nonpathogenic strains (*E. dispar*) are predominate among male homosexuals (10). The disease often is transmitted by asymptomatic carriers of *E. histolytica*.

### PERFORMANCE CHARACTERISTICS

In one study, a total of 62 specimens that had been tested by culture/zymodeme and frozen were evaluated in the *E. HISTOLYTICA II*<sup>TM</sup> test. The *E. HISTOLYTICA II*<sup>TM</sup> test detected 31 of 32 specimens that were confirmed positives by culture/zymodeme. The remaining 30 specimens were negative both by culture/zymodeme and by the *E. HISTOLYTICA II*<sup>TM</sup> test, giving a sensitivity and specificity of 96.9% and 100%, respectively. The predictive positive and negative values were 100% and 96.8%, respectively, and the correlation was 98.4%.

In a clinical evaluation performed in Mirpur, Bangladesh, a total of 757 specimens were tested with the *E. HISTOLYTICA II*<sup>TM</sup> test and the results were compared with culture results for *Entamoeba*. A total of 31 specimens were culture-positive. Zymodeme analysis identified 12 of these as positive for *E. histolytica* and 19 as *E. dispar*. The *E. HISTOLYTICA II*<sup>TM</sup> test was positive for all 12 *E. histolytica*-positive specimens and negative for 18 of the 19 *E. dispar*-positive specimens. There were 701 fecal specimens that were negative by culture/zymodeme and by the *E. HISTOLYTICA II*<sup>TM</sup> test. There were 25 specimens that were culture-negative but positive for fecal antigen with the *E. HISTOLYTICA II*<sup>TM</sup> test. Compared with culture/zymodeme analysis, the *E. HISTOLYTICA II*<sup>TM</sup> test exhibited a sensitivity and specificity of 100% and 94.7%. The predictive positive and negative values were 92.3% and 100%, and the correlation was 96.8%.

In a third study, a total of 55 specimens were tested in-house at TECHLAB. There were 47 specimens that were spiked with *E. histolytica* trophozoites and 8 that were not spiked. The specimens were assayed with the *E. histolytica* TEST, which is specific for *E. histolytica*, and with the *E. HISTOLYTICA II*<sup>TM</sup> test. The 47 specimens were positive in both tests and the 8 specimens were negative in both tests. The *E. HISTOLYTICA II*<sup>TM</sup> test exhibited higher A<sub>450</sub> readings than the *E. histolytica* TEST. Spectrophotometric readings from both tests correlated 100% with visual readings.

## ANALYTICAL SENSITIVITY

The *E. HISTOLYTICA II*<sup>TM</sup> test detects approximately 0.2 to 0.4 ng per well of adhesin from *E. histolytica*.

## REPRODUCIBILITY

The inter-assay variation of the *E. HISTOLYTICA II*<sup>TM</sup> test was determined by evaluating three positive specimens and three negative specimens. Each specimen was tested each day for five consecutive days. The average inter-assay coefficient of variation for positive specimens was 13.811%. The intra-assay variation was determined by analyzing eight different specimens consisting of four positive (PZ, pathogenic zymodeme) and four negative (NPZ, nonpathogenic zymodeme) specimens. Each specimen was tested in 12 wells. The average intra-assay coefficient of variation for positive specimens was 4.728%.

## CROSS-REACTIVITY

The *E. HISTOLYTICA II*<sup>TM</sup> test was evaluated using fecal specimens found to be positive for a variety of fecal parasites and bacteria. No cross-reactivity was observed with fecal specimens that contained any of the following enteric pathogens: *Ascaris lumbricoides*, *Blastocystis hominis*, *Clostridium difficile*, *Cryptosporidium* species, *Endolimax* species, *Entamoeba coli*, *Escherichia coli*, *Giardia lamblia*, Rotavirus, *Salmonella typhimurium*, *Shigella sonnei*, *Trichuris trichura*.

## LIMITATIONS OF THE *E. HISTOLYTICA II*<sup>TM</sup>

1. The *E. HISTOLYTICA II*<sup>TM</sup> test detects adhesin in fecal specimens. The test confirms the presence of the organism. This information should be considered by the physician along with the clinical history of the patient.
2. Fecal specimens are complex. Optimal results are obtained with specimens <24 hours old. If specimens are not assayed within this time period, they may be frozen and thawed. However, freezing and thawing, especially multiple times, may cause specimens to lose their activity due to degradation of adhesin.
3. Some specimens may give weak reactions that are inconclusive. This may be due to a number of factors such as the presence of binding substances or inactivating enzymes in the feces. Under these conditions, the specimen should be retested or a fresh specimen should be tested.
4. Additional tests that may be used in conjunction with the *E. HISTOLYTICA II*<sup>TM</sup> test include visual observation of organisms in preserved specimens. These additional tests, however, do not distinguish *E. histolytica* from *E. dispar*. The *E. HISTOLYTICA II*<sup>TM</sup> test or zymodeme analysis, which requires culturing and electrophoretic analysis of the isolate, is specific for *E. histolytica*.
5. The magnitude of the absorbance value does not correlate with organism load.
6. The *E. HISTOLYTICA II*<sup>TM</sup> test has only been validated in a culture positive population.

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For Informational Use Only

**BIBLIOGRAPHY**

1. Gathiram, V., and T. F. H. G. Jackson. 1987. A longitudinal study of asymptomatic carriers of pathogenic zymodemes of *Entamoeba histolytica*. South Afr. Med. J. 72:669-672.
2. Haque, R., K. Kress, S. Wood, T. Jackson, D. Lyerly, T. Wilkins, and W. A. Petri, Jr. 1993. Diagnosis of pathogenic *Entamoeba histolytica* infection using a stool ELISA based on monoclonal antibodies to the galactose-specific adhesin. J. Infect. Dis. 167:247-249.
3. Irušen, E. M., T. F. H. G. Jackson, and A. E. Simjee. 1992. Asymptomatic intestinal colonization by pathogenic *Entamoeba histolytica* in amebic liver abscess: prevalence, response to therapy, and pathogenic potential. Clin. Infect. Dis. 14:889-93.
4. Petri, W. A., Jr., T. F. H. G. Jackson, V. Gathiram, K. Kress, L.D. Saffer, T. L. Snodgrass, M. D. Chapman, Z. Keren, and D. Mirelman. 1990. Pathogenic and nonpathogenic strains of *Entamoeba histolytica* can be differentiated by monoclonal antibodies to the galactose-specific adherence lectin. Infect. Immun. 58:1802-1806.
5. Sargeant, P. G., J. E. Williams, and J. D. Grene. 1978. The differentiation of invasive and non-invasive *Entamoeba histolytica* by isoenzyme electrophoresis. Trans. R. Soc. Trop. Med. Hyg. 72:519-521.
6. Sargeant, P. G., and J. E. Williams. 1979. Electrophoretic isoenzyme patterns of the pathogenic and non-pathogenic intestinal amebae of man. Trans. R. Soc. Trop. Med. Hyg. 73:225-227.
7. Walsh, J. A. 1986. Problems in recognition and diagnosis of amebiasis: estimation of the global magnitude of morbidity and mortality. Rev. Infect. Dis. 8:228-238.
8. Bruckner, D. A. 1992. Amebiasis. Clin. Microbiol. Rev. 5:356-369.
9. Krogstad, D. J., H. C. Spencer, G. R. Healy, N. N. Gleason, D. J. Sexton, and C. A. Herron. 1978. Amebiasis: epidemiologic studies in the United States, 1971-1974. Ann. Intern. Med. 88:89-97.
10. Tannich, E., and G. D. Burchard. 1991. Differentiation of pathogenic from nonpathogenic *Entamoeba histolytica* by restriction fragment analysis of a single gene amplified in vitro. J. Clin. Microbiol. 29:250-255.

**Technical Support**

Further information can be obtained from contacting TECHLAB® Technical Support::

US	+1 800 TECHLAB
Phone	(540) 953-1664
Fax	(540) 953-1665
Email	ts@techlab.com

© 2021 TECHLAB, Inc. All rights reserved.

E. HISTOLYTICA II, the TECHLAB Logo, and TECHLAB are trademarks of TECHLAB, Inc. All other trademarks referenced are trademarks of their respective owners.