ASCA-CHEK™

An ELISA for the Detection of Anti-*Saccharomyces cerevisiae*
Antibodies (ASCA) as an aid in the diagnosis of Crohn's Disease
Catalog No. T5016 (96 Tests)

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ASCA-CHEK™

INTENDED USE

The ASCA-CHEK™ test is an enzyme-linked immunosorbent assay (ELISA) for the qualitative detection of human anti-*S. cerevisiae* antibodies (ASCA) in feces and serum. The test result is used as an aid in the diagnosis of Crohn's disease in combination with clinical and other laboratory findings. **FOR IN VITRO DIAGNOSTIC USE.**

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician.

EXPLANATION

An estimated 1 million Americans suffer from inflammatory bowel disease (IBD) that is comprised of Crohn's disease and ulcerative colitis (1). IBD is characterized by a chronic inflammatory response that results in histologic damage to the intestinal lining. Crohn's disease may involve the entire gastrointestinal tract and include inflammation extending into the transmural mucosa, whereas ulcerative colitis only affects the large bowel and includes inflammation of the mucosal lining. These two distinct diseases require a rapid differential diagnosis for optimal treatment. Conventional methods utilizing multiple endoscopy examinations and histological analysis may take years to confirm a diagnosis (2). Test methods for determining the presence of serum ASCA as a marker of Crohn's disease are currently available, with sensitivities ranging from 35 - 70% (3-6).

The ASCA-CHEK™ test is an ELISA for detecting ASCA in human feces and serum. The test provides a method utilizing antigens of *Saccharomyces cerevisiae* for measuring total fecal and serum ASCA as an aid to distinguish Crohn's disease from other gastrointestinal illnesses such as ulcerative colitis and irritable bowel syndrome.

PRINCIPLE OF THE TEST

The ASCA-CHEK™ test detects anti-*Saccharomyces cerevisiae* antibodies (ASCA). The microassay wells supplied with the kit contain immobilized antigens of *Saccharomyces cerevisiae*. The *Conjugate* consists of anti-human immunoglobulin antibody conjugated to horseradish peroxidase. In the assay, an aliquot of either fecal or serum specimen is emulsified in the *Diluent* and the diluted specimen is transferred to the microassay well. If ASCA are present in the specimen, they will bind to the immobilized antigens. After incubation, the wells are washed and the *Conjugate* is added. The *Conjugate* binds to the ASCA captured by the immobilized antigens. A second series of wash steps removes any unbound material. Following the addition of *Substrate*, a color is detected due to the enzyme-antibody-antigen complexes that form in the presence of ASCA.

REAGENTS

DIL 10X *10X Diluent*, 40 mL (10X concentrate of a phosphate buffered protein solution containing 0.2% thimerosal*). The 1X *Diluent* is also to be used as the negative control (see TEST PROCEDURE).

Signal Word: Warning

H373: May cause damage to organs through prolonged or repeated exposure

H411: Toxic to aquatic life with long lasting effects

P260, P273, P314, P391, P501

CONJ ENZ *Conjugate*, 7 mL (goat anti-human polyclonal immunoglobulin-HRP in a buffered solution containing 0.02% thimerosal*)

SUBS REAG *Substrate*, 14 mL (tetramethylbenzidine substrate and peroxide)

CONTROL + *Positive Control*, 3.5 mL (human ASCA in a buffered protein solution containing 0.02% thimerosal*)

WASHBUF 20X *Wash Buffer Concentrate*, 50 mL (20X concentrate containing phosphate-buffered saline, detergent and 0.2% thimerosal*)

Signal Word: Warning

H373: May cause damage to organs through prolonged or repeated exposure

H411: Toxic to aquatic life with long lasting effects

P260, P273, P314, P391, P501

H₂SO₄ 0.6N

Stop Solution, 7 mL (0.6 N sulfuric acid). CAUTION: Avoid contact with skin. Flush with water immediately if contact occurs.

Signal Word: Danger

H314: Causes severe skin burns and eye damage

P260, P264, P280, P301, P330, P331, P303, P361, P353, P363, P304, P340, P310, P321, P305, P351, P338, P501

MA PLT

Microassay Plate, 12 strips, 8 wells per strip, coated with antigens of *Saccharomyces cerevisiae* (stored with desiccant)

*contains mercury

Materials provided

2 Plastic adhesive sheets

100 Flared Transfer pipettes (flared section = 50 µL)

Materials and equipment required but not provided

Squirt bottle for 1X *Wash Solution*

Vortex mixer

Tubes for dilution of specimen

Discard container/absorbent paper

Refrigerator for storage

Deionized or distilled water

Bottle for 1X *Diluent*

37°C incubator

ELISA reader capable of reading 450 nm or 450/620 nm

PRECAUTIONS

1. Rx Only - Prescription Only
2. For professional use only.
3. Reagents from different kits should not be mixed. Do not use the kit past the expiration date.
4. Do not freeze the reagents. Store the kit between 2° and 8°C.
5. Reagents should be at room temperature before use.
6. Gently mix all reagents before use.
7. Caps and tips are color coded; do not mix!
8. When handling the microassay wells, avoid scratching the bottom of the wells because this may result in elevated absorbance readings.
9. Hold dropper bottles vertically when dispensing reagent to ensure proper drop size.
10. Unused microassay wells must be placed inside the resealable pouch immediately with the desiccant to protect them from moisture.
11. Perform the washing procedure as directed to avoid high background reactions.
12. Use fecal and serum specimens within 48 hours and 7 days, respectively, of collection to obtain optimal results. Even though a single freeze-thaw is acceptable, frozen specimens (-20°C or lower) may lose activity due to multiple freezing and thawing.
13. The *Substrate* is light sensitive and should be protected from direct sunlight or UV sources. If the *Substrate* is exposed to light and develops a color, it must be discarded.
14. Microbial contamination of reagents may decrease the accuracy of the assay. Avoid microbial contamination of reagents by using sterile disposable pipettes if removing aliquots from reagent bottles.
15. Optimal results are obtained by following the specified test procedure. The concentrations, incubation conditions, and processing specifications have been optimized for sensitivity and specificity. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test.
16. The *Positive Control* contains ASCA that is a human derived material. Material has been tested and found negative for antibody to HIV-1, HIV-2, HCV, and HbsAg. No known test method can offer complete assurance that infectious agents are absent. ALL HUMAN SOURCE PRODUCTS SHOULD BE HANDLED AS POTENTIALLY INFECTIOUS MATERIAL. A procedure for handling biohazards is published in the CDC/NIH Manual of Biosafety in Microbiology & Biomedical Laboratories.
17. Handle specimens and used microassay wells as if capable of transmitting infectious agents. Wear disposable gloves when doing the test.
18. The 10X *Diluent* and the 20X *Wash Buffer Concentrate* contain 0.2% thimerosal as a

preservative. Once diluted to normal use concentration these solutions are classified as non-hazardous. The *Stop Solution* contains 0.6 N sulfuric acid. Flush with water immediately if contact occurs. Take off contaminated clothing and wash it before reuse. Handle reagents according to existing regulations for laboratory safety and good laboratory practice. Safety Data Sheets for this product are available upon request, contact technical support.

19. Follow your national, regional, and local ordinances accordingly for waste disposal regulations. Do not place in trash, dispose of as hazardous waste.

PRELIMINARY PREPARATIONS

1. All reagents must be at room temperature prior to use in the assay.
2. **Prepare 1X Wash Solution.** The *Wash Solution* is supplied as a 20X concentrate (a precipitate may be noticed and is acceptable). Dilute to a total volume of 1 liter by adding 50 mL of the concentrate to 950 mL of deionized water. Label the bottle. Store any unused 1X *Wash Solution* between 2° and 8°C.
3. **Prepare 1X Diluent.** The *Diluent* is supplied as a 10X concentrate (a precipitate may be noticed and is acceptable). Dilute to a total volume of 400 mL by adding 40 mL of the concentrate to 360 mL of deionized water. Label the bottle. Store any unused 1X *Diluent* between 2° and 8°C.
4. **Microassay Plate preparation.** Each strip contains 8 microassay wells coated with antigens of *Saccharomyces cerevisiae*. Each specimen or control will require one of these coated wells. Avoid contact with the bottom of the wells because this is the optical window for ELISA readers. Microassay wells not used must be immediately returned to the plastic bag and carefully resealed with desiccant.
5. All reagents, with the exception of the *Wash Buffer Concentrate* and the *10X Diluent*, are supplied in ready-to-use bottles. Reagents can be dispensed directly from the dropper bottles or decanted for use with multichannel pipettes. If excess reagent has been decanted, the excess should be discarded. Do not pour back into the bottle. The *Substrate* should be stored in and used from the light protected bottle in which it is supplied. If an aliquot is removed from the original bottle for any reason, do not return unused *Substrate* to the original bottle.

COLLECTION OF **FECAL** SPECIMENS AND PREPARATION OF DILUTIONS

NOTE: Standard collection and handling procedures used in-house for fecal specimens for culture are appropriate. Do not use specimens that have been collected or stored in 10% Formalin, Merthiolate Formalin, Sodium Acetate Formalin, or Polyvinyl Alcohol. Store specimens at -20°C, or lower, if the test cannot be performed within 48 hours of collection. Thaw specimens at room temperature and mix well prior to sample preparation. Store diluted specimens between 2° and 8°C up to 48 hours. Whenever possible, test fecal specimens that are less than 48 hours old. **Mix (vortex) specimens thoroughly prior to performing the assay.** This includes complete mixing of the specimen prior to transfer to 1X *Diluent* as well as complete mixing of the diluted specimen prior to performing the assay.

1. Prepare Diluted Specimens:

For Liquid Fecal Specimens: Set up a plastic tube for each specimen to be tested. For each specimen, add 450 µL of 1X *Diluent* to each tube. Use a transfer pipette to add 50 µL (flared section) of liquid fecal specimen to the tube containing 1X *Diluent* for a 1:10 dilution and mix well using a vortex mixer.

For Formed/Solid Fecal Specimens: Set up a plastic tube for each specimen to be tested. For each specimen, add 450 µL of 1X *Diluent* to each tube. Add 0.05 g (flared section) or weigh 0.05 g of solid fecal specimen to the tube containing 1X *Diluent* for a 1:10 dilution and mix well using a vortex mixer.

2. Vortex the tubes for 10 seconds and store between 2° and 8°C until the test is performed. Vortex again before transferring diluted specimen to the microassay well. This ensures thorough mixing of the specimen.

COLLECTION OF SERUM SPECIMENS AND PREPARATION OF DILUTIONS

NOTE: Blood samples should be collected by venipuncture and allowed to clot naturally. Serum should be separated from the clot. Store undiluted specimens at -20°C , or lower, if the test cannot be performed within 7 days of collection. Thaw specimens at room temperature and mix well prior to sample preparation. Diluted specimens may be stored between 2° and 8°C for up to 7 days. Whenever possible, test serum specimens that are less than 7 days old. **Mix (vortex) specimens thoroughly prior to performing the assay.** This includes complete mixing of the specimen prior to transfer to 1X *Diluent* as well as complete mixing of the diluted specimen prior to performing the assay.

1. Prepare Diluted Specimens:

For Serum Specimens: Set up and label 3 plastic tubes for each serum specimen to be tested. For each specimen, add 450 μL of 1X *Diluent* to each of the 3 tubes.

Use a transfer pipette to add 50 μL (flared section) of serum specimen to the first tube containing 1X *Diluent* for a 1:10 dilution and mix well using a vortex mixer.

2. Next, transfer 50 μL (flared section) from the 1:10 dilution into the second tube containing *Diluent* for a 1:100 dilution and mix well using a vortex mixer.


3. Next, transfer 50 μL from the tube containing the 1:100 dilution into the third tube containing *Diluent* for a 1:1000 dilution and mix well using a vortex mixer.

Vortex the tubes for 10 seconds and store between 2° and 8°C until the test is performed. Vortex again before transferring diluted specimen to the microassay well. This ensures thorough mixing of the specimen.

TEST PROCEDURE

1. Determine the number of wells to be used. Add 1 drop of *Positive Control* (black cap) to a positive control well. Add 50 μL (flared section) of 1X *Diluent* to a negative control well. Add 100 μL (first mark past flared section) of diluted specimen to one well.

Transfer Pipette: Bulb 300 μL 200 μL 100 μL 50 μL Flared tip

2. Incubate the wells at $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$ for 30 minutes.

3. Shake out the contents of the microassay wells into a discard pan.

4. Wash each well using the 1X *Wash Solution* in a squirt bottle with a fine-tipped nozzle, directing the 1X *Wash Solution* to the bottom of the well with force. Fill the wells (approximately 400 μL when filled), then shake the 1X *Wash Solution* out of the well into a discard pan. Slap the inverted plate on a dry paper towel and repeat wash steps **four times** using a dry paper towel each time. If any particulate matter is seen in the wells, continue washing until all the matter is removed.

5. Add 1 drop of *Conjugate* (red cap) to each well. Incubate the wells at $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$ for 30 minutes.

6. Repeat Steps #3 and #4. Dispose of all paper towels and specimen containers properly.

7. Add 2 drops of *Substrate* (blue cap) to each well. Gently tap the wells to mix the contents. Incubate the wells at room temperature for 15 minutes. Gently tap the wells 1 or 2 times during the incubation period.

8. Add 1 drop of *Stop Solution* (yellow cap) to each well. Gently tap the wells and wait 2 minutes before reading. The addition of *Stop Solution* converts the blue color to a yellow color, which is detected by measuring the optical density at 450 nm on an ELISA reader. Wipe the underside of each well to remove moisture before measuring the optical density. If a dual reader is used, read at 450 nm and reference 620 nm. Read within two to ten minutes after adding *Stop Solution*.

9. Record absorbance values for the positive control, negative control, and each specimen tested.

QUALITY CONTROL

The positive and negative control must meet the following criteria or the test is invalid. The positive control well must be a visible yellow color. When read on a spectrophotometer, it must have an OD₄₅₀ and OD_{450/620} >0.500. The negative control must have an OD₄₅₀ <0.110 or an OD_{450/620} <0.080. To ensure that carryover did not occur, repeat testing if a sample gives a weak positive result (i.e., <0.400) and is adjacent to a strong positive well.

INTERPRETATION OF TEST RESULTS

Measurements should be determined at 450 nm on a single wavelength spectrophotometer or at 450/620 nm on a dual wavelength spectrophotometer.

	Spectrophotometric Reading	
	Single Wavelength at 450 nm	Dual Wavelength at 450/620 nm
Negative for fecal specimens	OD < 0.150	<0.110
Positive for fecal specimens	OD ≥ 0.150	≥0.110
Negative for serum specimens	OD < 0.110	<0.080
Positive for serum specimens	OD ≥ 0.110	≥0.080

A positive test result indicates the specimen contains ASCA.

A negative test result indicates the specimen does not contain detectable levels of ASCA.

EXPECTED VALUES

The prevalence of Crohn's disease is estimated from 10 to 198 people/100,000 and appears to be increasing. Crohn's disease is most common in northern Europe and North American populations (9). It affects men and women approximately equally.

In 2 separate groups of healthy subjects, a total of 98% (57/58) and 96% (90/94) tested negative for fecal and serum ASCA, respectively. There were 2 specimen panels consisting of feces (211) and serum (186) collected from non-Crohn's disease patients and healthy subjects from Virginia, Indiana, Minnesota, Massachusetts, and Germany that were tested with the ASCA-CHEK™ test. The age of the population ranged from 1 to 79 years old. For the 5 regional populations, the specificity for fecal specimens was 90% with a 95% Confidence Interval of 85 – 94%. For the 5 regional populations, the specificity for serum specimens was 93% with a 95% Confidence Interval of 88 – 96%.

SHELF-LIFE AND STORAGE

The expiration date of the kit is given on the outside label. Expiration dates for each component are listed on the individual labels. The kit containing the reagents with designated shelf life should be stored between 2° and 8°C and should be returned to the refrigerator as soon as possible after use.

PERFORMANCE CHARACTERISTICS

In the clinical evaluation of the ASCA-CHEK™ test using fecal specimens, there were 250 patients with IBD, 45 IBS/non-Crohn's disease patients and 58 healthy control persons enrolled at 5 separate study sites (4 clinical sites and 1 in-house site). In the IBD group, there were 142 patients with Crohn's disease and 108 with ulcerative colitis. The male to female ratio was approximately 1:1 and a total of 40% of patients was less than 18 years old. Fecal consistency ranged from solid to liquid and all specimens showed acceptable reactivity in the assay. A total of 57% of the patients with Crohn's disease were positive by the ASCA-CHEK™ test. The combined analysis for all 5 sites showed a sensitivity and specificity of 57% and 91%, respectively, with an overall correlation of 77% to clinically confirmed Crohn's disease. Table 1 shows the test results for the combined analysis of all 4 study sites and the single in-house site.

Table 1. ASCA-CHEK™ test results of 4 clinical study sites and 1 in-house site for distinguishing Crohn's disease from non-Crohn's disease using feces.

Feces N = 353	Crohn's disease	Ulcerative colitis, IBS and Non-IBD
ASCA-CHEK™ test Positive	81	20
ASCA-CHEK™ test Negative	61	191

		95% Confidence Intervals
Sensitivity	57%	49 - 65%
Specificity	91%	86 - 94%
Correlation	77%	73 - 81%

In a separate clinical evaluation of the ASCA-CHEK™ test using serum specimens, there were 218 patients with IBD, 39 IBS/non-Crohn's disease patients and 94 healthy control persons enrolled at 5 separate study sites (4 clinical sites and 1 in-house site). In the IBD group, there were 172 patients with Crohn's disease and 46 with ulcerative colitis. The male to female ratio was approximately 1:1 and a total of 40% of patients was less than 19 years old. A total of 62% of the patients with Crohn's disease was positive by the ASCA-CHEK™ test. The combined analysis for all 5 sites showed a sensitivity and specificity of 62% and 93%, respectively with an overall correlation of 78% to clinically confirmed Crohn's disease. Tables 2 and 3 show results with the ASCA-CHEK™ test and a commercial ASCA test compared to clinical assessments.

Table 2. ASCA-CHEK™ test results of 4 clinical study sites and 1 in-house site for distinguishing Crohn's disease from non-Crohn's disease using serum.

Serum N = 351	Crohn's disease	Ulcerative colitis, IBS and Non-IBD
ASCA-CHEK™ test Positive	107	13
ASCA-CHEK™ test Negative	65	166

		95% Confidence Intervals
Sensitivity	62%	55 - 69%
Specificity	93%	88 - 96%
Correlation	78%	73 - 82%

Table 3. Commercial ASCA test results of 4 clinical study sites and 1 in-house site for distinguishing Crohn's disease from non-Crohn's disease using serum.

Serum N = 274	Crohn's disease	Ulcerative colitis, IBS and Non-IBD
Commercial ASCA test Positive	89	7
Commercial ASCA test Negative	51	127

		95% Confidence Intervals
Sensitivity	64%	55 - 71%
Specificity	95%	89 - 98%
Correlation	79%	74 - 83%

In a comparative study between serum ASCA and fecal ASCA, a total of 82 IBD patients (45% adult and 55% pediatric) were enrolled and paired fecal and serum specimens were collected from each enrolled subject. The serum specimens were tested on a commercial serum ASCA IgG EIA and the fecal specimens were analyzed by the ASCA-CHEK™ test. There were a total of 5 fecal ASCA-positive specimens that were serum ASCA IgG-negative of which 4 were from patients with Crohn's disease. The comparative statistical analysis between the two assays showed a percent positive agreement of 70% and a percent negative agreement of 88% with an overall agreement of 79%. Table 4 shows the individual test results for the serum and fecal specimen comparative study.

Table 4. Comparative study done between the ASCA-CHEK™ test and a commercial serum ASCA IgG ELISA using paired fecal and serum specimens, respectively.

N = 82	Commercial Serum ASCA IgG ELISA Positive	Commercial Serum ASCA IgG ELISA Negative
Fecal ASCA-CHEK™ test Positive	28	5
Fecal ASCA-CHEK™ test Negative	12	37

		95% Confidence Intervals
Percent Positive Agreement	70%	53 - 83%
Percent Negative Agreement	88%	74 - 96%
Overall Percent Agreement	79%	70 - 86%

LIMITATIONS OF THE PROCEDURE

1. Fecal samples from the following patients should be excluded from use in the ASCA-CHEK™ test: patients with a history of HIV and/or Hepatitis B and C, patients with a history of infectious diarrhea (within 6 months), and patients having had a colostomy and/or ileostomy within 1 month.
2. The dilutions recommended for fecal and serum specimens in the Package Insert have been evaluated in clinical trials and have been found to be optimal for test results. Therefore, only the dilutions recommended in the Package Insert should be used.
3. Only fecal and serum specimens have been clinically evaluated using the ASCA-CHEK™ test.
4. Fecal specimens that have been preserved in 10% Formalin, Merthiolate Formalin, Sodium Acetate Formalin, Polyvinyl Alcohol or contain barium enema cannot be used. Serum specimens with microbial contamination or that have been heat treated should not be tested.
5. A positive ASCA-CHEK™ test indicates the presence of antibody to *S. cerevisiae* and without consideration of a patient's clinical history or physical examination does not necessarily indicate the presence of Crohn's disease.
6. Samples from patients with chronic liver disease and hypergammaglobulinemia were not evaluated by the ASCA-CHEK™ test for potential presence of fecal ASCA.
7. It is possible that the fecal and serum ASCA results may not agree.

CROSS-REACTIVITY FOR FECAL SPECIMENS

Various intestinal organisms were examined for cross-reactivity in the ASCA-CHEK™ test. For the analysis, media cultures were spiked into a ASCA-negative fecal specimen and tested in the ASCA-CHEK™ according to the Package Insert. None of the organisms below reacted in the ASCA-CHEK™ are listed in this Package Insert. The ASCA-CHEK™ test was not evaluated for potential cross-reactivity with enteric viruses.

Aeromonas hydrophila

Bacillus cereus

Bacillus subtilis

Bacteroides fragilis

Campylobacter coli

Campylobacter jejuni

Candida albicans

Clostridium butyricum

Clostridium difficile (non-toxigenic) VPI 11186

Clostridium difficile (toxigenic) VPI 10463

Clostridium difficile (ToxA-/ToxB+) 8864

Clostridium perfringens Type A

Clostridium septicum

Clostridium sordellii

Clostridium sporogenes

Enterobacter cloacae

Escherichia coli EIEC SD67

Escherichia coli ETEC E2348169

Escherichia coli O157:H7

Escherichia coli

Helicobacter pylori

Klebsiella pneumoniae

Porphyromonas asaccharolytica

Proteus vulgaris

Pseudomonas aeruginosa

Salmonella typhimurium

Serratia liquefaciens

Shigella dysenteriae

Shigella flexneri

Shigella sonnei

Staphylococcus aureus

Staphylococcus epidermidis

Streptococcus faecalis

Peptostreptococcus anaerobius

INTERFERING SUBSTANCES FOR SERUM

Potential interfering substances present in serum were tested in the ASCA-CHEK™ test. Each substance was tested for potential interference with ASCA-positive and ASCA-negative serum specimens. No interference was observed in the serum specimens tested. Table 5 shows the concentrations and methods of measurement for each substance tested.

Table 5. Interfering substances tested in the ASCA-CHEK™ test.

Test Substance	Method of analysis	Serum Concentration
Lipase	Hitachi 917	> 20 U/L
Cholesterol	Olympus AU5000	504 mg/dL
Hemolyzed blood	Visual red color	Hemolyzed
Rheumatoid Factor (RF)	Nephelometry	27.9 IU/mL
Antinuclear Antibodies (ANA)	Nephelometry	0.263 (ratio)
Bilirubin F (free)	Beckman LX20	19 mg/dL
Bilirubin C (conjugated)	Beckman LX20	18 mg/dL

INCIDENCE OF SERUM ASCA IN NON-CROHN'S DISEASE GROUPS

Sera from non-Crohn's disease groups were tested in the ASCA-CHEK™ test and showed a positivity range of 4% to 13%. Table 6 shows the individual disease groups with the corresponding percent positive rate. Serum ASCA may also be detected in patients with other autoimmune disorders (10).

Table 6. Incidence of serum ASCA in Non-Crohn's disease groups.

Non-Crohn's Disease Group	Number Tested	Number Positive	Percent Positive
Ulcerative colitis	46	6	13%
Irritable Bowel Syndrome	30	3	10%
Constipation	1	0	0%
Esophagitis	1	0	0%
<i>H. pylori</i>	1	0	0%
Alcoholic Cirrhosis	20	1	5%
ANA-positive	3	0	0%
Healthy	94	4	4%

REPRODUCIBILITY AND PRECISION

The inter-assay variation for fecal specimens was determined using 8 ASCA-positive and 4 ASCA-negative fecal specimens tested three times over a 6-day period. During the testing period all positive samples remained positive and all negative samples remained negative. The intra-assay variation was determined using 18 fecal specimens, 13 ASCA-positive and 5 ASCA-negative, tested in 8 replicates for each specimen for a single test run. During the testing period all positive samples remained positive and all negative samples remained negative.

The inter-assay variation for serum specimens was determined using 4 ASCA-positive and 4 ASCA-negative serum specimens tested eight times over a 2-day period. During the testing period all positive samples remained positive and all negative samples remained negative. The intra-assay variation for serum was determined using 4 ASCA-positive and 4 ASCA-negative serum specimens, tested in 8 replicates for each specimen for a single test run. During the testing period all positive samples remained positive and all negative samples remained negative.

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FOR INFORMATIONAL USE ONLY

REFERENCES CITED

1. Everhart, J.E. 1994. Digestive Diseases in the United States: Epidemiology and Impact. U.S. Department of Health and Human Services, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. U.S. Government Printing Office. NIH Publication no. 94-1447.
2. Sartor, R.B. 1995. Microbial Agents in Pathogenesis, Differential Diagnosis, and Complications of Inflammatory Bowel Disease. In M. Blaser, P. Smith, J. Ravdin, H. Greenberg, and R. Guerrant (ed.), Infections of the Gastrointestinal Tract. Raven Press, New York, NY.
3. Annese, V., A. Andreoli, A. Andriulli, R.D. Inca, P. Gionchetti, A. Latiano, G. Lombardi, A. Piepoli, D. Poulain, and B. Sendid. 2001. Familial Expression of Anti-*Saccharomyces cerevisiae* Mannan Antibodies in Crohn's Disease and Ulcerative Colitis: A GISC Study. *Am. J. Gastroenterology* 96:2407-2412.
4. Dubinsky, M.C., J.J. Ofman, M. Urman, S. Targan, and E.G. Seidman. 2001. Clinical Utility of Serodiagnosis Testing in Suspected Pediatric Inflammatory Bowel Disease. *Am. J. Gastroenterology* 96:758-765.
5. Peeters, M., S. Joosens, S. Vermeire, R. Vlietinck, X. Bossuyt and P. Rutgeerts. 2001. Diagnostic Value of Anti-*Saccharomyces cerevisiae* and Antineutrophil Cytoplasmic Autoantibodies in Inflammatory Bowel Disease. *Am. J. Gastroenterology* 96:730-734.
6. Sandborn, W.J., E.V. Loftus, J.F. Colombel, K.A. Fleming, F. Seibold, H.A. Homburger, B. Sendid, R.W. Chapman, W. J. Tremaine, D.K. Kaul, J. Wallace, W.S. Harmsen, A.R. Zinsmeister, and S. Targan. 2001. Evaluation of Serologic Disease Markers in a Population Based Cohort of Patients with Ulcerative Colitis and Crohn's Disease. *Inflamm. Bowel Dis.* 7:192-201.
7. Boone, J., Sandborn, W., Pelanne, L., and D. Lysterly. 2002. Measurement of Anti-*Saccharomyces Cervisiae* Antibodies in Human Feces as an Indicator of Crohn's disease. Supplement to *Am. J. Gastroenterol.* 97:S253:771. Scientific presentation at the Annual Scientific Meeting of the American College of Gastroenterology.
8. Buderus, S., Ringelmann, H., Lentze, M., Boone, J., and D. Lysterly. 2003. Fecal Anti-*Saccharomyces* Antibodies (ASCA) – Evaluation of a New Diagnostic Test in Pediatric Patients with Inflammatory Bowel Disease (IBD). Supplement to *Gastroenterology*. 124:A-196:S1327. Scientific presentation at Digestive Disease Week.
9. (www.niddk.nih.gov/health/digest/pubs/crohns/crohns.htm).
10. Muratori P, Muratori L, *et al.* Anti-*Saccharomyces cerevisiae* antibodies (ASCA) and autoimmune liver diseases. *Clin Exp Immunol* 2003; 132:473-476.

Technical Support

Further information can be obtained by contacting TECHLAB® Technical Support:

US	+ 1 800 TECHLAB
Phone	(540) 953-1664
Fax	(540) 953-1665
Email	ts@techlab.com

FOR INFORMATIONAL USE ONLY

© 2023 TECHLAB, Inc. All rights reserved.

ASCA-CHEK, the TECHLAB Logo and TECHLAB are trademarks of TECHLAB, Inc.