TECHLAB® C. DIFFICILE TOX A/B II™

An ELISA for the detection of
Clostridium difficile Toxins A and B.
Catalog No. 30397 (96 Tests) and 30397B (960 Tests)
IVD In Vitro Diagnostic Medical Device

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C. DIFFICILE TOX A/B II™

INTENDED USE

The *C. DIFFICILE TOX A/B II™* test is an enzyme immunoassay for the detection of toxins A and B produced by toxigenic strains of *Clostridium difficile*. It can be used to detect toxins A and B in fecal specimens from persons suspected of having *C. difficile* disease. The test is to be used as an aid in the diagnosis of *C. difficile* disease and results should be considered in conjunction with the patient history.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician.

EXPLANATION

After treatment with antibiotics, many patients develop gastrointestinal problems ranging from mild diarrhea to severe pseudomembranous colitis. Many cases of the milder forms of gastrointestinal illness and most cases of pseudomembranous colitis are caused by *Clostridium difficile* (1). This organism is an opportunistic anaerobic bacterium that grows in the intestine once the normal flora has been altered by the antibiotic. The disease results from the toxins that the organism produces. The clinical symptoms associated with the disease are believed to be primarily due to toxin A, which is a tissue-damaging enterotoxin (2,3). *Clostridium difficile* also produces a second toxin, designated toxin B. Toxin B, which has been referred to as the cytotoxin of the organism, is the toxin detected by the tissue culture assay currently used by many laboratories. Most strains either produce both toxins or neither toxin, although recently, toxin A negative/toxin B positive strains have been identified (4,5). These toxin A negative/toxin B positive strains test positive in the *C. DIFFICILE TOX A/B II™* test (6,7).

Clostridium difficile strains that produce high levels of toxin A also produce high levels of toxin B. Likewise, strains that produce low levels of toxin A produce low levels of toxin B, indicating that the toxin production may be regulated similarly. Tests that detect either or both toxins are being used as diagnostic aids. The genes for the toxins have been cloned and sequenced, and some properties of the toxins are now well-defined (8,9). Both toxins are large (M_r of toxin A, 308,000; M_r of toxin B, 279,000). Toxin A has a complex series of repeating units at the COOH-terminus of the molecule and these repeating units most likely serve as the binding portion that recognizes galactose-containing receptors (10). There is some evidence suggesting that toxins A and B act synergistically and that the initial tissue damage caused by toxin A allows toxin B to exert its toxicity. Therefore, although toxin A is believed to cause most of the clinical signs, toxin B also may play an important role in the disease.

PRINCIPLE OF THE TEST

The *C. DIFFICILE TOX A/B II™* test uses antibodies to *Clostridium difficile* toxins A and B. The microassay wells supplied with the kit contain immobilized affinity-purified polyclonal goat antibody against toxins A and B. The detecting antibody consists of a mixture of toxin A monoclonal mouse antibody conjugated to horseradish peroxidase and toxin B polyclonal goat antibody conjugated to horseradish peroxidase. In the assay, an aliquot of a fecal specimen is emulsified in the *Diluent* and the diluted specimen is then transferred to the microassay well containing the detecting antibody. If toxins A and B are present in the specimen, they will bind to the detecting antibody and to the immobilized polyclonal antibody during the incubation phase. Any unbound material is removed during the washing steps. Following the addition of substrate, a color is detected due to the enzyme-antibody-antigen complexes that form in the presence of toxin.

MATERIALS PROVIDED

DIL **SPE**

Diluent, 40 mL each, buffered protein solution + 0.02% thimerosal. The *Diluent* is also to be used as the negative control solution (see TEST PROCEDURE).*

H412 - Harmful to aquatic life with long lasting effects.

P273, P501

CONJ **ENZ**

Conjugate, 7 mL each, mouse monoclonal antibody specific for toxin A coupled to horseradish peroxidase and goat polyclonal antibody specific for

toxin B coupled to horseradish peroxidase in a buffered protein solution + 0.02% thimerosal*

SUBS REAG
CONTROL +

Substrate, 14 mL each, solution containing tetramethylbenzidine and peroxide
Positive Control, 3.5 mL each, inactivated toxins in a buffered protein solution containing 0.02% thimerosal*

H412 - Harmful to aquatic life with long lasting effects.

P273, P501

WASHBUF 20X

Wash Buffer Concentrate, 50 mL each, 20X concentrate containing phosphate-buffered saline, detergent, and 0.2% thimerosal*

Signal Word: Warning

H373: May cause damage to organs through prolonged or repeated exposure

H411: Toxic to aquatic life with long lasting effects

P260, P273, P314, P391, P501

H₂SO₄ 0.6N

Stop Solution, 7 mL each, 0.6N sulfuric acid. CAUTION: Avoid contact with skin; flush with water immediately if contact occurs.

Signal Word: Danger

H314: Causes severe skin burns and eye damage

P260, P264, P280, P301, P330, P331, P303, P361, P353, P363, P304, P340, P310, P321, P305, P351, P338, P501

MA PLT

Microassay Plate, 12 strips each, each strip consisting of 8 wells, coated with affinity purified goat antibodies specific for toxins A and B (stored with desiccant)

Wash Solution Label, 1 label

*contains mercury

ACCESSORIES

- 100/1000 Disposable plastic pipettes
- 2/20 Plastic adhesive sheets
- 50/500 Applicator sticks

MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

Wash bottle

Timer

Vortex mixer

Discard container

Distilled water

Tubes for dilution of specimen

Paper towels or absorbent sheets

Spectrophotometer capable of reading dual wavelength at 450/620 nm or single wavelength at 450 nm (a dual wavelength plate reader is recommended; absorbances should be measured at 450 nm and referenced at 620 nm)

Refrigerator set between 2° and 8°C

Incubator set at 37°C ± 2°C

SHELF LIFE AND STORAGE

The expiration date of the kit is given on the label. Expiration dates for each component are listed on the individual labels. The kit containing the reagents with designated shelf life should be stored between 2° and 8°C and should be returned to the refrigerator as soon as possible after use.

PRECAUTIONS

1. Rx - Prescription Only
2. For *in vitro* diagnostic use. For professional use only.
3. Reagents from different kits should not be mixed. Do not use a kit past the expiration date.
4. Reagents should be at room temperature before use.
5. Caps and tips are color coded; do not mix!
6. When handling assay wells, avoid scratching the bottom of the wells because this may result in elevated absorbance readings.
7. Hold dropper bottles vertically to ensure proper drop size.

8. *Unused microwells must be placed back inside of the resealable pouch with the desiccant to protect them from moisture.*
9. Perform the washing procedure as directed to avoid high background reactions.
10. Use fecal specimens within 24 hours of collection to obtain optimal results. Specimens that are frozen (-20°C or lower) may lose activity due to freezing and thawing.
11. The *Substrate* is light sensitive and should be protected from direct sunlight or UV sources.
12. Optimal results are obtained by following the specified test procedure. The concentrations, incubation conditions, and processing specifications have been optimized for sensitivity and specificity. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test.
13. Handle specimens and used microassay wells as if capable of transmitting infectious agents. Wear disposable gloves when doing the test.
14. The *20X Wash Buffer Concentrate* contains 0.2% Thimerosal as a preservative. Once diluted to normal use concentration this solution is classified as non-hazardous. The *Stop Solution* contains 0.6N sulfuric acid. Flush with water immediately if contact occurs. Take off contaminated clothing and wash it before reuse. Handle reagents according to existing regulations for laboratory safety and good laboratory practice. Safety Data Sheets for this product are available upon request; contact technical support.
15. Follow your national, regional, and local ordinances accordingly for waste disposal regulations. Do not place in trash, dispose of as hazardous waste.

COLLECTION AND HANDLING OF FECAL SPECIMENS

NOTE: *Standard collection and handling procedures used in-house for fecal specimens are appropriate. Specimens should be transported and diluted in the kit Diluent as soon as possible. Specimens should be stored between 2° and 8°C . Whenever possible, test samples which are less than 24 hours old. Store specimens at -20°C , or lower, if the test cannot be performed within 72 hours of collection (11). Freezing and thawing of the specimen, especially multiple times, may result in loss of activity due to degradation of the toxins. Fecal specimens that have been preserved in 10% Formalin, Merthiolate Formalin, Sodium Acetate Formalin, or Polyvinyl Alcohol cannot be used. Make sure that specimens are thoroughly mixed (vortexed) prior to performing the assay. This includes complete mixing of the specimen prior to transfer to Diluent as well as complete mixing of the diluted specimen prior to transfer to the microwell. The Diluent has been formulated to stabilize the toxins in fecal specimens and minimize degradation. Disposable pipettes which are included in the accessories are graduated at 50, 100, 200, and 300 μL .*

1. Set up one dilution tube for each specimen to be tested. Add 200 μL *Diluent* to each tube. *Label the tube directly on the side.*
2. For **formed fecal specimens**, use an applicator stick to transfer the fecal specimen to the tube. Transfer an amount equal to about 3 mm in diameter with the applicator stick to the *Diluent*. For **liquid fecal specimens**, use a plastic pipette to transfer 50 μL of specimen to tube. Make sure the liquid specimens are evenly suspended before transferring.
3. Vortex the tubes for 10 seconds and store between 2° and 8°C until the ELISA is performed (within 72 hours of collection). Vortex again before transferring diluted specimen to microassay well. This ensures thorough mixing of the specimen.
4. Semi-automated washing equipment may be used with specimens that have been centrifuged (5000 x g for 10 minutes) to remove particulate matter.

PRELIMINARY PREPARATIONS

1. **IMPORTANT:** All reagents must be at room temperature prior to use in the assay.
2. Prepare 1X *Wash Solution*. The *Wash Buffer Concentrate* is supplied as a 20X concentrate (*a precipitate may be noticed*). It should be diluted to a total volume of 1 liter by adding 50 mL of the concentrate to 950 mL of distilled water. Label the bottle. Store any unused 1X *Wash Solution* between 2° and 8°C .

3. **Assay Strip Preparation.** Each strip contains 8 wells coated with affinity purified polyclonal antibodies specific for toxins A and B. Each specimen or control will employ one of these coated wells. Determine the number of wells to be used. Avoid contact with the base of the wells. Assay wells not used must be returned to the plastic bag and carefully resealed with desiccant.

TEST PROCEDURE

1. Add 1 drop (50 μ L) of *Conjugate* (red cap) to each well. Be sure to hold each bottle vertically when adding the drops. Use 1 well for each fecal specimen, 1 well for the *Positive Control* and 1 well for the negative control (i.e., *Diluent*). Identification marks may be written directly on side of well.
2. Transfer 100 μ L (2 drops using a transfer pipette) of diluted specimen to the assay well. Add 1 drop (50 μ L) of the *Positive Control* (black cap) to the positive control well and 50 μ L (1 drop using a transfer pipette) of the negative control (i.e., *Diluent*) to the negative control well.
3. Cut the adhesive plastic sheet to the size necessary to cover the wells. Cover the wells and incubate them at $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$ for 50 minutes.
4. Shake out the contents of the assay wells into a discard pan.
5. Wash each well using the 1X *Wash Solution* in a squirt bottle with a fine-tipped nozzle, directing the *Wash Solution* to the bottom of the well with force. Fill the wells, then shake the *Wash Solution* out of the well into a discard pan. Slap the inverted plate on a dry paper towel and repeat step #5 **four times** using a dry paper towel each time. If any particulate matter is seen in the wells, continue washing until all the particulate matter is removed. **Note:** If using semi-automated or automated washing equipment, specimens must be centrifuged (5000 x g for 10 minutes) to remove any particulate matter. Add 350 μ L of 1X *Wash Solution* to each well. Wash for a total of 5 times. If any particulate matter is seen in the wells, continue washing until all the particulate matter is removed.
6. After washing, completely remove any residual liquid in the wells by striking the plate once again onto a dry paper towel until no liquid comes out. *Dispose of paper towels and specimen containers properly.*
7. Add 2 drops (100 μ L) of *Substrate* (blue cap) to each well. Gently tap the wells to mix the substrate. Incubate the wells at room temperature for 10 minutes. Gently tap the wells at 5 minutes.
8. Add 1 drop (50 μ L) of *Stop Solution* (yellow cap) to each well. Gently tap the wells and wait 2 minutes before reading. The addition of the *Stop Solution* converts the blue color to a yellow color which may be quantitated by measuring the optical density at 450 nm on a microplate ELISA reader. The instrument should be blanked against air. If a dual wavelength reader is used, blank against air at 620 and read at 450 nm. Wipe the underside of each well before measuring the optical density. If an ELISA reader is unavailable, the test may be read visually in good light against a white background. Read within ten minutes after adding *Stop Solution*.

ALTERNATE TEST PROCEDURE/RAPID FORMAT

Perform the regular test procedure according to the instructions provided above replacing the 50 minute incubation at $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$ with 20 minutes at 37°C using the *Stat Fax 2200 Incubator/Shaker* or equivalent incubator/shaker. If the *Stat Fax 2200 Incubator/Shaker* is used, set the shaker to speed 7 and the temperature to 37°C . If other shakers are used, a speed of 1500 rpm is recommended. Shaking should not cause spillage. If spillage occurs, reduce the speed accordingly (11).

QUALITY CONTROL

1. A positive and negative control must be run with each series of test specimens.
2. Positive and negative controls must fall within their respective ranges or the test is not valid.

- a. Positive Control must be a visible yellow color. If read on a spectrophotometer, the OD at 450 nm or using dual wavelength at 450/620 nm must be ≥ 0.500 .
- b. Negative Control must be visually clear. If read on a spectrophotometer, the OD at 450 nm must be < 0.120 . If read at 450/620 nm the absorbance must be < 0.080 .
3. Wells that are clear visually but give absorbance ≥ 0.120 should be wiped on the underside and remeasured.
4. Visual readings must be taken in good light against a white background.
5. A sample that yields a weak positive result (i.e., < 0.200) and is adjacent to a strong positive should be repeated to assure carryover did not occur.

INTERPRETATION OF RESULTS

1. Visual Reading

Negative = Colorless

Positive = Any yellow color

2. Spectrophotometric Single Wavelength at 450 nm

Negative = OD < 0.120

Positive = OD ≥ 0.120

3. Spectrophotometric Dual Wavelength at 450/620 nm

Negative = OD < 0.080

Positive = OD ≥ 0.080

A positive test result indicates that *C. difficile* toxin A and/or toxin B are present in the fecal specimen.

LIMITATIONS OF THE *C. DIFFICILE* TOX A/B II™ TEST

1. The *C. DIFFICILE* TOX A/B II™ test is used to detect *C. difficile* toxin in fecal specimens. The test confirms the presence of toxin in feces and this information should be taken under consideration by the physician in light of the clinical history of the patient. Inability to detect toxin A or B in fecal samples from patients suspected of having *C. difficile* disease may not preclude actual disease but may be caused by other factors (i.e., incorrect specimen collection, handling and/or storage, toxin levels lower than the kit detection limits). The *C. DIFFICILE* TOX A/B II™ test will detect Toxin A at levels ≥ 0.8 ng/mL and Toxin B at levels ≥ 2.5 ng/mL.
2. Fecal specimens represent an extremely complex clinical specimen. Optimal results with the *C. DIFFICILE* TOX A/B II™ test are obtained with specimens that are less than 24 hours old. Most specimens can be stored between 2° and 8°C for 72 hours before significant degradation of the toxin is noted. If specimens are not assayed within this time period, they may be frozen and thawed. However, freezing and thawing, especially multiple times, may cause specimens to lose their activity due to degradation of the toxin.
3. Some specimens may give weak reactions. This may be due to a number of factors such as the presence of a weakly toxigenic strain, low levels of toxin production *in vivo*, or the presence of binding substances or inactivating enzymes in the feces. *Under these conditions, the specimen should be retested or a fresh specimen should be tested.* Additional tests that may be used in conjunction with the *C. DIFFICILE* TOX A/B II™ test include isolation of the organism on selective media, latex agglutination assay for the detection of *C. difficile* (the organism), or tissue culture cytotoxicity assay for the detection of *C. difficile* toxin. The isolation of the organism does not confirm that it is toxigenic. This must be confirmed by additional testing of the isolate using either the ELISA or tissue culture assay to demonstrate toxigenicity. Likewise, the latex agglutination assay, which detects a nontoxic protein of *C. difficile*, does not demonstrate the presence of toxigenic *C. difficile*.
4. Certain toxigenic isolates of *Clostridium sordellii* produce toxins that are similar in their biologic, physicochemical and immunologic properties to the toxins of *C. difficile*. These isolates, however, have not been detected in patients with antibiotic-associated diarrhea and colitis.

5. Fecal specimens that have been preserved in 10% Formalin, Merthiolate Formalin, Sodium Acetate Formalin, or Polyvinyl Alcohol cannot be used.
6. The performance characteristics of the *C. DIFFICILE TOX A/B II™* test have not been thoroughly established in the pediatric population.

EXPECTED VALUES

The prevalence of a positive *C. DIFFICILE TOX A/B TEST* plus positive tissue culture and/or toxigenic culture was 5.4% in one study and 8.6% in another study. The prevalence will vary from location to location and hospitals may experience rates lower or higher than those observed at the sites used in the *C. DIFFICILE TOX A/B TEST* evaluation. *Clostridium difficile* disease is primarily a nosocomial disease of elderly patients, and hospitals that have higher numbers of elderly patients may experience higher rates. It is important to consider any test results in conjunction with clinical symptoms because some healthy adults and large numbers of healthy infants (up to 50%) will be positive for *C. difficile* toxin (toxin-positive by either tissue culture or ELISA). In addition, *C. difficile* carriage rates of 22% to 32% have been reported in cystic fibrosis patients (12-14).

PERFORMANCE CHARACTERISTICS

Clinical evaluation

The *C. DIFFICILE TOX A/B TEST* was compared with the tissue culture test at four U.S. hospitals and in-house at TECHLAB®, Inc. Specimens included in the evaluation were submitted to the clinical laboratory for routine testing. The tissue culture test was done according to the in-house procedure. Table 1 shows a comparison of the *C. DIFFICILE TOX A/B TEST* with tissue culture. Of the 1,152 specimens included in the evaluation, approximately 3.6% were from children ≤ 2 years of age. No *C. DIFFICILE TOX A/B TEST* (+)/tissue culture (-) specimens were identified. Overall, the sensitivity for the *C. DIFFICILE TOX A/B TEST* ranged from 83.3% to 96% (confidence interval of 87.4 to 95.0, p=0.05). The specificity was 100% in all studies. The predictive positive value was 100% in all studies and the predictive negative value ranged from 90% to 99.5% (confidence interval of 93.8 to 99.8, p=0.05). The correlation ranged from 94.9% to 99.5% (confidence interval of 96.6 to 99.4, p=0.05).

TABLE 1 Correlation of *C. DIFFICILE TOX A/B TEST* with Tissue Culture (n=1,152).

	Tissue Culture positive	Tissue Culture negative
<i>C. DIFFICILE TOX A/B TEST</i> positive	165	0
<i>C. DIFFICILE TOX A/B TEST</i> negative	14	973

Sensitivity	92.2%
Specificity	100%
Correlation	98.8%
Predictive Positive Value	100%
Predictive Negative Value	98.6%

In-house comparison

In a study performed in-house at TECHLAB, Inc., the *C. DIFFICILE TOX A/B II™* test was compared to the *C. DIFFICILE TOX A/B TEST*. In the study, 218 clinical specimens representative of those routinely submitted to the clinical laboratory for *C. difficile* testing were evaluated. The results obtained with both tests were comparable.

Centrifugation

A total of 337 fecal specimens, including 30 positives and 307 negatives, were evaluated to determine the effect of centrifugation on performance. For the analysis, specimens were diluted and vortexed as described in the package insert. The specimens were centrifuged (5,000 x g) to remove insoluble material and the supernatant fluid was assayed in the *C. DIFFICILE TOX A/B TEST*. Results were compared to results obtained with the same panel of diluted and vortexed specimens that had not been centrifuged. The results demonstrated a correlation of 100% between centrifuged and noncentrifuged specimens.

CROSS-REACTIVITY

Various organisms were examined for cross-reactivity in the *C. DIFFICILE TOX A/B TEST*. For the analysis, broth cultures mixed with *Diluent* were evaluated. Broth cultures at log phase containing $>10^8$ bacteria per mL were used. A listing of the organisms that did not react under any of the conditions is shown in Table 2. The only organisms that reacted were toxigenic *C. difficile* and a toxigenic strain of *C. sordellii* (VPI strain 9048) which produces toxins that crossreact significantly with *C. difficile* toxins. A nontoxigenic strain of *C. sordellii* that does not produce toxin HT was negative in the test. *C. sordellii* has not been implicated in pseudomembranous colitis or antibiotic-associated diarrhea. The strains of toxigenic *C. difficile* that were analyzed included six toxigenic strains that ranged from weakly to highly toxigenic and two nontoxigenic strains. The *C. DIFFICILE TOX A/B TEST* detected all six of the toxigenic strains and did not react with the two nontoxigenic strains.

TABLE 2 Organisms that do not react in the *C. DIFFICILE TOX A/B TEST*

<i>Acinetobacter lwoffii</i>	<i>Clostridium histolyticum</i>	<i>Pseudomonas aeruginosa</i>
<i>Aeromonas hydrophila</i>	<i>Clostridium novyi</i> (Types A, B, C)	<i>Salmonella choleraesuis</i>
<i>Bacillus cereus</i>	<i>Clostridium perfringens</i> (Types A-E)	<i>Salmonella enteritidis</i>
<i>Bacillus subtilis</i>	<i>Clostridium septicum</i>	<i>Salmonella typhimurium</i>
<i>Bacteroides fragilis</i> (toxigenic/nontoxigenic)	<i>Clostridium sordellii</i> (nontoxigenic)	<i>Shigella flexneri</i>
<i>Candida albicans</i>	<i>Clostridium sporiforme</i>	<i>Shigella sonnei</i>
<i>Candida krusei</i>	<i>Clostridium sporogenes</i>	<i>Staphylococcus aureus</i> (Protein A-negative)
<i>Candida tropicalis</i>	<i>Clostridium tetani</i>	<i>Staphylococcus aureus</i> (Protein A-positive @ $<10^8$)
<i>Clostridium bifermentans</i>	<i>Enterococcus faecalis</i>	<i>Streptococcus pyogenes</i>
<i>Clostridium botulinum</i> (Types A-G)	<i>Escherichia coli</i>	<i>Vibrio parahaemolyticus</i>
<i>Clostridium chauvoei</i>	<i>Escherichia coli</i> (enterohemorrhagic)	<i>Yersinia enterocolitica</i>
<i>Clostridium difficile</i> (nontoxigenic)	<i>Klebsiella pneumoniae</i>	
<i>Clostridium haemolyticum</i>	<i>Peptostreptococcus anaerobius</i>	
	<i>Proteus vulgaris</i>	

EFFECT OF FECAL SAMPLE CONSISTENCY

The reactions of fecal samples of varying consistencies in the *C. DIFFICILE TOX A/B TEST* and tissue culture assay are shown in Table 3. The rates of positive reactions were similar in all three types of fecal specimens. All of the specimens were submitted for *C. difficile* testing. The basis for submission is the clinical history of the patient and not the consistency of the specimen. In additional studies, highly purified toxin A and toxin B were used to spike liquid, semi-solid, and solid fecal specimens.

TABLE 3 Effect of sample consistency

Test	Liquid Samples	Semi-Solid Samples	Solid Samples
Number of specimens (n=435)	150	133	152
<i>C. DIFFICILE TOX A/B TEST</i> Positive	13 (8.7%)	11 (8.3%)	13 (8.6%)
Tissue Culture Positive	13 (8.7%)	14 (10.5%)	15 (9.9%)

The *C. DIFFICILE TOX A/B TEST* detected toxins A and B in liquid, semi-solid and solid fecal specimens at levels similar to those observed with toxin A and toxin B prepared in kit *Diluent*.

REPRODUCIBILITY AND PRECISION

Five fecal specimens (one negative specimen and four positive specimens) were sent to four independent laboratories for analysis using the *C. DIFFICILE TOX A/B TEST*. All specimens were kept frozen until the assay was performed. The results from each laboratory were compared with in-house results and found to be identical. The four positive specimens were confirmed positive and the negative specimen was confirmed negative at each site.

The intraassay coefficient of variation (CV) of the *C. DIFFICILE TOX A/B TEST* was determined by analyzing 32 positive and 32 negative control reactions along with 8 negative fecal specimens. Each fecal specimen was tested in 11 wells. The intraassay %CV was 7.190 with the positive control, 6.557 with the negative control, and 9.697 with the fecal specimens. The interassay CV was determined using four positive and one negative fecal specimens tested at time 0, 24, 48, and 72 hours. The %CV ranged from 9.9 to 29.6, with an average of 16.3.

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FOR INFORMATIONAL USE ONLY

REFERENCES

1. Lyerly, D.M., H.C. Krivan, and T.D. Wilkins. 1988. *Clostridium difficile*: its disease and toxins. Clin. Microbiol. Rev. **1**:1-18.
2. Lyerly, D.M., K.E. Saum, D.K. MacDonald, and T.D. Wilkins. 1985. Effects of *Clostridium difficile* toxins given intragastrically to animals. Infect. Immun. **47**:349-352.
3. Borriello, S.P., F.E. Barclay, A.R. Welch, J.M. Ketley, T.J. Mitchell, J. Stephen, and G.E. Griffin. 1985. Host and microbial determinants of the spectrum of *Clostridium difficile* mediated gastrointestinal disorders. Microecol. Ther. **15**:231-236.
4. Lyerly, D.M., N.M. Sullivan, and T.D. Wilkins. 1983. Enzyme-linked immunosorbent assay for *Clostridium difficile* toxin A. J. Clin. Microbiol. **17**:72-78.
5. Laughon, B.E., R.P. Viscidi, S.L. Gdovin, R.H. Yolken, and J.G. Bartlett. 1984. Enzyme immunoassays for detection of *Clostridium difficile* toxins A and B in fecal specimens. J. Infect. Dis. **149**:781-788.
6. Lyerly, D.M., L.A. Barroso, and T.D. Wilkins. 1992. Characterization of a toxin A-/toxin B+ isolate of *Clostridium difficile*. Infect. Immun. **60**:4633-4639.
7. Lyerly, D.M., L.M. Neville, D.T. Evans, J. Fill, S. Allen, W. Greene, R. Sautter, P. Hnatuck, D.J. Torpey, and R. Schwalbe. 1998. Multicenter Evaluation of the *Clostridium difficile* TOX A/B TEST. J. Clin. Micro. **36**:184-190.
8. Dove, C.H., S.-Z. Wang, S.B. Price, C.J. Phelps, D.M. Lyerly, T.D. Wilkins, and J.L. Johnson. 1990. Molecular characterization of the *Clostridium difficile* toxin A gene. Infect. Immun. **58**:480-488.
9. Barroso, L.A., S.-Z. Wang, C.J. Phelps, J.L. Johnson, and T.D. Wilkins. 1990. Nucleotide sequence of *Clostridium difficile* toxin B gene. Nucl. Acids Res. **18**:4004.
10. Krivan, H.C., G.F. Clark, D.F. Smith, and T.D. Wilkins. Cell surface binding site for *Clostridium difficile* enterotoxin: evidence for a glycoconjugate containing the sequence Gala1-3Gal β 1-4GlcNac. Infect. Immun. **53**:573-581.
11. Data on file.
12. Peach, S.L., S.P. Borriello, H. Gaya, F.E. Barclay, and A.R. Welch. 1986. Asymptomatic carriage of *Clostridium difficile* in patients with cystic fibrosis. J. Clin. Pathol. **39**:1013-1018.
13. Welkon, C.J., S.S. Long, C.M. Thompson, Jr., and P.H. Gilligan. 1985. *Clostridium difficile* in patients with cystic fibrosis. Amer. J. Dis. Child. **139**:805-808.
14. Wu, T.C., V.P. McCarthy, and V.J. Gill. 1983. Isolation rate and toxigenic potential of *Clostridium difficile* isolates from patients with cystic fibrosis. J. Infect. Dis. **148**:176.

Technical Support

Further information can be obtained by contacting TECHLAB® Technical Support:

US	+ 1 800 TECHLAB
Phone	(540) 953-1664
Fax	(540) 953-1665
Email	ts@techlab.com

© 2023 TECHLAB, Inc. All rights reserved.

C. DIFFICILE TOX A/B II, the TECHLAB Logo, and TECHLAB are trademarks of TECHLAB, Inc.