C. DIFFICILE TOX-B TEST

A Toxin/Antitoxin Kit for the Detection of
C. difficile Toxin B in Clinical Specimens
Catalog No. T5003 (96 Tests)

IVD *In Vitro* Diagnostic Medical Device

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C. DIFFICILE TOX-B TEST

INTENDED USE

The *C. DIFFICILE TOX-B TEST* is intended for use in conjunction with the tissue culture cytotoxicity assay for the detection of *C. difficile* toxin B in patient specimens. The test is to be used as an aid in the diagnosis of *C. difficile* disease and results should be considered in conjunction with patient history.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician.

EXPLANATION

After treatment with antibiotics, many patients develop gastrointestinal problems ranging from mild diarrhea to severe pseudomembranous colitis. Many cases of the milder forms of gastrointestinal illness and most cases of pseudomembranous colitis are caused by *Clostridium difficile* (1). This organism is an opportunistic anaerobic bacterium that grows in the intestine once the normal flora has been altered by the antibiotic. The disease results from the toxins that the organism produces. The clinical symptoms associated with the disease are believed to be primarily due to toxin A, which is a tissue-damaging enterotoxin (2,3).

Clostridium difficile also produces a second toxin, designated toxin B. Toxin B, which has been referred to as the cytotoxin of the organism, is the toxin detected by the tissue culture assay currently used by many laboratories. Most strains of *Clostridium difficile* either produce both toxins or neither toxin, although recently, toxin A negative/toxin B positive strains have been identified (4-7). Strains that produce high levels of toxin A also produce high levels of toxin B. Likewise, strains that produce low levels of toxin A produce low levels of toxin B, indicating that the toxin production may be similarly regulated. Tests that detect either toxin are being used as diagnostic aids in the identification of toxigenic strains. The genes for the toxins have been cloned and sequenced, and some properties of the toxins are now well-defined (8,9). Both toxins are large (molecular weight of toxin A, 308,000; molecular weight of toxin B, 279,000). Both toxins have a complex series of repeating units at the COOH-terminus of the molecule and these repeating units may serve as the binding portion that recognizes the receptor (10). There is some evidence suggesting that toxins A and B act synergistically and that the initial tissue damage caused by toxin A allows toxin B to exert its toxicity. Therefore, although toxin A is believed to cause most of the clinical signs, toxin B also may play an important role in the disease.

The disease caused by *Clostridium difficile* is complicated by a number of factors. Although the organism is responsible for almost all cases of pseudomembranous colitis, *C. difficile* is estimated to cause only about 25% of the cases of antibiotic-associated diarrhea. Thus, many cases of antibiotic-associated diarrhea are due to other causes. Toxigenic strains of *C. difficile* can range from weakly toxigenic to strains that produce extremely high levels of toxin. It is unclear if mild, self-limiting cases of diarrhea represent only milder forms of the disease caused by infection with weakly toxigenic strains or whether the presence of toxigenic *Clostridium difficile* occurs as a result of the diarrhea.

There are four methods currently used for the detection of *C. difficile* and its toxins in fecal specimens. These include (i) isolation of the organism, (ii) latex agglutination, (iii) tissue culture assay, and (iv) ELISA for detection of toxin A or both toxins A and B. Isolation of the organism and latex agglutination are useful for detecting the presence of the organism in fecal specimens. However, neither test demonstrates the presence of toxin. This has to be done using either the tissue culture test or an ELISA specific for the toxins. The tissue culture procedure has proved to be very useful because of the high activity of *C. difficile* toxin B in the assay. As little as one picogram of the toxin is sufficient to cause a positive reaction. Neutralization of the cytotoxic activity by *C. difficile* antitoxin confirms the presence of *C. difficile* toxin. The assay is to be used as an aid in the diagnosis of *C. difficile* disease and test results must be considered in conjunction with the patient history. This is due to the fact that some adults are asymptomatic carriers and can have detectable levels of toxin in their feces but have no clinical symptoms. In addition, up to 50% of healthy infants can have *C. difficile* toxin in their feces but show no sign of illness.

PRINCIPLE OF THE TEST

The *C. DIFFICILE TOX-B TEST* uses a tissue culture format to detect the presence of cytotoxic activity (noted by cell rounding) in fecal specimens and confirms the identification of *Clostridium difficile* toxin using specific antitoxin. Two reagents are provided in the *C. DIFFICILE TOX-B TEST*. One is the *Positive Toxin Control* reagent. The other is specific antitoxin to be used to confirm the presence of *C. difficile* toxin by neutralizing the cytotoxic activity. Also included is *Diluent* for preparing the fecal specimen. In the assay, an aliquot of a fecal specimen is emulsified in the *Diluent* and the diluted specimen is centrifuged and filtered. A sample of the filtrate is transferred to a tissue culture well containing phosphate-buffered saline (PBS) and to a well containing antitoxin. If *C. difficile* toxin is present, it will cause the cells in the well with PBS to become round, demonstrating the presence of the cytotoxic activity. The presence of *Clostridium difficile* toxin is confirmed if the cytotoxic activity is neutralized in the well containing antitoxin.

MATERIALS PROVIDED

AT	Goat Antitoxin against <i>C. difficile</i> , 3.0 mL (with preservative)
PBS	Phosphate-buffered saline (PBS) , 3.0 mL (with preservative)
CONTROL +	Positive Toxin Control , 2.0 mL active <i>C. difficile</i> toxin in a buffered protein solution containing preservative.
DIL SPE	Diluent , 100 mL (phosphate-buffered saline containing preservative and phenol red)

T5003 contains sufficient reagents for two 48-Well tissue culture plates.

MATERIALS REQUIRED BUT NOT PROVIDED

Microscope for observing cells	Incubator set at 37°C ± 2°C
Centrifuge and centrifuge tubes	Pipettes/pipettor
Test tubes for preparing dilutions	Syringe filters (0.45 to 0.8 µm)

SHELF LIFE AND STORAGE

The expiration date of the kit is given on the label. Expiration dates for each component are listed on the individual labels. Do not use reagents past the expiration date because the sensitivity decreases.

PRECAUTIONS

1. Rx Only - Prescription Only.
2. Use fecal specimens within 24 hours of collection to obtain optimal results. Specimens that are frozen (-20°C or lower) may lose activity due to freezing and thawing. **AVOID REPEATED CYCLES OF FREEZING AND THAWING THE SPECIMEN. THIS MAY INACTIVATE THE TOXIN.**
3. All specimens should be considered as potentially infectious and should be handled appropriately.
4. The *Positive Toxin Control* contains biologically active toxin and should be treated with caution. The antitoxin does not pose any known hazard.
5. Do not use the reagents beyond their expiration dates.
6. Store reagents according to the manufacturer's recommendations. This prolongs the shelf-life of the reagents.
7. Place reagents back in the refrigerator once the procedure has been completed.

COLLECTION AND HANDLING OF FECAL SPECIMENS

Standard collection and handling procedures used in-house for fecal specimens are appropriate. Specimens should be transported as soon as possible and stored between 2° and 8°C. Whenever possible, test fecal specimens which are less than 24 hours old. Store specimens at -20°C, or lower, if the test cannot be performed within 48 hours. **Freezing and thawing of the specimen, especially multiple times, may result in loss of activity due to degradation of the toxin.** Make sure that specimens are thoroughly mixed prior to performing the assay. This includes complete mixing of the specimen prior to preparing the fecal extract.

PROCESSING FECAL SPECIMENS

NOTE: Fecal specimens less than 24 hours old should be tested whenever possible. Specimens should be stored between 2° and 8°C until tested.

1. Set up one tube for each specimen to be tested. Add 1.8 mL *Diluent* to each tube. Label the tube directly on the side.
2. For formed fecal specimens, use a swab to transfer the fecal specimen to the tube. Coat the swab completely before transferring the specimen. Mix the swab in the diluent to remove as much sample as possible and squeeze the swab against the side of the tube to express any residual liquid. For liquid fecal specimens, use a plastic pipette to transfer 0.2 mL specimen to the tube (1/10 dilution). Make sure the liquid specimens are evenly suspended before transferring.
3. Vortex the tubes for 10 seconds. Transfer an aliquot of the specimen to a microcentrifuge tube (e.g., an Eppendorf tube) and centrifuge the specimen for 5 minutes at high speed (2,000 to 10,000 rpms) in a microcentrifuge.
4. Collect the supernatant fluid and filter it through a 0.45 µm membrane into a sterile tube. Store the filtrate between 2° and 8°C.

THE TISSUE CULTURE PROCEDURE

NOTE: The procedure described below may be used to detect *C. difficile* toxin in fecal specimens with the *C. DIFFICILE TOX-B TEST*. Alternatively, the kit may be used according to other modified in-house procedures. It is important, however, that the *Positive Toxin Control* and the antitoxin reagent be used according to the instructions.

1. Determine the number of wells that will be needed to perform the test. Two wells will be needed as control wells each time the test is performed. In addition, two wells will be needed for each specimen. Tissue culture cells that typically are used include human foreskin cells, MRC-5, WI-38, and Chinese Hamster Ovary cells.
2. Add 1 drop (50 µL) of PBS to one of the two wells in each set. Add 1 drop (50 µL) of *C. difficile* antitoxin to the other well in each set. Set up the wells as illustrated in table 1.

Table 1.

Well	PBS or Antigen	Test sample
Control Well 1 Control Well 2	1 drop (50 µL) PBS 1 drop (50 µL) antitoxin	1 drop (50 µL) positive control 1 drop (50 µL) positive control
Test Well 1 Test Well 2	1 drop (50 µL) PBS 1 drop (50 µL) antitoxin	1 drop (50 µL) fecal filtrate 1 drop (50 µL) fecal filtrate

The final dilution of fecal filtrate in the well is 1/50. Use new pipette tips between specimens. Please note that (i) one set of controls should be used each time the test is performed and (ii) the remaining fecal filtrate should be refrigerated between 2° and 8°C for up to 48 hours in case additional testing is required.

3. After the samples have been applied, incubate the wells at 37°C ± 2°C for 24 to 48 hours. Observe each well for cell rounding, indicating the presence of cytotoxic activity. The cytotoxic effect is considered positive if at least 50% of the cells in the well are rounded. Record the results for each well. The cytotoxic effect may be observed within 8 hours after applying the specimen if high amounts of toxin are present. However, specimens should not be considered negative until after 48 hours.

QUALITY CONTROL

A positive and negative control must be run with each series of test specimens. This includes the *Positive Toxin Control* added to a well containing buffer to demonstrate activity,

and *Positive Toxin Control* added to a well containing antitoxin to demonstrate neutralization of the cytotoxic activity. The positive control should show the typical rounding effect characteristic of *C. difficile* toxin. The test results are not valid unless these performance characteristics are met. If the performance characteristics are not met, please call Technical Services and do not report test results. Upon arrival, the kit should be inspected to ensure that components arrived in satisfactory condition. Test results along with control results should be recorded and reported according to in-house procedures and this information should be stored according to in-house procedures for future reference.

INTERPRETATION OF RESULTS

Wells Containing		INTERPRETATION OF RESULTS
PBS	<i>C. difficile</i> antitoxin	
-	-	Negative - no cytotoxic activity present
+	-	Cytotoxic activity present that is neutralized by <i>C. difficile</i> antitoxin, confirming the presence of <i>C. difficile</i> toxin in the specimen
+	+	Indeterminate results that may be caused by extremely high levels of <i>C. difficile</i> toxin. Fecal filtrate should be diluted an additional 1/10 in Diluent and the test should be repeated to rule out this possibility. Upon repeat testing under these conditions, if the cells in the well containing PBS and highly dilute fecal filtrate are rounded and the cells in the well containing antitoxin and highly diluted fecal filtrate are not, <i>C. difficile</i> toxin is present. If the cells in both wells are round, the specimen contains cytotoxic activity but it is not due to <i>C. difficile</i> toxin.

Control Results: The cells in the well that contains PBS and *Positive Toxin Control* should be round. The cells in the well that contains *C. difficile* antitoxin and *Positive Toxin Control* should appear normal, demonstrating the neutralization of the *Positive Toxin Control* by the antitoxin.

Test Sample Results: When specimen filtrates show evidence of specific *C. difficile* toxin, report "Patient specimen positive for *Clostridium difficile* toxin." When specimen filtrates show no evidence of cytotoxic activity at 48 hours, report "No *Clostridium difficile* toxin detected." When specimen filtrates show evidence of cytotoxic activity that is not neutralized by *C. difficile* antitoxin, report "Nonspecific reaction in patient specimen, not characteristic of *Clostridium difficile* toxin." A repeat specimen may be helpful to rule out patient disease caused by *Clostridium difficile* toxin.

LIMITATIONS OF THE *C. DIFFICILE TOX-B TEST*

1. The *C. DIFFICILE TOX-B TEST* is used to detect *C. difficile* toxin in fecal specimens. Because of the complex nature of feces, specimens should be <24 hours old.
2. Detection of toxin by this test is dependent on the biological activity of the toxin. Specimens must be handled properly and stored as recommended to minimize inactivation of the toxin.
3. Some specimens may give reactions that are inconclusive (e.g., partial rounding of the cells, stretching instead of rounding). Under these conditions, the specimen should be retested or a fresh specimen should be tested. Additional tests that may be used in conjunction with the *C. DIFFICILE TOX-B TEST* include the ELISA for toxin A, ELISA for toxins A and B, isolation of the organism on selective media, and latex agglutination assay for the detection of the organism.

- 6 4. Certain isolates of *C. sordellii* produce the same type of rounding on tissue culture cells as toxigenic *C. difficile*. This is due to the similarities of the *C. sordellii* toxins and the *C. difficile* toxins. The strains of *C. sordellii* that produce these related toxins have not been detected in patients with antibiotic-associated diarrhea and colitis.

EXPECTED VALUES

The prevalence of a positive *C. DIFFICILE TOX-B TEST* in our clinical studies ranged from 5.8% to 19%. The prevalence will vary from location to location and hospitals may experience rates lower or higher than those observed at the sites used in the *C. DIFFICILE TOX-B TEST* evaluation. *Clostridium difficile* disease is primarily a nosocomial disease of elderly patients and hospitals that have higher numbers of elderly patients may experience higher rates. It is important to consider any test results in conjunction with clinical symptoms because some healthy adults and large numbers of healthy infants (up to 50%) will be positive for *C. difficile* toxin A or toxin B. In addition, *C. difficile* carriage rates of 22% to 32% have been reported in cystic fibrosis patients (11-13).

PERFORMANCE CHARACTERISTICS

The *C. DIFFICILE TOX-B TEST* was evaluated for its performance at four different locations in the United States. For the analysis, the *C. DIFFICILE TOX-B TEST* was compared with another commercial Cytotoxicity Test for *C. difficile* toxins, which also is a tissue culture-based assay. All tissue culture assays were done using human foreskin cells. A summary of the results from the clinical investigations is shown in table 2. The range of overall relative sensitivity, specificity, and correlation for the evaluations is presented in parentheses. In studies performed at reference centers in the northeastern and southeastern U.S. and a medical center on the East Coast involving a total of 427 specimens, the *C. DIFFICILE TOX-B TEST* exhibited an overall relative correlation of 100% with the other commercial test. In a study performed at a medical center in the Mid West involving 517 specimens, the *C. DIFFICILE TOX-B TEST* exhibited an overall relative sensitivity and specificity of 96.7% and 98.5%, respectively. The overall relative correlation was 98.3%. When discrepant results were resolved by toxigenic culture, the overall relative sensitivity and specificity increased to 98.3% and 98.7%, with an overall relative correlation of 98.6%.

Table 2

	Commercial Cytotoxicity Assay/Toxigenic Culture	
	Positive	Negative
TOX-B TEST Positive	114	6
TOX-B TEST Negative	1	823

Overall relative sensitivity	99.1% (98.3 to 100%)
Overall relative specificity	99.3% (98.7 to 100%)
Overall relative correlation	99.3% (98.6% to 100%)
Number of specimens	944

CROSS-REACTIVITY

Toxic filtrates from *C. chauveoi*, *C. fallax*, *C. perfringens*, *C. septicum*, and *C. spiroforme* were not neutralized by the *C. difficile* antitoxin. In addition, *C. perfringens* alpha-toxin, *C. histolyticum* collagenase and clostripain, cholera toxin, diphtheria toxin, shiga-like toxin, and staphylococcal enterotoxin were not neutralized. The only organism known to cross-react in the *C. DIFFICILE TOX-B TEST* is *C. sordellii*. Certain toxigenic isolates of *C. sordellii* produce toxins that are very similar to the toxins of *C. difficile*. The hemorrhagic toxin (toxin HT) of *C. sordellii* is very similar to toxin A whereas the lethal toxin (toxin LT) is very similar to toxin B. Antibodies against the *C. difficile* toxins neutralize toxins HT and LT of *C. sordellii*. *Clostridium sordellii* has not been implicated in antibiotic-associated diarrhea and colitis.

REPRODUCIBILITY AND PRECISION

Identical sets of six fecal specimens were sent to three different clinical laboratories located in the Mid West, Southeast, and Northeast U.S. The *C. DIFFICILE TOX-B TEST* results from each laboratory were consistent with the results obtained in-house. In an in-house study, five positive fecal specimens were assayed at time 0, 24, 48, and 72 hours using the *C. DIFFICILE TOX-B TEST*. These specimens were shown to be consistently positive for *C. difficile* toxin at each of the time intervals, as determined by neutralization of the cytotoxic activity. The specimens were stored at 4°C during the testing. In additional in-house studies, three different lots of the positive control reagent were tested in duplicate and consistently gave cytotoxic titers of 10^5 . Further evaluation demonstrated that different lots of the *C. difficile* antitoxin neutralized the *Positive Toxin Control*.

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FOR INFORMATIONAL USE ONLY

BIBLIOGRAPHY [For details on *C. difficile* disease, see our 1988 review (1)]

1. Lyerly, D.M., H.C. Krivan and T.D. Wilkins. 1988. *Clostridium difficile*: its disease and toxins. Clin. Microbiol. Rev. 1: pp. 1-18.
2. D.M. Lyerly, K.E. Saum, D.K. MacDonald and T.D. Wilkins. 1985. Effects of *Clostridium difficile* toxins given intragastrically to animals. Infect. Immun. 47: 349-352.
3. Borriello, S.P., F.E. Barclay, A.R. Welch, J.M. Ketley, T.J. Mitchell, J. Stephen and G.E. Griffin. 1985. Host and microbial determinants of the spectrum of *Clostridium difficile* mediated gastrointestinal disorders. Microecol. Ther. 15:231-236.
4. Lyerly, D.M., N.M. Sullivan, and T.D. Wilkins. 1983. Enzyme-linked immunosorbent assay for *Clostridium difficile* toxin A. J. Clin. Microbiol. 17:72-78.
5. Laughon, B.E., R.P. Viscidi, S.L. Gdovin, R.H. Yolken, and J.G. Bartlett. 1984. Enzyme immunoassays for detection of *Clostridium difficile* toxins A and B in fecal specimens. J. Infect. Dis. 149:781-788.
6. Lyerly, D.M., L.A. Barroso, and T.D. Wilkins. 1992. Characterization of a toxin A-toxin B+ isolate of *Clostridium difficile*. Infect. Immun. 60:4633-4639.
7. Lyerly, D.M., L.M. Neville, D.T. Evans, J. Fill, S. Allen, W. Greene, R. Sautter, P. Hnatuck, D.J. Torpey, and R. Schwalbe. 1998. Multicenter Evaluation of the *Clostridium difficile* TOX A/B TEST. J. Clin. Micro. 36:184-190.
8. Dove, C.H. S.Z. Wang, S.B. Price, C.J. Phelps, D.M. Lyerly, T.D. Wilkins and J.L. Johnson. 1990. Molecular characterization of the *Clostridium difficile* toxin A gene. Infect. Immun. 58: 480-488.
9. Barroso, L.A., S.Z. Wang, C.J. Phelps, J.L. Johnson and T.D. Wilkins. 1990. Nucleotide sequence of *Clostridium difficile* toxin B gene. Nucl. Acids. Res. 18:4004.
10. Krivan, H.C., G.F. Clark, D.F. Smith, and T.D. Wilkins. Cell surface binding site for *Clostridium difficile* enterotoxin: evidence for a glycoconjugate containing the sequence Gal α 1-3Gal β 1-4GlcNac. Infect. Immun. 53: 573-581.
11. Peach, S.L., S.P. Borriello, H. Gaya, F.E. Barclay and A.R. Welch. 1986. Asymptomatic carriage of *Clostridium difficile* in cystic fibrosis patients. J. Clin. Path. 39: 1013-1018.
12. Welkon, C.J., S.S. Long, C.M. Thompson, Jr. and P.H. Gilligan. 1985. *Clostridium difficile* in patients with cystic fibrosis. Amer. J. Dis. Child. 139: 805-808.
13. Wu, T.C., V.P. McCarthy and V.J. Gill. 1983. Isolation rate and toxigenic potential of *Clostridium difficile* isolates from patients with cystic fibrosis. J. Infect. Dis. 148:176.

Technical Support

Further information can be obtained from contacting
TECHLAB® Technical Support:

US + 1 800 TECHLAB
Phone (540) 953-1664
FAX (540) 953-1665
email ts@techlab.com
www.techlab.com

FOR INFORMATIONAL USE ONLY

The TECHLAB Logo and TECHLAB are trademarks of TECHLAB, Inc.

© 2023 TECHLAB, Inc. All rights reserved.