# GIARDIA II™

A Monoclonal ELISA for Detecting *Giardia lamblia* Antigen  
in Fecal Specimens

Catalog No. PT5012 (96 Tests)

**IVD** *In Vitro* Diagnostic Medical Device

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## GIARDIA II™

### INTENDED USE

The *GIARDIA II*™ test is an enzyme immunoassay for the qualitative detection of *Giardia lamblia* cyst antigen in human fecal specimens. It is indicated for use with fecal specimens from patients with diarrhea to determine the presence of *G. lamblia* gastrointestinal infection. This test can be used for fecal specimens submitted for routine clinical testing from adults or children.

**Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician**

### EXPLANATION

*Giardia* is a binucleated flagellated protozoan parasite which exists in two forms: a noninfectious, pear-shaped trophozoite (9 to 20 µm) inhabiting the small intestine and the highly infectious cyst form which is elliptical in shape and ranges in size from 8 to 12 µm (1). Survival outside its host varies greatly between the two forms: the trophozoite which is extremely labile, lasts only a matter of hours outside the body, while the cyst form may survive for several days in an external environment (1). The parasite is responsible for infections due to water contamination and travelers have been found to contract giardiasis from endemic areas (2-5). Transmission also occurs by direct contact especially by asymptomatic carriers and by food contamination (6,7). High-risk categories include young children, immunocompromised patients and those without previous exposure (8). More recently, giardiasis has become a common sexually transmitted disease (9).

*Giardia* has been found in all animal hosts studied and animal fecal contamination especially of water is another route of transmission in humans (4,10,11). Clinical manifestations of giardiasis range from asymptomatic carriage with the passing of cysts to chronic debilitating diarrhea, weight loss and malabsorption (8,12,13). Diagnosis by microscopy has been the most common method used for giardiasis. However, this procedure requires extensive experience and the presence of intact cysts in the feces. An alternative test is the ELISA. This procedure is fairly simple to perform and exhibits increased sensitivity compared to microscopic examination. Fecal specimens are examined using various techniques (14) and accuracy of results are dependent on the skill of the technician. The success rate of fecal examination varies between 50 and 70% (15,16) and several specimens are usually necessary to establish a diagnosis. Giardiasis also may be present in the absence of whole organisms and can be confused with other illnesses such as Crohn's disease and ulcerative colitis (17). When infection is present but parasites are not detected, sampling and testing of duodenal fluid may detect trophozoites, but this method is invasive and expensive (14,16). Detection of the organism and antigens by ELISA provides an alternative method of establishing a diagnosis and is sensitive and specific (18-20). Large numbers of specimens can be tested rapidly and objectively, and the procedure is less labor intensive than some microscopy methods.

### BIOLOGICAL PRINCIPLES

The *GIARDIA II*™ test uses monoclonal and polyclonal antibodies to a cell-surface antigen of the organism. The *Microassay Plate* in the kit contains immobilized monoclonal antibody and the *Conjugate* consists of polyclonal antibody, both of which are specific for the cell-surface antigen. In the assay, an aliquot of a diluted fecal specimen is transferred to a microassay well. The immobilized monoclonal antibody binds the *Giardia* antigen if the antigen is present. Upon addition, *Conjugate* then binds to the antigen/antibody complex. Any unbound materials are removed during the washing steps. Following the addition of substrate, a color is detected due to the enzyme-antibody-antigen complexes that formed in the presence of *Giardia* antigen and *Conjugate*.

### MATERIALS PROVIDED

CONJ	ENZ
------	-----

**Conjugate**, 7 mL (Rabbit polyclonal antibody to a cell-surface antigen of *Giardia* in a protein buffered solution containing 0.02% Thimerosal)\*

**DIL SPE** **Diluent**, 50 mL (Buffered protein solution containing 0.02% Thimerosal)\*. The *Diluent* is also to be used as the negative control solution (see TEST PROCEDURE).

H412 - Harmful to aquatic life with long lasting effects.

P273, P501

**H<sub>2</sub>SO<sub>4</sub> 0.6N** **Stop Solution**, 7 mL (0.6N sulfuric acid). Caution: Avoid contact with skin. Flush with water immediately if contact occurs.

Signal Word: Danger

H314: Causes severe skin burns and eye damage

P260, P264, P280, P301, P330, P331, P303, P361, P353, P363, P304, P340, P310, P321, P305, P351, P338, P501

**CONTROL +** **Positive Control**, 3.5 mL (*Giardia* antigen in a protein buffered solution containing 0.02% thimerosal)\*

**SUBS REAG** **Substrate**, 14 mL (solution containing tetramethylbenzidine and peroxide)

**WASHBUF 20X** **Wash Buffer Concentrate**, 50 mL (20X concentrate containing phosphate-buffered saline, detergent and 0.2% thimerosal)\*

Signal Word: Warning

H373: May cause damage to organs through prolonged or repeated exposure

H411: Toxic to aquatic life with long lasting effects

P260, P273, P314, P391, P501

**MA PLT** **Microassay Plate**, 12 strips, each consisting of 8 wells coated with monoclonal antibody to *Giardia* cell-surface antigen (stored with desiccant)

**2 plastic adhesive sheets**

**100 graduated disposable pipettes**

\*contains mercury

## MATERIALS REQUIRED BUT NOT PROVIDED

Squirt bottle for wash reagent

950 mL distilled water for diluting wash reagent

Absorbent paper

Applicator sticks

ELISA reader capable of reading at 450 nm or 450/620 nm

Vortex mixer

Discard container

Microcentrifuge tubes

## PRECAUTIONS

1. Rx Only – Prescription Only
2. For *in vitro* diagnostic use. For professional use only.
3. Reagents from different kits should not be mixed. Do not use a kit past the assigned expiration date.
4. Reagents should be at room temperature before use.
5. Caps and tips are color-coded; do not mix.
6. Each component in the kit should be inspected for any signs of leakage. Upon arrival, the kit should be inspected to ensure that components are not frozen or warm-to-the touch due to improper shipping conditions.
7. When handling assay wells, avoid scratching the bottom of the wells because this may result in elevated absorbance readings.
8. Hold dropper bottles vertically to ensure proper drop size.
9. Unused microassay wells must be placed back inside of the resealable pouch with the desiccant to protect them from moisture.
10. Perform the washing procedure as directed to avoid high background reactions.
11. The *Substrate* is light sensitive and should be protected from direct sunlight or UV sources.
12. Optimal results are obtained by following the specified test procedure. The concentrations, incubation conditions, and processing specifications have been optimized for sensitivity and specificity. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test.

- 4
13. Handle specimens and used microassay wells as if capable of transmitting infectious agents. Wear disposable gloves when doing the test.
  14. The *20X Wash Buffer Concentrate* contains 0.2% Thimerosal as a preservative. Once diluted to normal use concentration this solution is classified as non-hazardous. The *Stop Solution* contains 0.6N sulfuric acid. Flush with water immediately if contact occurs. Take off contaminated clothing and wash it before reuse. Handle reagents according to existing regulations for laboratory safety and good laboratory practice. Safety Data Sheets for this product are available upon request, contact technical support.
  15. Follow your national, regional, and local ordinances accordingly for waste disposal regulations. Do not place in trash, dispose of as hazardous waste.

### **SHELF LIFE AND STORAGE**

The expiration date of the kit is given on the label. Expiration dates for each component are listed on individual labels. The kit should be stored between 2° and 8°C and should be returned to the refrigerator as soon as possible after use.

### **SPECIMEN HANDLING**

1. Standard collection and handling procedures used in-house for fecal specimens for culture are appropriate. No modification of collection methods used for standard microscopic O&P examinations is needed. Fecal specimens may be used as unpreserved or frozen, or in preservation media of 10% buffered formalin or Sodium Acetate Formalin (SAF).
2. Unpreserved specimens should be kept between 2° and 8°C and tested within 24 hours of collection. Specimens that cannot be tested within this time should be stored at -20°C or less until tested.
3. Preserved specimens may be kept at room temperature and tested within 18 months of collection.
4. Concentration steps are not needed (or recommended) for fecal specimens.
5. Make sure that specimens are thoroughly mixed (vortexed) prior to performing the test. This includes complete mixing of the specimen prior to transfer to the *Diluent* and/or microassay well.
6. All dilutions must be made with *Diluent*.

### **REAGENT PREPARATION**

1. The contents of the kit should be brought to room temperature before use.
2. Prepare *1X Wash Solution*. The *Wash Buffer Concentrate* is supplied as a *20X* concentrate (a precipitate may be noticed). It should be diluted to a total volume of 1 liter by adding 50 mL of the concentrate to 950 mL of distilled water. The *1X Wash Solution* can be stored between 2° and 8°C.

### **SPECIMEN PREPARATION**

1. **Fresh/Frozen Fecals:** Frozen fecal specimens should be thawed. Add 400 µL of *Diluent* to a microcentrifuge tube (one per sample), then add 100 µL (2<sup>nd</sup> graduation mark on pipette) of sample to the tube and mix well. This is a 1:5 dilution. If the specimen cannot be pipetted, use an applicator stick to transfer approximately 0.1 gram of feces. This is about the size of a small pea (about 4 mm in diameter). Perform final dilution in microassay wells as directed in the PROCEDURE.
2. **Preserved Fecals:** Mix (vortex) contents of container thoroughly before transferring specimen. No further processing or dilution is necessary. Perform final dilution in microassay wells as directed in PROCEDURE.

### **PROCEDURE**

1. Two control wells must be used each time the test is performed. These wells serve as positive and negative controls. One well is needed for each test run. Hold the dropper bottles vertically when adding the drops. Identification marks may be written directly on side of well. Add 1 drop (50 µL) of the *Positive Control* (black cap) to the positive control well and 2 drops (100 µL, 2<sup>nd</sup> graduation mark on transfer pipette) of the negative control (i.e., *Diluent*) to the negative control well.
2. If fresh or frozen samples (i.e. unpreserved) are to be used, make a 1:5 dilution of the sample as stated in SPECIMEN PREPARATION above.

3. Transfer 100  $\mu\text{L}$  of *Diluent* to each test well of the *Microassay Plate*. Using plastic pipettes, transfer 1 drop (50  $\mu\text{L}$ , 1<sup>st</sup> graduation mark on pipette) of sample (preserved or diluted as above) to each test well already containing *Diluent* and gently tap the wells to mix. Seal with a plate sealer and incubate 1 hour at room temperature.
4. Shake out contents of assay wells into a discard pan. Wash each well using the 1X *Wash Solution* in squirt bottle with fine-tipped nozzle, directing the *Wash Solution* to the bottom of the well with force. Fill the wells, shake the wash solution out of the well into a discard pan and slap the plate hard onto a dry paper towel. Repeat the washing step 3X (for a total of 4 washes). If any fecal material remains in the wells wash the plate until it appears clean.
5. After washing, completely remove any residual liquid in the wells by striking the plate onto a dry paper towel until no liquid comes out. Dispose of paper towels and specimen containers properly.
6. Add 1 drop (50  $\mu\text{L}$ ) of *Conjugate* (red cap) to each well and gently tap to mix. Seal with a plastic adhesive sheet. Incubate the wells for 30 minutes at room temperature.
7. Repeat the washing procedure (Steps 4 and 5).
8. Add 2 drops (100  $\mu\text{L}$ ) of *Substrate* (blue cap) to each well. Gently tap the wells to mix. Incubate the wells at room temperature for 10 minutes.
9. Add 1 drop (50  $\mu\text{L}$ ) of *Stop Solution* (yellow cap) to each well. Gently tap the wells and wait 2 minutes before reading. The addition of the *Stop Solution* converts the blue color to a yellow color, which may be quantitated by measuring the absorbance at 450 nm on a microplate ELISA reader. The instrument should be blanked against air. Wipe the underside of each well before measuring the absorbance. If a dual reader is used, blank against air at 620 nm and read at 450 nm. Visual readings should be noted. Read within ten minutes after adding *Stop Solution*.

### QUALITY CONTROL

1. A positive and a negative control must be run with each series of test specimens.
2. Each positive control well should be an easily visible yellow color and should give an absorbance of 0.500 or higher. Any well that gives a positive reading without visible color should be repositioned, wiped on the underside of the well, and read again.
3. Negative control wells should be colorless or they may have a faint yellow color but they must have an absorbance value of  $<0.150$  when read on single wavelength ( $\text{OD}_{450\text{nm}}$ ) or  $<0.090$  when read on dual wavelength ( $\text{OD}_{450/620\text{nm}}$ ).
4. Test results are not valid unless the performance characteristics of the positive and negative controls are met. If these results are not observed, call Technical Services.
5. Test results along with control absorbance values should be recorded and reported according to in-house procedures and should be stored according to in-house procedures for future reference.

### INTERPRETATION OF RESULTS

Absorbance		Visual Color	Interpretation
450 nm	450/620 nm		
$<0.150$	$<0.090$	Clear to slight yellow	Negative - Below detectable limits of assay
$\geq 0.150$	$\geq 0.090$	Pale yellow to strong yellow	Positive - Specimen contains <i>Giardia</i> antigen

#### Visual Interpretation

**Negative:** Any sample that is colorless or resembles the negative control well in intensity of color.

**Positive:** Any sample that is obviously more yellow than the negative control well.

NOTE: The negative control, as well as some test wells, may show some slight yellow color. A sample well must be obviously more yellow than the negative control well to be called a positive result.

**6**

**Spectrophotometric Interpretation**

1. Determine the absorbance value of the negative control. The negative control reading should be  $< 0.150 \text{ OD}_{450}$  or  $< 0.090 \text{ OD}_{450/620}$ . If not, the test is invalid and should be repeated, paying attention to the wash procedure.
2. The reading for the *Positive Control* should be  $\geq 0.500$ .
3. Test results  
 Negative:  $< 0.150$  (absorbance at 450 nm) or  $< 0.090$  (absorbance at 450/620 nm)  
 Positive:  $\geq 0.150$  (absorbance at 450 nm) or  $\geq 0.090$  (absorbance at 450/620 nm)

**LIMITATIONS OF THE PROCEDURE**

1. The *GIARDIA II™* test detects the presence of *Giardia* antigen in fecal specimens. This antigen is present primarily in the cyst form of the parasite.
2. The test results should be interpreted by a physician in consideration of other laboratory results and clinical history.
3. Concentrated fecal specimens should not be tested in the *GIARDIA II™* test and will not give accurate results.
4. The predictive value of a positive result decreases when testing in a low prevalence population.
5. The *GIARDIA II™* test is for the qualitative detection of *Giardia* antigen in fecal specimens. It has not been evaluated for quantitative determinations of organism load, and the magnitude of the absorbance value does not correlate with organism load.

**EXPECTED VALUES**

Normal healthy individuals should not be infected with *Giardia* and should test negative in the *GIARDIA II™* test. A positive test result in the *GIARDIA II™* test indicates that the person is shedding detectable amounts of *Giardia* antigen. The incidence of *Giardia* infection varies significantly between populations and geographic regions. Children in daycare settings have exhibited higher rates of infection with *Giardia* than the normal population (21). In addition, homosexual men have shown higher rates of infection (22).

**PERFORMANCE CHARACTERISTICS**

In a study performed at TECHLAB Inc. the *GIARDIA II™* test was evaluated using fresh fecal specimens and fecal specimens preserved in 10% Formalin or SAF. Specimens used in the study were identified as *Giardia*-positive by microscopy at outside laboratories. The positive specimens were from patients ranging from age three to eighty-nine. Of 203 specimens analyzed, 110 were determined to be positive by microscopy. The *GIARDIA II™* test detected all 110 of the positive specimens. A total of 93 *Giardia* negative specimens were tested. All of the 93 *Giardia* negative specimens were negative in the *GIARDIA II™* test. The results are shown in Table 1.

**Table 1. Comparison of the *GIARDIA II™* test with microscopy (n = 203).**

		Microscopy	
		+	-
<i>GIARDIA II™</i>	+	110	0
	-	0	93

Sensitivity	100%
Specificity	100%
Predictive Positive Value	100%
Predictive Negative Value	100%
Correlation	100%

In an additional in-house study, the *GIARDIA II*<sup>TM</sup> test was compared to another commercially available ELISA using fresh fecal specimens and fecal specimens preserved in 10% formalin or SAF. Specimens used in the study were identified as *Giardia*-positive by microscopy by outside laboratories. Of 128 specimens analyzed, 88 were determined to be positive by microscopy. The *GIARDIA II*<sup>TM</sup> test detected all 88 of the positive specimens while the other commercially available ELISA detected 86 of the 88 positive specimens. A total of 40 negative specimens were tested. All 40 of the *Giardia*-negative specimens were negative in both the *GIARDIA II*<sup>TM</sup> test and the other commercially available ELISA. The results are shown in Table 2.

**Table 2. Comparison of another commercial ELISA and *GIARDIA II*<sup>TM</sup> test with microscopy (n = 128).**

		Microscopy				Microscopy	
		+	-			+	-
Other Commercial ELISA	+	86	0	<i>GIARDIA II</i> <sup>TM</sup>	+	88	0
	-	2	40		-	0	40

	Other Commercial ELISA	<i>GIARDIA II</i> <sup>TM</sup>
Sensitivity	97.7%	100%
Specificity	100%	100%
Predictive Positive Value	100%	100%
Predictive Negative Value	95.2%	100%
Correlation	98.4%	100%

### CROSS-REACTIVITY

The *GIARDIA II*<sup>TM</sup> test was evaluated using fecal specimens found to be positive for a variety of intestinal pathogens. No cross-reactivity was observed with fecal specimens that contained any of the pathogens listed below. The number of specimens tested with each organism is shown in parenthesis.

<i>Ascaris lumbricoides</i> (2)	<i>Blastocystis hominis</i> (8)	<i>Chilomastix mesnili</i> (2)
<i>Clonorchis</i> (1)	<i>Cryptosporidium parvum</i> (8)	<i>Dientamoeba fragilis</i> (3)
<i>Endolimax nana</i> (3)	<i>Entamoeba coli</i> (4)	<i>Entamoeba hartmanni</i> (2)
<i>Entamoeba histolytica</i> (2)	Hookworm (2)	<i>Iodamoeba butchlii</i> (2)
<i>Isospora belli</i> (2)	Rotavirus (11)	<i>Schistosoma mansoni</i> (2)
<i>Strongyloides stercoralis</i> (3)	<i>Trichuris trichiura</i> (2)	

### PRECISION - INTER-ASSAY

For the determination of inter-assay performance, four positive and four negative fecal specimens were assayed a total of five times over a five day period in the *GIARDIA II*<sup>TM</sup> test. The % Coefficient of Variation (CV) of positive specimens ranged from 3.5 to 8.0, with an average of 5.2. The % Coefficient of Variation (CV) of negative specimens ranged from 3.4 to 14.9, with an average of 8.3.

### PRECISION - INTRA-ASSAY

The intra-assay % Coefficient of Variation (CV) was determined by analyzing four positive and four negative fecal specimens. Each specimen was assayed in twelve test wells. The % Coefficient of Variation (CV) of positive specimens ranged from 2.9 to 11.4, with an average of 7.1. The % Coefficient of Variation (CV) of negative specimens ranged from 7.9 to 36.9, with an average of 22.8.

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© 2023 TECHLAB, Inc. All rights reserved.
GIARDIA II, the TECHLAB Logo and TECHLAB are trademarks of TECHLAB, Inc.