TRI-COMBO PARASITE SCREEN

A Monoclonal ELISA for Detecting *Giardia*, *Cryptosporidium* and
E. histolytica Antigens in Fecal Specimens

Catalog No. T30408 (96 Tests)

IVD In Vitro Diagnostic Medical Device

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TRI-COMBO PARASITE SCREEN

INTENDED USE

The TECHLAB® TRI-COMBO PARASITE SCREEN test is an enzyme immunoassay for the simultaneous qualitative detection of *Giardia* spp., *Cryptosporidium* spp. and/ or *E. histolytica* antigen in human fecal specimens. The test is indicated as an aid in the diagnosis of gastrointestinal infection when giardiasis, cryptosporidiosis and amebiasis is suspected. The test does not differentiate between the three parasites and follow-up testing is required for all positive results to confirm the specific diagnosis.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician

EXPLANATION

Giardia is a binucleated, flagellated protozoan parasite which exists in two forms: a non-infectious, pear-shaped trophozoite (9 to 21 µm) inhabiting the small intestine and the highly-infectious cyst form, which is elliptical in shape and ranges in size from 8 to 12 µm (reviewed in 1). Survival outside its host varies greatly between the two forms: the trophozoite, which is extremely labile, lasts only a matter of hours outside the body, while the cyst form may survive for several days in an external environment (1). *Giardia* often is responsible for infections due to water contamination, and travelers have been found to contract giardiasis from endemic areas. Transmission also occurs by direct contact, especially with asymptomatic carriers and by food contamination. Giardiasis is known to be a common sexually-transmitted disease (2). Animal fecal contamination, especially of water, is another route of transmission in humans (1).

Clinical manifestations of giardiasis range from asymptomatic carriage with the passing of cysts to chronic debilitating diarrhea, weight loss and malabsorption. High-risk categories include young children, immunocompromised patients and those without previous exposure (1).

Cryptosporidium is a protozoan parasite of vertebrates originally thought to cause diarrhea only in animals (3). In 1976, the first human infection was reported (4). Since then, *Cryptosporidium* spp. has been found to be associated with diarrheal illness in most parts of the world and is a frequent cause of traveler's diarrhea. The disease is transmitted by the thick-walled oocyst form, 2-6 µm in diameter, which is remarkably resistant to common disinfectants and routine chlorination of drinking water. Person-to-person transmission, especially among children, is common. *Cryptosporidium* has little or no host specificity, and animals such as rodents, cattle and domestic pets can serve as a reservoir for zoonotic transmission to humans, transmission occurring either by direct contact or by contamination of water supplies with fecal matter. Cryptosporidiosis is a serious opportunistic infection in patients with acquired immunodeficiency syndrome (AIDS) and is potentially a sexually-transmitted disease (5,6).

Clinical manifestations of cryptosporidiosis include cholera-like diarrhea, abdominal pain, nausea, vomiting and weight loss. In normally healthy persons, the infection is usually self-limiting and of short term. In AIDS and other immunocompromised patients, cryptosporidiosis can result in prolonged and life-threatening illness due to excessive fluid loss (7). In these patients the infection may also spread to the respiratory and biliary tracts (5).

The life cycle of *Entamoeba histolytica* is similar to that observed with other amoebae. The organism exists either as a trophozoite, ranging from 12-60 µm in diameter, or as a cyst (10-15 µm) (8). Humans serve as the primary reservoir, with the organism being spread through ingestion of contaminated food and water or through sexual contact. The cyst is very stable. Once it enters the intestine, it begins to divide into trophozoites that bind to the intestinal mucosa via a galactose or N-acetyl-D-galactosamine-binding lectin referred to as the galactose adhesin. Once the trophozoites have attached, they release tissue-damaging enzymes and proteins that lyse the mucosal cells. The galactose adhesins of *E. histolytica* and *E. dispar* cross-react serologically but contain distinct epitopes. The adhesin is antigenically conserved, but monoclonal antibodies are used to distinguish between the adhesin from pathogenic *E. histolytica* and non-pathogenic *E. dispar* (reviewed in 9).

E. histolytica and *E. dispar* are intestinal parasites that infect approximately half a billion people worldwide annually (10). Only malaria and schistosomiasis are believed to be more prevalent parasitic causes of morbidity and mortality. Of the huge number of persons infected, most are infected with *E. dispar*, which has not been associated with disease. Infection with *E. dispar* rather than *E. histolytica* is believed to explain, at least in part, the low prevalence of disease considering the high rate of infection. Even so, it is estimated that approximately 10% of the half billion people infected each year are infected with *E. histolytica*. These individuals can become symptomatic and develop colitis and liver abscesses, resulting in a mortality rate estimated between 40,000 and 120,000 persons annually. It is necessary to distinguish between the two species because *E. histolytica* is pathogenic, whereas *E. dispar* is not associated with colitis or liver abscess. Inaccurate diagnosis may result in unwarranted and unnecessary treatment with harsh drugs with side effects.

Patients with pathogenic *E. histolytica* may show a wide range of conditions. Some are completely asymptomatic. These persons may shed millions of cysts daily, representing a potential reservoir for dissemination. Some patients may show mild diarrhea that develops into bloody diarrhea with abdominal cramps, eventually resulting in fulminant colitis. Because of the tissue damage that can occur, perforation of the intestine may result and the amoebae may disseminate to other parts of the body. Approximately 10% of persons with invasive amebiasis develop liver abscesses (10).

The most common method used to detect giardiasis, cryptosporidiosis and amebiasis has been diagnosis by ova and parasite (O&P) light microscopy. However, accuracy of O&P results depends upon the skill of the technician and also relies on the presence of intact cysts in the feces, which may not be present in all samples. In addition, it is only rarely possible to distinguish between the pathogenic and non-pathogenic species of *Entamoeba* using microscopy. The success rate of fecal examination by microscopy varies between 50 and 70%, and multiple specimens are usually necessary to establish a diagnosis. When infection is present but parasites are not detected by microscopy, sampling and testing of duodenal fluid may detect trophozoites, but this method is invasive and expensive. Detection of the organism and antigens by ELISA provides an alternative method of diagnosis and is sensitive and specific (9). The ELISA procedure is straightforward to perform and exhibits increased sensitivity compared to microscopic examination. Large numbers of specimens can be tested rapidly and objectively, and the procedure is less labor-intensive than most microscopy methods.

PRINCIPLE OF THE TEST

The **TRI-COMBO PARASITE SCREEN** test uses monoclonal and polyclonal antibodies to cell-surface antigens of *Giardia*, *Cryptosporidium* and *E. histolytica*. The **Microassay Plate** in the kit contains immobilized monoclonal antibodies against the antigens, and the **Conjugate** consists of polyclonal antibodies against the antigens. In the assay, an aliquot of a diluted fecal specimen is transferred to a microassay well. The immobilized monoclonal antibodies bind the *Giardia*, *Cryptosporidium* and/or *E. histolytica* antigens if they are present. Upon addition, **Conjugate** then binds to the antigen/antibody complex. Any unbound materials are removed during the washing steps. Following the addition of **Substrate**, a color is detected due to the enzyme-antibody-antigen complexes that formed in the presence of antigens and conjugate.

MATERIALS PROVIDED

CONJ | ENZ

Conjugate (7 mL) – Antibodies against antigens of *Giardia* spp., *Cryptosporidium* spp., and *E. histolytica* in a buffered protein solution containing 0.05% ProClin® 300

Signal Word: Warning

H317: May cause an allergic skin reaction

P261, P272, P280, P302, P352, P333, P313, P321, P362, P363, P364, P501

DIL | SPE

Diluent (50 mL) – Buffered protein solution containing 0.02% thimerosal.

The Diluent is also to be used as the negative control solution (see TEST PROCEDURE).*

H₂SO₄	0.6N	Stop Solution (7 mL) – 0.6 N sulfuric acid. CAUTION: Avoid contact with skin and eyes; flush with water immediately if contact occurs Signal Word: Danger H314: Causes severe skin burns and eye damage P260, P264, P280, P301, P330, P331, P303, P361, P353, P363, P304, P340, P310, P321, P305, P351, P338, P501	
CONTROL	+	Giardia Positive Control (3.5 mL) – <i>Giardia</i> antigen in a buffered protein solution with 0.02% thimerosal*	
CONTROL	+	Cryptosporidium Positive Control (3.5 mL) – <i>Cryptosporidium</i> antigen in a buffered protein solution with 0.02% thimerosal*	
CONTROL	+	E. histolytica Positive Control (3.5 mL) – <i>E. histolytica</i> antigen in a buffered protein solution with 0.02% thimerosal*	
SUBS/REAG		Substrate (14 mL) – solution containing tetramethylbenzidine and peroxide	
WASHBUF	20X	Wash Buffer Concentrate (50 mL) – 20X concentrate containing phosphate buffered saline, detergent, and 0.2% thimerosal* Signal Word: Warning H373: May cause damage to organs through prolonged or repeated exposure H411: Toxic to aquatic life with long lasting effects P260, P273, P314, P391, P501	 
MA	PLT	Microassay Plate – 12 strips, each consisting of 8 wells coated with antibodies to <i>Giardia</i> spp., <i>Cryptosporidium</i> spp., and <i>E. histolytica</i> antigen (stored with desiccant)	

*contains mercury

ACCESSORIES

100 Disposable plastic transfer pipettes

2 Plastic adhesive sheets

1 Wash Solution Label

MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

Squirt bottle for wash reagent

Vortex mixer

950 mL distilled water for diluting wash reagent

Discard container

ELISA reader capable of reading at 450 nm or 450/620 nm

Absorbent paper

Small tubes for dilution of fecal specimens (e.g., microcentrifuge tubes)

Applicator sticks

SHELF LIFE AND STORAGE

The expiration date of this kit is given on the box label. Expiration dates for each component are listed on the individual labels. The kit should be stored between 2°C and 8°C and should be returned to the refrigerator as soon as possible after use.

PRECAUTIONS

1. Rx Only – Prescription use only.
2. For *in vitro* diagnostic use. For professional use only.
3. Each component in the kit should be inspected for any signs of leakage. Upon arrival, the kit should be inspected to ensure that components are not frozen or warm to the touch due to improper shipping conditions.
4. Reagents from different kits should not be mixed. Do not use a kit past the assigned expiration date.
5. Bring all components to ROOM TEMPERATURE BEFORE USE!
6. Do not freeze the reagents. The kit should be stored between 2°C and 8°C.
7. Caps and tips are color-coded; do NOT mix or interchange!
8. When handling assay wells, avoid scratching the bottom of the wells because this may result in elevated absorbance readings.
9. Unused microwells must be placed back inside of the resealable pouch with the desiccant and tightly sealed to protect them from moisture.
10. Hold reagent bottles vertically when dispensing to ensure proper drop size and correct volume.

11. Microbial contamination of reagents may decrease the accuracy of the assay. Avoid microbial contamination of reagents by using sterile disposable pipettes if removing aliquots from reagent bottles.
12. All reagents, with the exception of the *Wash Buffer Concentrate*, are supplied in ready-to-use bottles. Reagents can be dispensed directly from the dropper bottles or decanted for use with multichannel pipettes. If excess reagent has been decanted, the excess should be discarded. Do not pour back into the bottle. The *Substrate* should be stored in and used from the light-protected bottle in which it is supplied. If an aliquot is removed from the original bottle for any reason, do not return unused *Substrate* to the original bottle.
13. Perform the washing procedure as directed to avoid high background reactions.
14. The *Substrate* is light sensitive and should be protected from direct sunlight or UV sources.
15. The test has been optimized for sensitivity and specificity. Do not deviate from the specified procedure. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test.
16. Fecal specimens may contain potentially infectious agents and should be handled at "Biosafety Level 2" as recommended in the CDC/NIH Manual, "Biosafety in Microbiological and Biomedical Laboratories."
17. Handle specimens and used microassay wells as if capable of transmitting infectious agents. Wear disposable gloves when performing the test.
18. Reagents contain 0.02% thimerosal as a preservative and should be handled with normal laboratory caution.
19. The *Conjugate* reagent contains 0.05% ProClin® 300 as a preservative. Although the concentration is low, ProClin® 300 is known to be harmful. If skin irritation or rash occurs, get medical advice/attention. The *20X Wash Buffer Concentrate* contains 0.2% thimerosal as a preservative. Once diluted to normal use concentration this solution is classified as non-hazardous. The *Stop Solution* contains 0.6 N sulfuric acid. Flush with water immediately if contact with skin or eyes occurs. Take off contaminated clothing and wash it before reuse. Handle reagents according to existing regulations for laboratory safety and good laboratory practice. Safety Data Sheets for this product are available upon request, contact technical support.
20. Follow your national, regional, and local ordinances accordingly for waste disposal regulations. Do not place in trash, dispose of as hazardous waste.

PRELIMINARY PREPARATIONS

1. **All reagents must be at room temperature prior to use in the assay.**
2. **Prepare 1X Wash Solution.** The *Wash Buffer Concentrate* is supplied as a 20X concentrate (a precipitate may be noticed). It should be diluted to a total volume of 1 liter by adding 50 mL of the concentrate to 950 mL of distilled water. The 1X *Wash Solution* can be stored between 2° and 8°C.
3. **Assay Strip Preparation.** Each strip contains 8 wells coated with antibodies to *Giardia* spp., *Cryptosporidium* spp., and *E. histolytica* antigen. Each specimen or control will use one of these coated wells. Determine the number of wells to be used. Avoid contact with the base of the wells. Assay wells not used must be returned to the foil pouch and carefully resealed with desiccant.

COLLECTION, HANDLING, AND STORAGE OF FECAL SPECIMENS

Standard collection and handling procedures used in-house for fecal specimens are appropriate. No modification of collection methods used for standard microscopic O&P examinations is needed.

Acceptable Sample Types
Fresh Fecal Specimens
Frozen Fecal Specimens (frozen undiluted)
Fecal specimens in transport media (e.g., Cary Blair, C&S)


Do Not Use
Fecal specimens in Formalin-based fixative (e.g., Sodium Acetate Formalin, 10% formalin)
Fecal specimens in alcohol-based fixative (e.g., polyvinyl alcohol)
Concentrated Fecal Specimens

Sample Storage Temperature	Acceptable length of storage	Comments
Refrigerated (2°C – 8°C)	10 days	Fresh samples that will be tested within 24 hours may remain at room temperature (18°C – 25°C). If testing is not planned in <24 hours, refrigerate (2°C – 8°C) as soon as possible following collection.
Frozen ≤ -10°C	8 weeks	Freeze specimens and store at ≤ -10°C if testing cannot be performed within 10 days of collection. Thaw at room temperature. Freezing and thawing multiple times may result in loss of specimen reactivity due to antigen degradation.

Samples in Transport Media	
Storage Condition	Recommended storage time
Samples stored in Cary Blair media between 20°C and 25°C	96 hours
Samples stored in C&S media between 2°C and 8°C	96 hours
Samples stored in C&S media between 20°C and 25°C	48 hours

- Use standard in-house collection and handling procedures for fecal specimens. Collect fecal specimens in clean, leak-proof containers.
- Do not store fecal specimens in the *Diluent*.
- Set up and label one test tube for each sample as necessary. **Add 400 µL *Diluent* to each tube.**
- Obtain one disposable plastic transfer pipette (supplied with the kit) for each sample.

Transfer Pipette

- Mix all specimens thoroughly regardless of consistency – it is essential that the samples be evenly suspended before sampling.**
 - For **fresh or frozen/thawed specimens**, using the disposable plastic transfer pipette, add 100 µL (second graduation mark) of fecal specimen to the tube containing *Diluent* and mix well. If the specimen cannot be pipetted, use an applicator stick to transfer approximately 0.1 g of feces. This is about the size of a small pea (about 4 mm in diameter).

- b. For **specimens in Transport Medium** (e.g., Cary Blair, C&S), add 400 μL of the specimen to the tube containing *Diluent*. This can be achieved by pipetting 200 μL (third graduation mark in above diagram) of specimen TWICE for a total of 400 μL .
6. Thoroughly mix each diluted specimen. Proper mixing can be achieved by vortexing or inverting the tube several times.
7. If using semi-automated or automated washing equipment, once diluted, specimens must be centrifuged (5000 x g for 10 minutes) to remove any particulate matter from the supernatant before transfer to assay wells.

TEST PROCEDURE

1. Bring all reagents and the required number of test strips to room temperature before use.
2. **Prepare Controls**
 - a. Positive Controls – shake each *Positive Control* bottle for several seconds, then add one drop to each positive control well.
 - b. Negative Controls – add 100 μL of *Diluent* to an additional well to serve as the negative control.
3. **Using a separate transfer pipette for each specimen, transfer 100 μL of diluted specimen to the assay well of the *Microassay Plate*.**
4. Cut the adhesive plastic sheet to the size necessary to cover the wells. **Cover the wells and incubate them at for 1 hour at room temperature.**
5. Shake out contents of assay wells into a discard pan.
6. **Wash each well using the 1X Wash Solution in a squirt bottle with a fine-tipped nozzle**, directing the *Wash Solution* to the bottom of the well with force. Fill the wells, and then shake the *Wash Solution* out of the well into a discard pan. Slap the inverted plate on a dry paper towel.
 Note: If using semi-automated or automated washing equipment, add 350 μL of 1X *Wash Solution* to each well. Wash for a total of 5 times.
7. Repeat step 6 four additional times using a dry paper towel each time. If any particulate matter is seen in the wells, continue washing until all the particulate matter is removed.
8. After washing, completely remove any residual liquid in the wells by striking the plate onto a dry paper towel until no liquid comes out. Dispose of paper towels and specimen containers properly.
9. **Add 1 drop (50 μL) of Conjugate (red cap) to each well.**
10. Seal with a plastic adhesive sheet and gently tap to mix. **Incubate the wells for 30 minutes at room temperature.**
11. After incubating, **repeat the washing procedure described in steps 5-8.**
12. **Add 2 drops (100 μL) of Substrate (blue cap) to each well.** Gently tap the wells to mix the substrate. Incubate the wells at room temperature for 10 minutes. Gently tap the wells at 5 minutes.
13. **Add 1 drop (50 μL) of Stop Solution (yellow cap) to each well.** Gently tap the wells and wait 2 minutes before reading. The addition of the *Stop Solution* converts the blue color to a yellow color which may be quantitated by measuring the optical density at 450 nm on a microplate ELISA reader. The instrument should be blanked against air. If a dual wavelength reader is used, blank against air at 620 and read at 450 nm. Wipe the underside of each well before measuring the optical density. If an ELISA reader is unavailable, the test may be read visually in good light against a white background. Read within ten minutes after adding *Stop Solution*.

QUALITY CONTROL

1. Positive and negative controls must be run with each series of test specimens. The positive controls demonstrate that the assay is functioning properly for the detection of each parasite antigen in fecal specimens. The negative control demonstrates that the assay is not reacting nonspecifically.

2. Positive and negative controls must fall within their respective ranges or the test results are not valid.
 - a. **Positive Control must be a visible yellow color.** If read on a spectrophotometer, the OD at 450 nm or using dual wavelength at 450/620 nm must be ≥ 0.500 . Any well that gives a positive reading without visible color should be repositioned, wiped on the underside of the well, and read again.
 - b. **Negative Control must be visually clear or have a faint yellow color.** If read on a spectrophotometer, the OD at 450 nm must be < 0.200 . If read at 450/620 nm the absorbance must be < 0.160 . If not, the test is invalid and should be repeated, paying attention to the wash procedure.
3. Visual readings must be taken in good light against a white background.
4. Test results are not valid unless the performance characteristics of the positive and negative controls are met. If these results are not observed, call Technical Services.

INTERPRETATION OF RESULTS

	Spectrophotometric Reading	
	Single Wavelength at 450 nm	Dual Wavelength at 450/620 nm
Negative	OD < 0.200	OD < 0.160
Positive	OD ≥ 0.200	OD ≥ 0.160

Visual Interpretation

The negative control well should be colorless or have only a faint yellow color. The positive control well should give a visible yellow color. If these results are not observed, call Technical Services. A test sample is considered positive if it has an obvious yellow color when compared to the negative control well. It may be less yellow or more yellow than the color observed in the positive control well. A test sample is considered negative if the reaction is colorless or less yellow than the negative control well.

A positive result indicates that antigen from *Giardia* spp., *Cryptosporidium* spp. and/or *E. histolytica* is present in the specimen. A negative result indicates that *Giardia* spp., *Cryptosporidium* spp. and/or *E. histolytica* antigens are absent or the level is below the detection limit of the test.

LIMITATIONS OF THE PROCEDURE

1. The *TRI-COMBO PARASITE SCREEN* test detects the presence of *Giardia*, *Cryptosporidium* and/or *E. histolytica* antigens in fecal specimens. Test results should be taken under consideration by the physician in light of the clinical history and physical examination of the patient.
2. The *TRI-COMBO PARASITE SCREEN* test is intended for the qualitative detection of *Giardia* spp., *Cryptosporidium* spp., and/or *E. histolytica* antigens in fecal specimens. The intensity of the color should not be interpreted quantitatively.
3. Concentrated fecal specimens should not be used in the test and may not give accurate results.
4. Due to the small number of positive specimens collected during the prospective clinical study, performance characteristics for *E. histolytica* and *Cryptosporidium* spp. were also established with retrospective clinical specimens.

EXPECTED VALUES

Normal healthy individuals should not be infected with *Giardia*, *Cryptosporidium* or *E. histolytica* and should test negative in the *TRI-COMBO PARASITE SCREEN* test. A positive test result in the *TRI-COMBO PARASITE SCREEN* test indicates that the person is shedding detectable amounts of *Giardia*, *Cryptosporidium* and/or *E. histolytica* antigen. The incidence of *Giardia*, *Cryptosporidium* and *E. histolytica* infection varies significantly between populations and geographic regions. Children in daycare settings

have exhibited higher rates of infection with *Giardia* than the general population (11). In addition, homosexual men have shown higher rates of infection (9,12). In general, laboratory-confirmed incidence of cryptosporidiosis in developed countries ranges from 1 to 2% overall with a higher incidence in children. It is estimated that *Entamoeba histolytica* infects about 50 million people around the world (9). Roughly 90% of these persons remain asymptomatic, whereas about 10% develop clinical symptoms ranging from gastrointestinal disease to liver abscesses. High risk groups include persons who have traveled abroad, immigrants, immunocompromised persons, migrant workers, and active male homosexuals (9,14). Nonpathogenic strains (*E. dispar*) are predominant among male homosexuals (15). The disease often is transmitted by asymptomatic carriers of *E. histolytica*.

PERFORMANCE CHARACTERISTICS

Prospective Study

The performance of the *TRI-COMBO PARASITE SCREEN* test was evaluated at three independent sites. The three sites yield a total of 14 microscopy positive samples (13 *Giardia* and one *E. histolytica*). The remaining 740 samples were negative. Of the 740 negative samples, four were *E. histolytica* positive by molecular comparison only but were negative by the *TRI-COMBO PARASITE SCREEN* test, negative by microscopy, and negative by FDA-cleared antigen test for *E. histolytica*. Table 1 summarizes the performance observed, which is primarily a study to evaluate specificity due to the low number of positive specimens (see retrospective study for evaluation of sensitivity). Prospective testing was also read visually and performance was not significantly different from spectrophotometric readings.

Table 1. Summary of prospective clinical performance comparing the *TRI-COMBO PARASITE SCREEN* test to microscopy for *Giardia*, *Cryptosporidium* and *E. histolytica*

<i>TRI-COMBO PARASITE SCREEN</i> (N = 754)	Microscopy	
	Positive	Negative
Positive	13	14*
Negative	1	726
95% Confidence Limits		
Sensitivity	92.9%	68.5% - 98.7%
Specificity	98.1%	96.9% - 98.9%
Specimens Positive for <i>Giardia</i> by microscopy/ Specimens Detected by the <i>TRI-COMBO PARASITE SCREEN</i>		12/13
Specimens Positive for <i>E. histolytica</i> by microscopy/ Specimens Detected by the <i>TRI-COMBO PARASITE SCREEN</i>		1/1

* The fourteen *TRI-COMBO PARASITE SCREEN* positives that were microscopy negative were confirmed to be positive for *Giardia* with an alternate FDA-cleared antigen test or by PCR with sequencing.

Retrospective Study

Testing consisted of 96 archived specimens previously collected and frozen from one clinical site. The frozen samples are included in the bank based on being characterized as microscopy and PCR positive. The specimens were collected from an *E. histolytica* endemic area and contained specimens also positives for *Giardia* and *Cryptosporidium*. Table 2 summarizes the performance observed. Retrospective testing was also read visually and performance was not significantly different from spectrophotometric readings.

Table 2. Summary of retrospective clinical performance comparing the TRI-COMBO PARASITE SCREEN test to Microscopy and PCR

TRI-COMBO PARASITE SCREEN (N = 96)	Microscopy and PCR	
	Positive	Negative
Positive	85	0
Negative	5	6
		95% Confidence Limits
Sensitivity	94.4%	87.7% - 97.9%
Specificity	100%	61.0% - 100%
Specimens Positive for <i>Giardia</i> / Specimens Detected by the TRI-COMBO PARASITE SCREEN		41/41
Specimens Positive for <i>Cryptosporidium</i> / Specimens Detected by the TRI-COMBO PARASITE SCREEN		27/30
Specimens Positive for <i>E. histolytica</i> / Specimens Detected by the TRI-COMBO PARASITE SCREEN		28/30

Note: Eight specimens were dual positive for *Giardia* and *E. histolytica* by the Microscopy and PCR and tested positive in the TRI-COMBO PARASITE SCREEN test. Three specimens were dual positive for *Giardia* and *Cryptosporidium* by Microscopy and PCR and tested positive in the TRI-COMBO PARASITE SCREEN test.

The prospective study results were analyzed by considering composite results from multiple tests that consisted of light microscopy, molecular testing consisting of a commercial FDA-cleared device and PCR with sequencing for the identification of *Giardia* spp., *Cryptosporidium* spp., in addition to identification and subspeciation of *E. histolytica*. This testing was mainly done because identification of *E. histolytica* organisms cannot be determined solely by microscopy because it is morphologically indistinguishable from the non-pathogenic *E. dispar*. Use of an alternate molecular testing is needed to confirm Entamoeba speciation. The molecular testing algorithm used provides a comparator method that is highly sensitive at the detection of *Giardia* spp., *Cryptosporidium* spp. and *E. histolytica*. The performance is summarized in table 3 and is presented as positive percent agreement and negative percent agreement.

Table 3. Summary of prospective clinical performance comparing the TRI-COMBO PARASITE SCREEN test to microscopy and molecular testing

TRI-COMBO PARASITE SCREEN (N = 754)	Microscopy and Molecular testing	
	Positive	Negative
Positive	18	9*
Negative	11**	716
		95% Confidence Limits
Positive Percent Agreement	62.1%	44.0% - 77.3%
Negative Percent Agreement	98.8%	97.7% - 99.4%

Specimens Positive for <i>Giardia</i> / Specimens Detected by the TRI-COMBO PARASITE SCREEN	17/24
Specimens Positive for <i>E. histolytica</i> / Specimens Detected by the TRI-COMBO PARASITE SCREEN	1/5
Specimens Positive for <i>Cryptosporidium</i> / Specimens Detected by the TRI-COMBO PARASITE SCREEN	0/0

* These nine specimens were tested with an alternate FDA cleared antigen test resulting in 9/9 *Giardia* determined to be antigen positive.

** These eleven specimens were tested with an alternate FDA cleared antigen test resulting in 6/7 *Giardia* and 4/4 *E. histolytica* were determined to be antigen negative.

REPRODUCIBILITY

The reproducibility of the TRI-COMBO PARASITE SCREEN test was determined using 20 human fecal specimens coded to prevent their identification during testing. Testing was performed at 2 independent laboratories and on-site at TECHLAB, Inc. The samples were tested twice a day over a 5-day period by multiple technicians at each site using 2 different kit lots. Positive and negative controls were run with each panel of the masked samples. The results from each laboratory were submitted to TECHLAB, Inc. and compared with in-house results. The results were consistent among the different locations and exhibited a correlation of 99.9%. The samples produced the expected results 99.9% of the time.

CROSS REACTIVITY

The TRI-COMBO PARASITE SCREEN test was evaluated for cross reactivity with the bacterial and viral strains listed below. None of the strains were shown to interfere with the performance of the TRI-COMBO PARASITE SCREEN test.

<i>Aeromonas hydrophila</i>	<i>Enterococcus faecalis</i>	<i>Shigella sonnei</i>
<i>Bacillus cereus</i>	<i>Escherichia coli</i>	<i>Staphylococcus aureus</i>
<i>Bacillus subtilis</i>	<i>Escherichia coli</i> O157:H7	<i>Staphylococcus aureus</i>
<i>Bacteroides fragilis</i>	<i>Escherichia coli</i> EIEC	(Cowan's)
<i>Campylobacter coli</i>	<i>Escherichia coli</i> EPEC	<i>Staphylococcus epidermidis</i>
<i>Campylobacter fetus</i>	<i>Escherichia coli</i> ETEC	<i>Vibrio parahaemolyticus</i>
<i>Campylobacter jejuni</i>	<i>Klebsiella pneumoniae</i>	<i>Yersinia enterocolitica</i>
<i>Candida albicans</i>	<i>Salmonella typhimurium</i>	
<i>Clostridium bifermentans</i>	<i>Shigella dysenteriae</i>	
<i>Clostridium difficile</i>	<i>Shigella flexneri</i>	

Calicivirus	Human Coxsackievirus B2, B3, B4, B5
Cytomegalovirus	Human Echovirus 9
Echovirus 11, 18, 33	Human Enterovirus 68, 69, 70, 71
Human Adenovirus 1, 2, 3, 5, 40, 41	Human parechovirus 1 [Echovirus 22]
Human Coronavirus	Human Rotavirus

Additionally, the TRI-COMBO PARASITE SCREEN test was run on fecal specimens documented to be positive for other parasites by microscopy. The number in parentheses is the quantity of each organism found in the clinical specimens. No cross reactivity was seen with the following organisms.

<i>Ascaris lumbricoides</i> and with eggs (22)	<i>Entamoeba moshkovskii</i> (3)
<i>Blastocystis hominis</i> (11)	<i>Iodamoeba bütschlii</i> (10)
<i>Entamoeba bangladeshi</i> (3)	<i>Trichuris trichiura</i> eggs (12)
<i>Entamoeba coli</i> (16)	

Cross reactivity with Norovirus is unknown because it was not tested in analytical studies. However, Norovirus GI/GII was identified in 34 clinical specimens and Enterotoxigenic *E. Coli* - ETEC LT/ST was identified in 107 clinical specimens, using an FDA cleared multiplex NAAT assay during clinical testing and no cross reactivity was found using the *TRI-COMBO PARASITE SCREEN* in those samples.

Strain Specific Study

Due to the similarity in morphology between pathogenic and non-pathogenic *Entamoeba* species, 3 specimens identified by PCR as positive for non-pathogenic *Entamoeba moshkovskii* and 3 positive for non-pathogenic *Entamoeba bangladeshi* were evaluated using the *TRI-COMBO PARASITE SCREEN* test. These 6 specimens tested negative in the *TRI-COMBO PARASITE SCREEN* test.

INTERFERING SUBSTANCES (U.S. FORMULATION)

The following substances had no effect on positive or negative *TRI-COMBO PARASITE SCREEN* test results analyzed at the concentrations indicated: Barium sulfate (5% w/v), Benzalkonium Chloride (1% w/v), Ciprofloxacin (0.25% w/v), Ethanol (1% w/v), Hog gastric mucin (3.5% w/v), Human blood (40% v/v), Hydrocortisone (1% w/v), Imodium® (5% v/v), Kaopectate® (5% v/v), Leukocytes (0.05% w/v), Maalox® Advanced (5% v/v), Mesalazine (10% w/v), Metronidazole (0.25% w/v), Mineral Oil (10% w/v), Mylanta® (4.2 mg/mL), Naproxen Sodium (5% w/v), Nonoxonyl-9 (1% w/v), Nystatin (1% w/v), Palmitic Acid/Fecal Fat (40% w/v), Pepto-Bismol® (5% v/v), Phenylephrine (1% w/v), Polyethylene glycol 3350 (10% w/v), Prilosec OTC® (5 µg/mL), Sennosides (1% w/v), Simethicone (10% w/v), Stearic Acid/Fecal Fat (40% w/v), Tagamet® (5 µg/mL), TUMS® (50 µg/mL), Human Urine (5% v/v), and Vancomycin (0.25% w/v).

PRECISION – INTRA-ASSAY

For the determination of intra-assay performance, 24 fecal samples were analyzed by the *TRI-COMBO PARASITE SCREEN* test. The samples included 2 negative, 2 high negative, 2 low positive, and 2 moderate positive samples for each analyte. Each specimen was assayed a total of five times using two different kit lots. Positive specimens consistently tested positive and negative specimens consistently tested negative. High negative samples tested within 95% agreement of each other for all three analytes.

PRECISION – INTER-ASSAY

For the determination of inter-assay performance, 24 fecal samples were analyzed by the *TRI-COMBO PARASITE SCREEN* test. The samples included 2 negative, 2 high negative, 2 low positive, and 2 moderate positive samples for each analyte. The samples were tested twice a day by multiple technicians over a 12-day period using 2 different kit lots. All positive samples remained positive and all negative samples remained negative.

ANALYTICAL SENSITIVITY

The analytical sensitivity of the test was determined by using purified *Giardia* cysts, *Cryptosporidium* oocysts, and *E. histolytica* pathogenic zymodemes in a sample matrix. The concentration of *Giardia* cysts, *Cryptosporidium* oocysts, and *E. histolytica* pathogenic zymodemes in fecal matrix at which specimens were positive by the *TRI-COMBO PARASITE SCREEN* test 95% of the time was the assay limit-of-detection (LoD). Test results determined the LoD for the assay to be 8450 cysts/mL of feces for *Giardia* (equivalent to 169 cysts detected per test), 47,962 oocysts/mL of feces for *Cryptosporidium* (equivalent to 959 oocysts detected per test), and 1,676 PZs/mL (equivalent to 34 PZs detected per test).

For fecal matrix/Cary Blair, the LoD for the assay was determined to be 34,155 cysts/mL of feces for *Giardia* (equivalent to 427 cysts detected per test), 99,456 oocysts/mL of feces for *Cryptosporidium* (equivalent to 1243 oocysts detected per test), and 4655 PZs/mL (equivalent to 58 PZs detected per test). For fecal matrix/C&S, the LoD for the assay was determined to be 37,095 cysts/mL of feces for *Giardia* (equivalent to 464 cysts detected

per test), 122,299 oocysts/mL of feces for *Cryptosporidium* (equivalent to 1529 oocysts detected per test), and 3948 PZs/mL (equivalent to 49 PZs detected per test).

Because the *TRI-COMBO PARASITE SCREEN* test detects soluble antigen in fecal specimens in addition to cysts, oocysts, and trophozoites, this LoD study represents an estimate of analytical sensitivity.

FRESH VERSUS FROZEN SAMPLES

The effect of long term frozen specimen storage on antigen stability was evaluated. For the analysis, a total of 39 fecal samples were tested with the *TRI-COMBO PARASITE SCREEN* test. The samples consisted of negative fecal samples, high negative fecal samples, low positive fecal samples, moderate positive fecal samples, and high positive fecal samples. Samples were prepared by spiking a negative fecal pool with *Giardia* cysts, *Cryptosporidium* oocysts, or *E. histolytica* pathogenic zymodemes at respective concentrations and stored at $\leq -10^{\circ}\text{C}$. These samples were tested at 0, 1, 4 and 8 weeks. No conversion of positive-to-negative or negative-to-positive was observed in any of the samples at the specified time points.

PROZONE

To ensure that a high concentration of analyte does not interfere with a positive reaction in *TRI-COMBO PARASITE SCREEN* test, high samples were prepared by spiking a negative fecal pool with *Giardia* cysts, *Cryptosporidium* oocysts, or *E. histolytica* pathogenic zymodemes, and then tested. A total of 5 different dilutions of each, up to and including the clinically observed high concentration, were prepared and tested in triplicate. The results demonstrated that there was no overall prozone affect, that elevated levels of analyte did not affect the detection of each organism.

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© 2023 TECHLAB, Inc. All rights reserved.
The TECHLAB Logo and TECHLAB are trademarks of TECHLAB, Inc.