# ***Clostridium difficile* Toxin/Antitoxin Kit**

A Toxin/Antitoxin Kit for the Detection of  
*Clostridium difficile* Toxin in Clinical Specimens  
Catalog No. T5000

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International Symbol Key:

**REF** Catalog Number

**IVD** *In Vitro* Diagnostic Medical Device

**LOT** Lot Information

Contains sufficient reagents

Temperature Limitation

Use By/Expiration Date

**CE** CE Symbol

Caution, consult

## **C. difficile Toxin/Antitoxin Kit**

### **INTENDED USE**

The *Clostridium difficile* Toxin/Antitoxin Kit is intended for use in conjunction with the tissue culture cytotoxicity assay for the confirmation of *Clostridium difficile* toxin in patient specimens. The kit includes a toxin control reagent and a specific *C. difficile* antitoxin for use in a tissue culture assay. The test is to be used as an aid in the diagnosis of *C. difficile* disease and results should be considered in conjunction with the patient history.

FOR *IN VITRO* DIAGNOSTIC USE.

### **EXPLANATION**

*Clostridium difficile* is a major nosocomial pathogen that causes virtually all cases of pseudomembranous colitis and many cases of antibiotic-associated diarrhea (1). Its role as a major enteropathogen became apparent in the late 1970's when the results from a number of laboratories showed that (i) *C. difficile* could be isolated from patients with pseudomembranous colitis, (ii) the isolates produced the same cytotoxic (i.e., rounding) effect as that seen in fecal specimens from the patients, and (iii) the isolates produced a similar type of gastrointestinal disease in hamsters treated with antibiotics. The disease caused by *C. difficile* is due to the toxins that it produces. There are two toxins, toxin A and toxin B, that are involved in the disease. Toxin A is a potent enterotoxin and is lethal for experimental animals and cytotoxic for tissue culture cells. Toxin B is also lethal and cytotoxic for cells. Toxin B is often referred to as the cytotoxin because it is much more cytotoxic than toxin A on most tissue culture cells. It is this potent cytotoxic or rounding activity that has resulted in the use of toxin B as a marker for *C. difficile* toxin in fecal specimens. However, because of the enterotoxic activity of toxin A in the intestine, it is believed that toxin A causes most of the symptoms associated with the disease (2,3). This does not rule out a role for toxin B since there is evidence suggesting a synergistic action between the toxins.

The actual discovery of *C. difficile* toxin in fecal specimens from patients with pseudomembranous colitis (PMC) was due to the use of the tissue culture assay. In the mid-to-late 1970's, clinical investigators noted the presence of a cytotoxic or rounding activity in patients with PMC and subsequently discovered that the cytotoxic effect was neutralized by gas gangrene antitoxin. Further research led to the discovery that *C. sordellii* antitoxin neutralized the activity, leading scientists to suspect *C. sordellii* as the etiologic agent. However, carefully performed studies that followed these discoveries identified *C. difficile* as the cause of the disease. The neutralizing activity of the *C. sordellii* antitoxin was ultimately shown to occur because certain toxigenic strains of this organism produce toxins (toxin HT and toxin LT) that are very similar to toxins A and B of *C. difficile*. Thus, *C. sordellii* antitoxin contains cross reacting antibodies that neutralize the toxins of *C. difficile*.

The toxins of *C. difficile* are the largest bacterial toxins known. Toxins A and B have molecular weights of 308,000 and 279,000, respectively. Both toxins have a series of repeating units at the COOH-terminus of the molecule. The repeating units on the toxin A molecule are the binding portion that recognizes the carbohydrate receptors. The toxins have about 45% overall similarity in their amino acid sequence (4,5). This observation and the fact that the toxins elicit the same effect on tissue culture cells suggests that the toxins have the same mechanism of action. The genes for the toxins are located approximately 1 kb apart with a small open-reading frame located between them. Most toxigenic strains of *C. difficile* either produce both toxins or neither toxin, although recently, toxin A negative/toxin B positive strains have been identified (6-8).

There are four methods currently used for the detection of *C. difficile* and its toxins in fecal specimens. These include (i) isolation of the organism, (ii) antigen detection, (iii) tissue culture assay, and (iv) ELISA for the detection of toxin. Isolation of the organism and antigen detection are useful for detecting the presence of the organism in fecal

specimens. However, neither of these tests demonstrates the presence of toxin. This has to be done using either the tissue culture test or an ELISA specific for the toxins. The tissue culture procedure has proved to be very useful because of the high activity of *C. difficile* toxin B in the assay. As little as one picogram of the toxin is sufficient to cause a positive reaction. Neutralization of the cytotoxic activity by *C. difficile* antitoxin confirms the presence of *C. difficile* toxin. The assay is to be used as an aid in the diagnosis of *C. difficile* disease and test results must be considered in conjunction with the patient history. This is due to the fact that some adults are asymptomatic carriers and can have detectable levels of toxin in their feces but have no clinical symptoms. In addition, up to 50% of healthy infants can have *C. difficile* toxin in their feces but show no signs of illness. Therefore, it is important that test results be considered in light of the patient history.

### PRINCIPLE OF THE TEST

Two reagents are provided in the *Clostridium difficile* Toxin/Antitoxin Kit. One is the toxin control reagent to be used as the positive control material. The other is specific *C. difficile* antitoxin and is used to confirm the presence of *C. difficile* toxin by neutralizing the cytotoxic activity. For the assay, dilutions of the toxin control reagent are prepared and mixed either with buffer or diluted antitoxin. In addition, an extract of the clinical specimen is prepared and dilutions of the extract are prepared. The extract dilutions are then mixed either with buffer or diluted antitoxin. After an incubation period, an aliquot of each mixture is transferred to a microtiter well or tube containing a monolayer of tissue culture cells. In the presence of *C. difficile* toxin, the cells show a characteristic rounding effect. The toxin control reagent is used to demonstrate this effect and to show that the effect is neutralized by the *C. difficile* antitoxin. If the dilutions of specimen extract show a rounding effect and the effect is neutralized by the *C. difficile* antitoxin, the presence of *C. difficile* toxin is confirmed.

### MATERIALS PROVIDED

1. **T3000** Lyophilized *C. difficile* toxin control reagent - 1 mL (cytotoxic titer of  $>10^4$ )
2. **T1000** Lyophilized *C. difficile* antitoxin - 3 mL (prepared in goats)

### MATERIALS NOT PROVIDED

Test tubes for preparing dilutions	Pipettes/pipettor
Centrifuge	Sterile filters for preparing extracts
Tissue culture cells	Phosphate buffered saline (PBS)
Incubator	Microscope for observing cells

### PRECAUTIONS

1. Use fecal specimens within 24 hours of collection to obtain optimal results. Specimens that are frozen ( $-20^{\circ}\text{C}$  or lower) may lose activity due to freezing and thawing. AVOID REPEATED CYCLES OF FREEZING AND THAWING THE SPECIMEN. THIS MAY INACTIVATE THE TOXIN.
2. All specimens should be considered as potentially infectious and should be handled appropriately.
3. The toxin control reagent contains biologically active toxin and should be treated with caution. The antitoxin does not pose any known hazard.
4. Do not use the reagents beyond their expiration dates.
5. Keep reagents stored according to the manufacturer's recommendations. This prolongs the shelf-life of the reagents.
6. Place reagents back in refrigerator once the procedure has been completed.

### COLLECTION AND HANDLING OF FECAL SPECIMENS

**NOTE:** Standard collection and handling procedures used in-house are appropriate. Specimens should be transported as soon as possible and stored between  $2^{\circ}$  and  $8^{\circ}\text{C}$ . Whenever possible, test fecal specimens that are less than 24 hours old. Store

specimens at -20°C, or lower, if the test cannot be performed within 48 hours. FREEZING AND THAWING OF THE SPECIMEN, ESPECIALLY MULTIPLE TIMES, MAY RESULT IN LOSS OF ACTIVITY DUE TO DEGRADATION OF THE TOXIN. Make sure that specimens are thoroughly mixed (vortexed) prior to performing the assay. This includes complete mixing of the specimen prior to preparing the fecal extract.

#### PROCEDURE

**NOTE:** The procedure described below may be used for detecting *C. difficile* toxin in fecal specimens with the *Clostridium difficile* Toxin/Antitoxin Kit (8). Alternatively, the kit may be used according to other modified in-house procedures. It is important, however, that the toxin control reagent and the *C. difficile* antitoxin be used according to the instructions. Once the toxin control reagent and the antitoxin reagent have been reconstituted, aliquot the reagents in small volumes (50 to 100 µL) and freeze the undiluted aliquots. Airtight vials should be used to prevent dehydration of the aliquotted toxin and antitoxin. Do not freeze the **diluted** toxin control reagent or the diluted antitoxin reagent. Thaw and use undiluted aliquots of each reagent as needed.

#### PREPARATION OF CONTROLS AND FECAL FILTRATES

1. TOXIN CONTROL REAGENT. Add 1 mL sterile distilled water to lyophilized toxin control reagent. Aliquot small volumes (50 to 100 µL) of the undiluted reagent and store the unused aliquots at -20°C or lower in a freezer that is not self-defrosting. To further minimize dehydration of the undiluted aliquots, airtight vials should be used. To use the toxin control reagent, prepare tenfold dilutions ( $10^1$  to  $10^5$ ) in sterile PBS.
2. ANTITOXIN REAGENT. Add 3 mL sterile distilled water to the lyophilized antitoxin reagent. Aliquot small volumes (50 to 100 µL) of the undiluted reagent and store the unused aliquots at -20°C or lower. To use the antitoxin reagent, prepare a 1:25 dilution with sterile PBS.
3. FECAL FILTRATE. Use undiluted fecal specimen for watery samples or dilute 1:1 with sterile PBS. Centrifuge the specimen at 10,000 x g for 10 minutes and filter the supernatant fluid through a 0.45 µm membrane. Prepare tenfold dilutions ( $10^1$  to  $10^5$ ) in sterile PBS.

#### PERFORMING THE ASSAY IN MICROTITER WELLS

(Mix all reagents thoroughly before using or dispensing. Mixtures may be prepared in sterile microtiter wells or test tubes)

1. Mix 1 part (0.1 mL) of each dilution of the toxin control reagent with 1 part (0.1 mL) of sterile PBS. This series represents the positive control and demonstrates the characteristic cell rounding caused by *C. difficile* toxin.
2. Mix 1 part (0.1 mL) of each dilution of toxin control with 1 part (0.1 mL) of the diluted antitoxin. This series represents the antitoxin control and demonstrates the neutralization of *C. difficile* toxin by the antitoxin.
3. Mix 1 part (0.1 mL) of each dilution of fecal filtrate with 1 part (0.1 mL) of sterile PBS. This series will detect the presence of cytotoxic activity (i.e., cell rounding) in the fecal specimen.
4. Mix 1 part (0.1 mL) of each dilution of fecal filtrate with 1 part (0.1 mL) of diluted antitoxin. This series will confirm the presence of *C. difficile* toxin if the cytotoxic activity is neutralized by the antitoxin.
5. Incubate the mixtures for 30 minutes at room temperature.
6. Using a micropipette and new tips for each dilution, add 20 µL of each test mixture to microtiter wells containing cell monolayers. Incubate the wells overnight (18 hours) at  $37^\circ \pm 2^\circ\text{C}$  under conditions that are appropriate for the tissue cultured cells used in the assay. Many cell lines require  $\text{CO}_2$  and a humid atmosphere. Other cell lines perform adequately in HEPES-buffered media and do not require  $\text{CO}_2$ . Observe the condition of the cells. A positive reaction is noted by the rounding of the cells.

## PERFORMING THE ASSAY IN CULTURE TUBES

(Mix all reagents thoroughly before using or dispensing. Mixtures should be prepared in sterile test tubes)

1. Mix 1 part (0.2 mL) of each dilution of the toxin control reagent with 1 part (0.2 mL) of sterile PBS. This series represents the positive control and demonstrates the characteristic cell rounding caused by *C. difficile* toxin.
2. Mix 1 part (0.2 mL) of each dilution of toxin control with 1 part (0.2 mL) of the diluted antitoxin. This series represents the antitoxin control and demonstrates the neutralization of *C. difficile* toxin by the antitoxin.
3. Mix 1 part (0.2 mL) of each dilution of fecal filtrate with 1 part (0.2 mL) of sterile PBS. This series will detect the presence of cytotoxic activity (i.e., cell rounding) in the fecal specimen.
4. Mix 1 part (0.2 mL) of each dilution of fecal filtrate with 1 part (0.2 mL) of diluted antitoxin. This series will confirm the presence of *C. difficile* toxin if the cytotoxic activity is neutralized by the antitoxin.
5. Incubate the mixtures for 30 minutes at room temperature.
6. Using a sterile pipette for each control and specimen filtrate, add sufficient volumes of each test mixture to culture tubes containing cell monolayers to achieve a 1/10 dilution. For example, add 0.2 mL of mixture to a tube containing 1.8 mL. Replace caps on inoculated culture tubes. Incubate the wells overnight (18 hours) at  $37^{\circ} \pm 2^{\circ}\text{C}$  under conditions that are appropriate for the tissue cultured cells used in the assay. Many cell lines require  $\text{CO}_2$  and a humid atmosphere. Other cell lines perform adequately in HEPES-buffered media and do not require  $\text{CO}_2$ . Observe the condition of the cells. A positive reaction is noted by the rounding of the cells.

## QUALITY CONTROL

A positive and negative control must be run with each series of test specimens. This includes the toxin control reagent mixed with buffer as a positive control and the toxin control reagent mixed with antitoxin to demonstrate neutralization of the cytotoxic activity. The positive control should show the typical rounding effect characteristic of *C. difficile* toxin. The test results are not valid unless these performance characteristics are met. If the performance characteristics are not met, please call Technical Services, and do not report test results.

Upon arrival, the kit should be inspected to ensure that components arrived in satisfactory condition. Each vial should be examined for the presence of a lyophilized powder. Test results along with control results should be recorded and reported according to in-house procedures and this information should be stored according to in-house procedures for future reference.

## INTERPRETATION AND REPORTING RESULTS

Dilutions of toxin control reagent mixed with buffer should show cell rounding out to the last or next-to-last dilution. Dilutions of toxin control reagent mixed with antitoxin should be neutralized through the  $10^4$  dilution +/- one dilution. In mixtures containing dilutions of fecal filtrate plus buffer, cell rounding indicates the presence of cytotoxic activity. If this activity is neutralized in dilutions of fecal filtrate plus antitoxin, the presence of *C. difficile* toxin is confirmed.

Specimens that contain cytotoxic activity which is neutralized by *C. difficile* antitoxin should be reported as "Positive for *C. difficile* toxin". Specimens that do not contain cytotoxic activity should be reported as "Negative for *C. difficile* toxin". Specimens that cause cell rounding which is not neutralized by *C. difficile* antitoxin should be reported as "Cytotoxic activity is present but is not characteristic of *C. difficile* toxin".

## LIMITATIONS OF THE PROCEDURE

1. The *Clostridium difficile* Toxin/Antitoxin Kit is used to detect *C. difficile* toxin in fecal specimens. Detection of the toxin is dependent on the biological activity of the

toxin. Specimens must be handled properly and stored as recommended to minimize inactivation of the toxin.

2. Some specimens may give reactions that are inconclusive (e.g., partial rounding of the cells, stretching instead of rounding). Under these conditions, the specimen should be retested or a fresh specimen should be tested. Additional tests that may be used in conjunction with the *Clostridium difficile* Toxin/Antitoxin Kit include the ELISA for toxin, isolation of the organism on selective media, or antigen detection.
3. Certain isolates of *Clostridium sordellii* produce the same type of rounding on tissue culture cells as toxigenic *Clostridium difficile*, due to the similarities of the *C. sordellii* toxins and the *C. difficile* toxins. *Clostridium sordellii* has not been detected in patients with antibiotic associated-diarrhea and colitis. It is unlikely that *C. sordellii* will be present in human fecal specimens.

## PERFORMANCE CHARACTERISTICS

In a study performed in-house at TECHLAB<sup>®</sup>, Inc., of 100 specimens submitted from a previous *C. difficile* study, 87 were positive for cytotoxic activity. All 87 were confirmed positive for TECHLAB<sup>®</sup> *C. difficile* toxin using *C. sordellii* antitoxin from the U.S. Bureau of Biologics and by the *C. difficile* antitoxin. Therefore, the *C. difficile* antitoxin showed a sensitivity and specificity of 100% compared to *C. sordellii* antitoxin. The predictive positive and negative values and the correlation with the use of *C. sordellii* antitoxin also were 100%.

In another study at a Midwestern VA medical center, of 85 specimens submitted for *C. difficile* testing, 20 were confirmed positive for *C. difficile* toxin using *C. sordellii* antitoxin whereas 19 were confirmed positive by the *Clostridium difficile* Toxin/Antitoxin Kit. The one specimen that was positive with the *C. sordellii* antitoxin and negative with the *Clostridium difficile* Toxin/Antitoxin Kit was only weakly positive. In this study, the sensitivity and specificity were 95% and 100%, respectively. The predictive positive and negative values were 100% and 98.5%, respectively, and the correlation was 98.8%.

Overall, the sensitivity and specificity of the *Clostridium difficile* Toxin/Antitoxin Kit was 99.2% and 100%, respectively, compared with *C. sordellii* antitoxin. The predictive positive and negative values were 100% and 98.8%, respectively, and the overall correlation was 99.5%.

## EXPECTED VALUES

In a study in which specimens were analyzed in-house at TECHLAB<sup>®</sup>, Inc., the prevalence of a positive tissue culture assay in specimens actually submitted to the clinical laboratory was about 20%. In the study at the Midwestern VA medical center, the prevalence was 23.5%. The prevalence will vary from location to location and hospitals may experience rates lower or higher than those observed at this particular site. *Clostridium difficile* disease is primarily a nosocomial disease of elderly patients and hospitals that have higher numbers of elderly patients may experience higher rates. It is important to consider any test results in conjunction with clinical symptoms since some healthy adults and large numbers of healthy infants (up to 50%) will be positive for *C. difficile* toxin. In addition, *C. difficile* carriage rates of 22% to 32% have been reported in cystic fibrosis patients (10-12).

## CROSS-REACTIVITY

Toxic filtrates from *C. chauveoi*, *C. fallax*, *C. perfringens*, *C. septicum*, and *C. spiroforme* were not neutralized by the TECHLAB<sup>®</sup> *C. difficile* antitoxin. In addition, alpha-toxin from *C. perfringens*, collagenase and clostripain from *C. histolyticum*, cholera toxin, diphtheria toxin, shiga-like toxin, and staphylococcal enterotoxin were not neutralized by the antitoxin. The only organism that is known to cross-react in the *Clostridium difficile* Toxin/Antitoxin Kit is *C. sordellii*. Certain toxigenic isolates of *C. sordellii* produce toxins that are very similar to the toxins of *C. difficile*. The hemorrhagic toxin (toxin HT) of *C. sordellii* is very similar to toxin A whereas the lethal toxin

(toxin LT) is very similar to toxin B. Antibodies against the *C. difficile* toxins neutralize the toxins HT and LT of *C. sordellii*, and antibodies against toxins HT and LT neutralize toxins A and B. It should be emphasized that *C. sordellii* has not been implicated in antibiotic-associated diarrhea and colitis.

#### REPRODUCIBILITY AND PRECISION

Identical set of five fecal specimens were sent from TECHLAB®, Inc. to three different laboratories for analysis with the *Clostridium difficile* Toxin/Antitoxin Kit. One laboratory was located on the east coast and the other two laboratories were located in the Midwest. The test results from each laboratory were consistent with the test results obtained in-house at TECHLAB®, Inc. In an in-house study, ten positive fecal specimens were assayed at time 0, 24, 48, and 72 hours using the *Clostridium difficile* Toxin/Antitoxin Kit. These specimens were shown to be consistently positive for *C. difficile* toxin at each of the time intervals, as determined by neutralization of the cytotoxic activity. It should be noted that the specimens were stored at 4°C during the testing period and that in some fecal specimens, the cytotoxic activity is not stable and will decrease upon incubation at 4°C. In additional in-house studies, three different lots of the positive control reagent were tested in duplicate and consistently gave cytotoxic titers of 10<sup>4</sup> or higher. Further evaluation demonstrated that five different lots of the *C. difficile* antitoxin consistently neutralized the dilutions of positive control reagent.

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